Document Description:	ABPM-7100 PROTOCOL DOES NOT RUN	Document Number: Version:	80021507 A
	Welch Allyn™		
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.			

Customer Service Bulletin						
Product:	ABPM-7100 when used with CPWS 1.6.6		Date: 2016-05-06			
Subject:	ABPM-7100; Programmed Protocol Does Not Run					
HW Version(s) Affected:	All ABPM-7100 Devices	SW Version(s) Affected:	CPWS 1.6.6			
Serial Numbers Affected:	ABPM-7100 Devices: All CPWS 1.6.6: All	Lot or Date Code Affected:	N/A			
Classification: Informational Only						
Distribution:	□ Customer Care	□ Product Service □	⊠ Field Service			
		□ Customers	☐ Company Confidential			
Training Require	ed: □ Yes ⊠ No					
It has been discovered that, when creating a protocol for the ABPM-7100, if the protocol contains a single <i>Time interval</i> and the <i>Starting time</i> is adjusted from the default of midnight (0:00), the protocol does not run on the device and no automatic blood pressures will be taken. Note: When creating a protocol that contains a single <i>Time interval</i> , the <i>Starting time</i> must not be changed from the default of 0:00. • Examples of protocols that will not run on the ABPM-7100 (Single <i>Time interval</i> with a <i>Starting time</i> set to something other than 0:00): Protocol Designer Protocol 20:090 Protocol 30:090 P						


• Examples of protocols that run properly as programmed on the ABPM-7100 (Single *Time interval* with a *Starting time* set to 0:00, or multiple *Time intervals* with customized *Starting times*):

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.

Protocol 5min Protocol 6min

FCD-1482 Rev K Parent: WI-17206803

Document Description: ABPM-7100 PROTOCOL DOES NOT RUN	Document Number: 80021507 Version: A		
WelchAllyn ™			
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.			



Note: Be sure that all non-working protocols in the Protocol Designer have been edited or removed. While in the Protocol Designer, inspect each protocol to be sure that this instruction has been applied to all affected protocols. Each single *Time interval* protocol needs to be checked to confirm that the *Starting time* is set to 0:00. Single *Time interval* protocols that have a *Starting Time* other than 0:00 need to be edited so that the *Starting time* is set to 0:00, or the protocols deleted and recreated using the default *Starting time* of 0:00.

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2016-05-06	SMB

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.