

## **Patient Monitoring for COVID-19 Outside the ICU**

Role of Hillrom's Patient Vital Signs Monitoring and Diagnostic Cardiology Devices

### **Background and Purpose**

With the overwhelming demand on intensive care units (ICU) due to the COVID-19 pandemic, there is an immediate need to effectively risk stratify patients and identify patients most appropriate for the ICU. For patients treated outside the ICU, appropriate monitoring must be deployed to identify those with increasing severity early.

COVID-19 disease due to the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) presents with a spectrum of clinical features from mild to severe life-threatening disease. Major complications include severe pneumonia, acute respiratory distress syndrome (ARDS), and septic shock.<sup>1</sup> Patients with chronic conditions or comorbidities such as hypertension, diabetes mellitus, and/or cardiovascular disease are at greatest risk of deterioration and death.<sup>2</sup>

The majority of COVID-19 cases are non-critical. In a study of all reported cases in China through February 11, 2020 (N = 44,672), 80.9% of cases were mild.<sup>2</sup> Among hospitalized patients in China, studies have shown 73.9% to 93.9% of cases did not require admission to ICU.<sup>1,3</sup> Data from the US CDC shows 21.7 – 31.4% of COVID-19 cases in the US were hospitalized with only 4.9 – 11.5% admitted to the ICU.<sup>4</sup>

With such a large proportion of patients treated outside the ICU, determining the appropriate monitoring for clinical deterioration is important. The purpose of this document is to show how Hillrom vital signs and diagnostic cardiology products support the monitoring needs of COVID-19 patients outside the ICU.

### **Clinical Application**

#### **Vital Signs Monitoring for COVID-19 Patients Outside the ICU**

Vital signs are important identifiers of early deterioration of COVID-19 patients. According to the World Health Organization (WHO) guidance on COVID-19 pneumonia, "early identification of those with severe illness, such as severe pneumonia, allows for optimized supportive care treatments and safe, rapid referral" to appropriate COVID-19 treatment settings.<sup>5</sup>

A retrospective study looking at clinical outcomes and risk factors for COVID-19 showed the importance of vital sign derangement on mortality. Specifically, admitting vitals signs such as respiratory rate > 24 was present in 63% of non survivors

vs. 16% in survivors and pulse rate  $\geq 125$  /min in 4% of non-survivors vs 0% in survivors. Zhou et al., 2020 found an 8.89 (95% confidence interval 4.34 – 18.19) odds ratio of death in patients with a respiratory rate  $> 24$ .<sup>6</sup>

Sun et al. 2020 demonstrated vital signs play a significant role in mortality reduction. Using a protocol with vital signs-based triggers to escalate care, researchers noted a COVID-19 cure rate of 96.67%: a better rate than data at the national level in China. The mortality rate in this population was 3.3% as compared to 4.34% in Hubei province. This example demonstrates the benefit of establishing an early identification system of deteriorating patients before they progress to critical illness severity.

For non-critical COVID-19 patients, in the Sun et al. study, the determining factors to transfer patients to critical care management were these vital signs:

- Respiratory rate  $> 30$  breaths per minute
- SpO<sub>2</sub>  $< 93\%$
- Heart rate  $> 120$  beats per minute.<sup>7</sup>

Interestingly, WHO guidance on the clinical management of respiratory infection in COVID-19 patients uses the same respiratory rate and oxygen saturation cutoffs to identify severe pneumonia.<sup>5</sup> The 2019 American Thoracic Society guideline on community-acquired pneumonia also identifies a respiratory rate  $\geq 30$  as a criterion for severe disease.<sup>8</sup>

High risk patients had their vital signs monitored continuously whereas low risk patients were monitored twice per day. Differentiation between high risk and low risk patients was determined based on several demographic, lab and radiographic findings (see Figure 1 in appendix).<sup>7</sup>


On **Centrella® Smart Bed with Contact-Free Continuous Monitoring**, respiratory rate and heart rate are automatically and continuously monitored when a patient lies in the bed and provides trending data for these vital signs. These trending data provide clinicians with easy access to respiratory rate and heart rate changes over time. The **Welch Allyn® Connex® Vital Signs Monitor** continuously monitors end-tidal CO<sub>2</sub> and oxygen saturation (SpO<sub>2</sub>). These technologies, with alarm thresholds, alert clinicians when there is a significant vital sign change.

For lower risk patients, the **Welch Allyn® Connex® Spot Monitor** can measure vital signs at preset time intervals. (see Figure 3 in Appendix)

**Connex® Monitors** display patients' vitals remotely to decrease infection risk. Continuous vitals monitoring or automated interval measurements can be used with a **Welch Allyn Connex® Central Station** and/or the EMR to display the values. Similarly, **Centrella® Smart Bed with Contact-Free Continuous Monitoring** sends alarms to the

nurse call system to alert clinicians when thresholds have been crossed. These methods of remote monitoring may decrease the frequency of exposure and risk of transmission for healthcare providers and preserve personal protective equipment (PPE). (see Figure 2 see Appendix)

In addition, the **Connex® Monitors** may be configured for automated early warning scoring (e.g. MEWS, NEWS, PEWS) or a facility preferred score. These scores can be visualized at the bedside, integrated into the EMR, as well as the **Connex® Central Station** to support remote monitoring of patients outside the ICU.



Product	Heart Rate	Resp. Rate	Temp.	SpO <sub>2</sub>	Blood Pressure*	ECG	EtCO <sub>2</sub>	EWS
<b>Centrella®</b>	✓	✓						
<b>Connex® Vital Signs Monitor</b>	✓	✓	✓	✓	✓	3/5-Lead	✓	✓
<b>Connex® Central Station</b>	✓	✓	✓	✓	✓	✓	✓	✓
<b>Connex® Spot Monitor</b>	✓	✓	✓	✓	✓			✓
<b>ELI® 280 ECG / ELI® 380 ECG</b>						12-Lead		

\* Non-invasive blood pressure can be automated on a pre-set interval  
 Highlighted cells can be continuously monitored

## Clinical Application

### ECG Monitoring for Cardiac Abnormalities

Cardiac injury is another manifestation of COVID-19 including arrhythmia, myocarditis, and acute coronary syndrome.<sup>9</sup> Arrhythmias, as a complication of COVID-19, occurred in 16.7% of cases in one study.<sup>1</sup> Acute myopericarditis can occur as a complication associated with the disease as seen from a case study in Italy. This presented on ECG analysis minimal diffuse ST segment elevation and associated ST segment depression with T wave inversion in different leads on a 12-lead ECG.<sup>10</sup> A meta-analysis showed that at least 8.0% patients with COVID-19 suffered acute cardiac injury.<sup>11</sup> Another meta-analysis from COVID-19 literature shows an increase in cardiac biomarkers. The authors of the study hypothesize that identification of patients with possible cardiac injury could predict the progression of COVID-19 towards a worse clinical picture.<sup>12</sup>

Additionally, in those admitted with multiple comorbidities, the effect of medication interaction and risk of QT prolongation is of concern. Hospitalized patients come under higher risk for developing Torsades de Pointes than outpatients with the same QT prolonging drugs. Hospitalized patients may include the elderly with underlying heart disease who may also have renal or hepatic dysfunction, electrolyte abnormalities, or bradycardia and to whom drugs may be administered rapidly intravenously. Concomitant administration of drugs inhibiting cytochrome P450 especially imidazole antifungals, macrolide antibiotics, and/or those that can prolong the QT interval, and those that cause electrolyte disturbance increase risk. For any patients on QT-prolonging drugs, surveillance EKGs are recommended before and after initiation of these drugs.<sup>13</sup>

The **Connex® Vital Signs Monitor** with continuous 3- and 5-lead ECG option alarms for critical arrhythmias (ventricular tachycardia, ventricular fibrillation, and asystole) will alert and notify the clinician to immediately assess the patient.

In addition to arrhythmia monitoring, routine 12-lead resting ECG with the **ELI® 280 or ELI® 380** may help identify cardiac injury in COVID-19 patients. For patients on arrhythmogenic drugs (including those that prolong the QT-interval), routine baseline and post-administration 12-lead ECGs are recommended.<sup>13</sup> Furthermore, if signs of acute cardiac disease are identified while cardiac monitoring using a 3 or 5-lead ECG, a resting 12-lead ECG will provide a more detailed view to guide diagnosis and treatment.

The ELI 280 and 380 have benefits in infection prevention as well. With the Wireless Acquisition Module (WAM), there are no cables from the cart to the patient; therefore, there is less opportunity for pathogens to transfer. Additionally, the ELI 280




touch screen and ELI 380 glass keyboard make cleaning easier (see Figure 4 and 5 in Appendix).

## Conclusion

The majority of hospitalized COVID-19 patients will not require treatment in the ICU. For those patients, it is important to monitor vital signs to help identify clinical deterioration. In higher risk patients, continuous monitoring of vital signs outside the ICU can help identify deterioration earlier. Use of smart technology can help overburdened staff closely monitor patients and potentially decrease exposure to the SARS-CoV-2 virus.


## Appendix



**Figure 1:** Risk Stratification to initiate continuous or intermittent monitoring on the general ward




## Continuous Monitoring for Isolation




**Continuous HR & RR monitoring in the bed**

- Centrella® Smart+ Bed with EarlySense® – designed for med-surg
- Continuous HR/RR monitoring with alerts
- Special program for immediate activation\*



**Bedside continuous monitoring**

- SpO2
- etCO2 & respiration rate
- Heart rate
- Interval NIBPs



**Centralized remote monitoring of isolation patients**

- Real-time visibility to all monitored patients
- Alerts at the station

**Benefits:**

- Aligns to WHO's recommendation of dedicated monitor-to-patient for isolation infection prevention and control.<sup>14</sup>
- Limits clinician exposure to infected patients
- SpO2 and etCO2 for accurate respiratory monitoring
- Helps safely and effectively handle increased patient volume – monitor 48 patients per central station
- Customers with existing vital signs monitors can leverage their current devices
- Simple IT ramp-up – software-only central station, and operates on hospital's standard network

**Figure 2:** Continuous Monitoring for isolation with Centrella beds, Connex Monitors and Connex Central Station

## Centralized Interval Monitoring

### Overview:

- Clinician sign-on via active directory
- Patient identification via ADT
- Interval measurements displayed at the Central Station – no clinical interaction required
- To chart vitals in the EMR, clinician must select desired records in the device review tab and save / send.



1. Scan clinician badge to login



2. Scan wristband to ID patient



3. Apply sensor to Patient


4. Begin interval monitoring program on device



4. Interval measurements displayed on central station

**Figure 3:** Connex Monitors in an interval workflow for intermittent monitoring





**Figure 4: Resting ECG**

(A) ELI® 280 ECG and ELI® 380 Resting ECG with the Wireless Acquisition Module WAM™

(B) Typical keyboard versus ELI 380 after cleaning



## References:

1. Wang D, Hu B, Hu C, et al. Clinical Characteristics of 138 Hospitalized Patients With 2019 Novel Coronavirus–Infected Pneumonia in Wuhan, China. *JAMA*. 2020;323(11):1061–1069. doi:10.1001/jama.2020.1585.
2. Novel Coronavirus Pneumonia Emergency Response Epidemiology Team. Vital surveillances: the epidemiological characteristics of an outbreak of 2019 novel coronavirus diseases (COVID-19)—China, 2020. *China CDC Weekly*, 2(8), 113-122.
3. Guan, W. J., Ni, Z. Y., Hu, Y., Liang, W. H., Ou, C. Q., He, J. X., ... & Du, B. (2020). Clinical characteristics of coronavirus disease 2019 in China. *New England Journal of Medicine*.
4. CDC COVID-19 Response Team. *Severe Outcomes Among Patients with Coronavirus Disease 2019 (COVID-19) — United States, February 12–March 16, 2020*. MMWR. Accessed March 30, 2020. [https://www.cdc.gov/mmwr/volumes/69/wr/pdfs/mm6912e2-H.pdf?deliveryName=USCDC\\_921-DM23064](https://www.cdc.gov/mmwr/volumes/69/wr/pdfs/mm6912e2-H.pdf?deliveryName=USCDC_921-DM23064)
5. World Health Organization. (2020). Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: interim guidance, 13 March 2020 (No. WHO/2019-nCoV/clinical/2020.4). World Health Organization.
6. Zhou.F *et.al*. Clinical course and risk factors for mortality of adult inpatients with COVID-19 in Wuhan, China: a retrospective cohort study. *Lancet*. 2020 Mar 28;395(10229):1054-1062. doi: 10.1016/S0140-6736(20)30566-3.
7. Sun Q, Qiu H, Huang M, Yang Y. Lower mortality of COVID-19 by early recognition and intervention: experience from Jiangsu Province. *Ann Intensive Care*. 2020;10(1):33. Published 2020 Mar 18. doi:10.1186/s13613-020-00650-2.
8. Metlay JP, Waterer GW, Long AC, et al. Diagnosis and Treatment of Adults with Community-acquired Pneumonia. An Official Clinical Practice Guideline of the American Thoracic Society and Infectious Diseases Society of America. *Am J Respir Crit Care Med*. 2019;200(7): e45-e67.
9. Driggin E, Madhavan MV, Bikdeli B, et al. Cardiovascular Considerations for Patients, Health Care Workers, and Health Systems During the Coronavirus Disease 2019 (COVID-19) Pandemic. *J Am Coll Cardiol*. 2020.
10. Inciardi RM, Lupi L, Zacccone G, et al. Cardiac Involvement in a Patient with Coronavirus Disease 2019 (COVID-19). *JAMA Cardiol*. Published online March 27, 2020. doi:10.1001/jamacardio.2020.1096.
11. Li, B., Yang, J., Zhao, F. *et al*. Prevalence and impact of cardiovascular metabolic diseases on COVID-19 in China. *Clin Res Cardiol* (2020). <https://doi.org/10.1007/s00392-020-01626-9>.



12. Lippi, C.J. Lavie and F. Sanchis-Gomar, Cardiac troponin I in patients with coronavirus disease 2019 (COVID-19): Evidence from a meta-analysis, *Progress in Cardiovascular Diseases*, <https://doi.org/10.1016/j.pcad.2020.03.001>.
13. Nachimuthu S, Assar MD, Schussler JM. Drug-induced QT interval prolongation: mechanisms and clinical management. *Ther Adv Drug Saf*. 2012;3(5):241–253. doi:10.1177/2042098612454283
14. Infection Prevention and Control of Epidemic- and Pandemic-Prone Acute Respiratory Infections in Health Care. Geneva: World Health Organization; 2014. Annex E, Isolation rooms or areas. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK214341/>