

Competency Skills Checklist

Staff Signature: _____

Trainer's Signature: _____

Date: _____



COMPLETED

Feature/Function Completed**Siderails**

1. Demonstrate brake off audible indicator
2. Identify two different ways a patient can move themselves up in bed
3. Demonstrate use of the care grip positioning handles
4. Identify trend angle indicator
5. Demonstrate the integrated line managers
6. Demonstrate use of patient storage

Controls

7. Demonstrate head up/down
8. Demonstrate SlideGuard® mechanism from patient siderail
9. Demonstrate knee up/down
10. Demonstrate high/low
11. Demonstrate Trendelenburg in a stat situation and how to get bed back in normal position
12. Demonstrate use of enable key
13. Demonstrate use of the bed exit system
14. Demonstrate electric vascular foot mechanism
15. Identify "Bed Not Down" indicators
16. Demonstrate select control on patient pendant
17. Demonstrate use of volume up/down on pendant
18. Demonstrate use of siderail lockout controls
19. Demonstrate chair position
20. Identify battery backup and how to activate it
21. Demonstrate SafeView® alerts (optional feature)
22. Discuss location of digital head-of-bed angle indicator (optional feature)
23. Demonstrate 30° head-of-bed alarm

Scale

24. Identify how to zero the scale
25. Identify what the scale should be zeroed with (per hospital protocol)
26. Demonstrate how to turn scale on/off
27. Demonstrate how to properly weigh a patient
28. Demonstrate how to add or remove items from the bed
29. Demonstrate how to change from lbs to kgs

**Feature/Function Completed****COMPLETED****Patient Position Monitor (Bed Exit System)**

30. Demonstrate on/off of PPM bed exit system

31. Demonstrate volume change on bed exit (scale version)

32. Demonstrate tone change on bed exit (scale version)

33. Demonstrate how to set each mode of bed exit

34. Demonstrate how to reset bed exit system

35. Demonstrate how to turn off bed exit alarm

Other

36. Demonstrate headboard removal

37. Demonstrate utilizing the headboard for CPR

38. Demonstrate IV pole placement

39. Identify restraint strap holders

40. Discuss traction

41. Demonstrate moving/steering bed

42. Demonstrate patient ingress/egress

43. Demonstrate bed cleaning

44. Demonstrate use of brake, neutral and steer functions

45. Demonstrate use of CPR release

46. Demonstrate bed extender

47. Demonstrate IntelliDrive® powered transport (optional feature)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2013 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

ORDER NUMBER 182897 rev 1 E 08/28/2013 ENG – US

Zinc Number US-AC-0513-0208

www.hill-rom.com

USA 800-445-3730

CANADA 800-267-2337

Enhancing outcomes for
patients and their caregivers:

Hill-Rom