

VOLARA SYSTEM



MAXIMIZE PATIENT OUTCOMES

In Acute Care Settings

EXCITING ADVANCEMENTS IN A THERAPY YOU TRUST

Oscillation and Lung Expansion (OLE) therapy has been widely embraced as a proven airway clearance therapy for hospital patients.¹ Now, with the **Volara** System, industry-leading OLE therapy is easier to deliver, more consistently, for a broader range of patients, throughout the hospital—and even at home.

SHORTER CYCLE. FASTER RESULTS.¹

Alternating cycles of Continuous Positive Expiratory Pressure (CPEP) and Continuous High Frequency Oscillation (CHFO) combine with aerosol delivery to help maximize therapy effectiveness. In fact, OLE therapy with the **Volara** System takes just 10 minutes to deliver.

CPEP

Continuous Positive Expiratory Pressure treats and helps prevent pulmonary atelectasis by combining continuous positive pressure with aerosol medication, to help expand and hold the airways open.

CHFO

Continuous High Frequency Oscillation, a pneumatic form of chest therapy, delivers continuous pulses of positive pressure combined with aerosol medication to mobilize and treat retained secretions.

NEBULIZER

Integrated with CPEP and CHFO for maximum treatment efficiency, medication is delivered during therapy to help loosen secretions.

The **Volara** System gives you the versatility to deliver important benefits to patients in multiple hospital settings, with features including:

- Wide therapeutic pressure range: Allows you to treat a broad range of patients
- Programmable therapy presets: Facilitates standardization of therapy across patient groups
- Enhanced device mobility: Enables use throughout the hospital

CLINICALLY PROVEN THERAPY FOR PATIENTS HOSPITAL-WIDE^{2,3}

Tested and proven in acute care settings, OLE therapy has been shown to provide outstanding clinical outcomes.

1.6 DAYS¹ REDUCTION IN HOSPITAL LENGTH OF STAY

2 DAYS¹ REDUCTION IN ICU LENGTH OF STAY

68%¹ REDUCTION IN TIME ON VENTILATOR

12%¹-0%⁴ ASSOCIATED WITH A REDUCTION IN ICU BOUNCE-BACK

VERSATILE DESIGN PRECISION DELIVERY

INTUITIVE INTERFACE FOR SIMPLIFIED OPERATION

Interface shows therapy parameters at a glance

Digital touchscreen allows precise setting of pressures and programs

Preset protocols support precision and continuity of care

One-button operation provides quick access to frequently used presets

Send therapy data to EHR with a standard HL7 interface





MADE FOR PORTABILITY

Lightweight, compact design with ergonomically designed stand for easy portability throughout the hospital

TREAT A BROAD RANGE OF PATIENTS, WITH SEAMLESSLY-INTEGRATED CIRCUIT ADAPTORS.

Patented **Smart-Filter** technology with integrated RFID tag tracks single-patient-use circuit for up to 90 treatments





Dual power source with battery operation and builtin compressor, for operation in a variety of hospital settings



Give your team more control by programming multiple protocols into the device



At-a-glance therapy summary shows you treatment details when session is complete

HELP MORE PATIENTS BREATHE EASIER

The **Volara** System makes OLE therapy practical and accessible to more patients in acute care settings. And when they're discharged, they can continue their OLE therapy at home, using your settings and protocols.

PATIENTS

CONDITIONS

DEPARTMENTS

Acute Chronic Trauma Post-op Geriatric Bariatric Pediatric Pneumonia Bronchiectasis Post-Op Cystic fibrosis COPD Spinal cord inju Trauma Stroke Neuro Transplant Pediatrics Burn

Bronchiectasis Emergency Cystic fibrosis Department Spinal cord injury Gen Med/Surgery Stroke Pediatric

















FOR FURTHER INFORMATION ABOUT PRODUCTS AND SERVICES

Please contact your Baxter sales representative or call Baxter customer service at **1-800-426-4224**.

References

- 1. Huynh TT, Liesching TN, Cereda M, Lei Y, Frazer MJ, Nahouraii MR, et al. Efficacy of Oscillation and Lung Expansion in Reducing Postoperative Pulmonary Complication, *JACS* (2019)
- 2. Clinical Network Survey. September 10, 2018. n=44. Data on file.
- 3. Hillrom. (2022). OLE: A Guide to Clinical Outcomes (206762 rev 3). St. Paul, MN.
- 4. Nyland B, Spilman S, Halub M, et al. A preventative respiratory protocol to identify trauma subjects at risk for respiratory compromise on a general in-patient ward. *Respir Care*. 2016. doi:10.4187/respcare.04729.

Baxter.com

Baxter International Inc. One Baxter Parkway / Deerfield, Illinois 60015

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Hillrom, Smart-Filter and Volara are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC189-220076 v3 08/23