



Staff Signature: \_\_\_\_\_

Trainer's Signature: \_\_\_\_\_

Date: \_\_\_\_\_



**BED OPERATIONS**

COMPLETED

1. Identify appropriate patient population (400 lbs – 1,000 lbs).	
2. Demonstrate how to lower and raise both head and foot siderails.	
3. Demonstrate the width and height adjustment of siderails.	
4. Demonstrate how to activate emergency CPR position.	
5. Demonstrate independent head and knee positioning and raising/lowering bed.	
6. Demonstrate Trendelenburg and Reverse Trendelenburg.	
7. Demonstrate cardiac chair position.	
8. Demonstrate brake and steer mechanism.	
9. Demonstrate patient transport procedures.	
10. Demonstrate how to zero the scale.	
11. Demonstrate how to weigh a patient.	
12. Demonstrate how to set up the weight change mode.	
13. Demonstrate how to switch between lbs and kgs.	
14. Demonstrate how to expand the sleep deck.	
15. Demonstrate how to expand surface with accompanying interchangeable mattress bolsters.	
16. Demonstrate egress from side of the bed.	
17. Identify location of angle indicators.	
18. Demonstrate how to activate battery backup.	



**OPTIONAL SYNERGY® AIR ELITE SURFACE**

COMPLETED

19. Demonstrate how to adjust comfort control feature.	
20. Initiate Alternating Therapy mode and demonstrate how to adjust cycle timer.	
21. Initiate Static mode.	
22. Demonstrate Auto Firm mode.	
23. Demonstrate how to activate CPR release on control unit.	

For more information, please contact your local distributor or Hillrom sales representative at 1-800-445-3730.

**hillrom.com**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.