Document Description:	CSB-CVSM/CIWS DISPLAY REPLACEMENT	Document Number: Version:	80020168 B		
	Welch Allyn™	Page 1	of 3		
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only					

Customer Service Bulletin

Product:	CVSM – Connex Vital Signs Monitor all models	Date:	2017-10-18
ı ı Ouuct.	O V OIVI OOI II CX VII AI OIGI IS IVIOI II IO I AII II IO I CIS	Daic.	_0.7.10

CIWS - Connex Integrated Wall Systems all models

CSB-CVSM/CIWS DISPLAY REPLACEMENT Subject:

CVSM Models

6300/6400/6500/6700/6800

HW Version(s) CIWS Models Affected:

8400/8500

SW Version(s) Affected:

All 1.x versions 1.71.03 and earlier. All 2.x versions 2.20.00 and earlier

Serial Numbers

Affected:

Lot or Date NA

Code Affected:

NA

Classification: Informational Only

Distribution: □ Customer Care □ Product Service

> ☐ Company Confidential ⊠ ASPs □ Distributors

Training Required: ☐ Yes \bowtie No

Summary: All CVSM and CIWS models are now being built with a new display (P/N 411161)

offering a more durable touchscreen paired with a longer life backlight LCD display. The old display (P/N 405883) will continue to be used in service kit P/N 103351 until inventory is exhausted. The new display will be available in service kit P/N 106546.

Issue: When the new display is installed in a CVSM or CIWS running older software versions an

artifact may be seen for about 300ms during power up. The artifact seen briefly on the display does not impact device performance or functionality. An updated version of the

host software prevents the artifact from being displayed.

Use software version 2.20.03 or later to update CVSM and CIWS models running version 2.xx

Use software version 1.71.04 to update CVSM and CIWS models running 1.xx

CVSM Models 6700/6800 and CIWS models 8400/8500 are currently shipping with version 2.30.00

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.

Document Description: CSB-CVSM/CIWS DISPLAY REPLACEMENT	Document Number: 80020168 Version: B			
Welch Allyn ™	Page 2 of 3			
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.				

CVSM Models 6300/6400/6500 ship with version 1.71.03.

Customers outside of China and Singapore may obtain the updated software via PartnerConnect by contacting Technical Support. The software update is not currently registered for use in Singapore or China.

To receive and install software upgrades for CVSM and CIWS you must have the Welch Allyn service tool installed on a Windows PC with an active internet connection.

Some customers may notice a reduced viewing angle. Increasing the display brightness and or positioning the device to minimize reflected light may improve readability.

The original display (P/N 405883) and the new display (P/N 411161) use a different bezel. The display service kits ship with a matching bezel. When replacing the display use the bezel included in the service kit.

Compatibility

CVSM/CIWS Compatibility with Support Software and Systems

	Welch Allyn								Custom Da	nta Schema
Device SW	Service Tool (WAST)	Service Monitor ¹	Service Dash board	Connex VM ^{2,3}	Connex CS ^{4,5}	Data Catcher	Connect ivity SDK ⁴	CDIS- NCE	Custom modifiers	Custom Scores
2.30.00	1.8	1.1, 1.2, 1.3, 1.3.1.0	2.0	2.0.1	1.5, 1.6, 1.7, 1.7.1	4.0.0.2 4.0.1.0	2.4 3.00.19	5.0.0.3/ 5.1.0.16	2.0/2.01.00	2.0/2.01.00
2.20.03	1.8 /1.7	1.1, 1.2, 1.3	2.0	2.0.1	1.5, 1.6, 1.7, 1.7.1	4.0.0.2	2.3/2.4/ 3.0.0.19	5.0.0.3	2.0	2.0
1.71.04	1.8 /1.7/1.6	N/A	N/A	2.0/2.0.1	1.1, 1.2, 1.5/1.6/ 1.7/ 1.7.1	4.0.0.2	2.3		NA	NA

Compatibility table foot notes:

- 1. The following CVSM 2.20.0x features are not supported in Service Monitor 1.1: Remote Configuration, manual 'Sync to server'
- 2. Does not support Custom Modifiers or Custom Scores from device.
- 3. Connex VM 1.1 does not support CVSM match to save features in host SW 1.71.02, 1.71.03 or 1.71.04.
- 4. Connex CS supports device SW 1.71.03 and 1.71.04 episodic monitoring only, continuous monitoring is not supported
- 5. See DIR 80018292 for detailed information on host software compatibility with Connex CS.

Related documents: 20014010 TSB-CVSM/CIWS Display replacement

20012631 CVSM/CIWS Minimum Hardware/Software level Requirement

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.

Document Description: CSB-CVSM/CIWS DISPLAY REPLACEMENT	Document Number: 80020168 Version: B			
Welch Allyn™	Page 3 of 3			
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only				

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	D*	GS
В	All	Updated to include CIWS. Replaced version 2.20.01 with 2.20.03	D*	GS
	Issue	1st bullet changed 2.20.01 to 2.20.03 or later for all 2.xx CVSM/CIWS 2nd bullet add CIWS for all 1.xx versions Delete notes 1 & 2 Change "CVSM models 6700 and 6800" to add CIWS models 8400/8500 Change 2nd paragraph from the end replace "provide some relief", with "improve read ability"		
	Compatibilit y table	Added row for SW version 2.30.00, added column for CDIS-NCE, Changed SW version 2.20.01 to 2.20.03 and updated compatibility		

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.