Hillrom™

Bariatric Nissen Strap/Thigh Retractors

Instructions for Use

Product No. A-80400, A-80400E, A-80400UK

80028174

Version D

INSTRUCTIONS FOR USE

ENGLISH	3
中文简体	19
中文繁體	36
ČESKÝ	52
DANSK	67
FRANÇAIS	82
DEUTSCH	99
ΕΛΛΗΝΙΚΑ	115
BAHASA INDONESIA	131
ITALIANO	147
日本語版	163
한국어	178
LATVISKI	193
PORTUGUÊS	208
ROMÂNESC	225
SRPSKI	242
SLOVENSKY	258
SLOVENŠČINA	274
ESPAÑOL.....	290
SVENSKA	306
TIẾNG VIỆT	320

INSTRUCTIONS FOR USE

IMPORTANT NOTICES

Prior to using this or any other type of medical apparatus with a patient, it is recommended that you read the **Instructions for Use** and familiarize yourself with the product.

- Read and understand all warnings in this manual and on the device itself prior to use with a patient.
- The  symbol is intended to alert the user to important procedures or safety instructions regarding the use of this device.

The  symbol on the labels is intended to show when the IFU should be referenced for use.

- The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this device remains with the attending physician.
- Device function should be checked prior to each usage.
- This device should only be operated by trained personnel.
- All modifications, upgrades, or repairs must be performed by an authorized specialist.
- Keep this manual available for future reference.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority listed in this document.

INSTRUCTIONS FOR USE

Table of Contents

Bariatric Nissen Strap/Thigh Retractors (A-80400, A-80400E, A-80400UK)

1. General Information:	6
1.1 Copyright Notice:	6
1.2 Trademarks:	6
1.3 Contact Details:	7
1.4 Safety Considerations:	7
1.4.1 Safety hazard symbol notice:	7
1.4.2 Equipment misuse notice:	7
1.4.3 Notice to users and/or patients:	7
1.4.4 Safe disposal:	8
1.5 Operating the system:	8
1.5.1 Applicable Symbols:	8
1.5.2 Intended User and Patient Population:	9
1.5.3 Compliance with medical device regulations:	10
1.6 EMC considerations:	10
1.7 EC authorized representative:	10
1.8 Manufacturing Information:	10
1.9 EU Importer Information:	10
1.10 Authorised Australian sponsor:	10
2 System	11
2.1 System components Identification:	11
2.2 Product Code and Description:	11
2.3 List of Accessories and Consumable Components Table:	11
2.4 Indication for use:	12
2.5 Intended use:	12
3 Equipment Setup and Use:	12
3.1 Prior to use:	12
3.2 Setup:	12
3.3 Device controls and indicators:	13
3.4 Storage, Handling and Removal Instructions:	13
3.4.1 Storage and Handling:	13

INSTRUCTIONS FOR USE

3.4.2 Removal Instruction:13

3.5 Troubleshooting Guide:14

3.6 Device Maintenance:14

4 Safety Precautions and General Information:14

4.1 General Safety Warnings and Cautions:14

4.2 Product Specifications:15

4.3 Sterilization Instruction:.....15

4.4 Cleaning and Disinfection Instruction:15

5 List of Applicable Standards:16

INSTRUCTIONS FOR USE

1. General Information:

Allen Medical Systems, Inc. is a subsidiary of Hill-Rom, Inc., (NYSE: HRC), a leading worldwide manufacturer and provider of medical technologies and related services for the health care industry. As an industry leader in patient positioning, our passion is improving patient outcomes and caregiver safety, while enhancing our customers' efficiency. Our inspiration comes from providing innovative solutions to address our customers' most pressing needs. We immerse ourselves in our customers' world, to better address these needs and the daily challenges of their environment. Whether developing a solution to address patient positioning challenges or creating a system to offer safe and effective surgical site access for the surgical team, we are committed to providing products of exceptional value and quality.

Allen products are backed by responsive and reliable service and complimentary on-site product demonstrations.

1.1 Copyright Notice:

Revision

© 2019 Allen Medical Systems Inc. ALL RIGHTS RESERVED.

No part of this text shall be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or by any information or retrieval system without written permission from Allen Medical Systems, Inc. (Allen Medical).

The information in this manual is confidential and may not be disclosed to third parties without the prior written consent of Allen Medical.

1.2 Trademarks:

Trademark information can be found at [Allenmedical.com/pages/terms-conditions](https://allenmedical.com/pages/terms-conditions).

Products may be covered by one or more patents. Please consult listing at [Hill-rom.com/patents](https://hill-rom.com/patents) for any patent(s).

INSTRUCTIONS FOR USE

1.3 Contact Details:

For ordering information please see catalog.

Allen Customer Service Contact Information:

 North America	 International
 (800) 433-5774	 +1 978 266 4200 ext.4286
 (978) 263-8846	 +1 978 266 4426

1.4 Safety Considerations:

1.4.1 Safety hazard symbol notice:

DO NOT USE IF PRODUCT SHOWS VISIBLE DAMAGE OR MATERIAL DEGRADATION.

1.4.2 Equipment misuse notice:

Do not use the product if package is damaged or unintentionally opened before use.

All modifications, upgrades, or repairs must be performed by an authorized specialist.

1.4.3 Notice to users and/or patients:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Note: Refer to the surgical table manufacturer's user guide for instructions on use. Always refer to the surgical table manufacturer's weight limits.

NEVER EXCEED THE WEIGHT CAPACITY OF THE OPERATING ROOM TABLE

INSTRUCTIONS FOR USE

1.4.4 Safe disposal:

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories. If in doubt, the user of the device shall first contact Hill-Rom Technical Support for guidance on safe disposal protocols.

1.5 Operating the system:

1.5.1 Applicable Symbols:

Symbol used	Description	Reference
	Indicates the device is a medical device	MDR 2017/745
	Indicates the medical device manufacturer	EN ISO 15223-1
	Indicates the manufacturer's serial number. The device serial number is encoded as 1YYWWSSSSSS . YY indicates the year of manufacture. i.e. 18WWSSSSSS where 18 represents the year 2018. WW indicates the number of the manufacturing week per a standard shop calendar. (Leading zeros included.) SSSSSS is a sequential unique number.	EN ISO 15223-1
	Indicates the medical device Global Trade Item Number	21 CFR 830 MDR 2017/745
	Indicates the manufacturer's lot code using the Julian Date yyddd, where yy indicates the last two digits of the year and ddd indicates the day of the year. i.e. April 4 th , 2019 would be represented as 19094.	EN ISO 15223-1

INSTRUCTIONS FOR USE

	Indicates the date when the medical device was manufactured	EN ISO 15223-1
	Indicates the manufacturer's catalogue number	EN ISO 15223-1
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions.	EN ISO 15223-1
	Indicates the device do not contain natural rubber or dry natural rubber latex	EN ISO 15223-1
	Indicates the authorized representative in the European Community	EN ISO 15223-1
	Indicates the Medical Device complies to REGULATION (EU) 2017/745	MDR 2017/745
	Indicates a Warning	IEC 60601-1
	Indicates the need for the user to consult the instruction for use	EN ISO 15223-1

1.5.2 Intended User and Patient Population:

Intended User: Surgeons, Nurses, Doctors, Physicians and OR healthcare professionals involved in the device intended procedure. Not intended for Lay persons.

Intended Populations:

This device is intended to be used with patients that do not exceed the weight in the safe working load field specified in the product specification section 4.2

INSTRUCTIONS FOR USE

1.5.3 Compliance with medical device regulations:

This Product is a non-invasive, Class I Medical Device. This system is CE-marked according to Annex VIII, Rule 1, of the Medical Device Regulations (REGULATION (EU) 2017/745)

1.6 EMC considerations:

This is not an electromechanical device. Therefore, EMC Declarations are not applicable.

1.7 EC authorized representative:

HILL-ROM SAS
B.P. 14 - Z.I. DU TALHOUE
56330 PLUVIGNER
FRANCE
TEL: +33 (0)2 97 50 92 12

1.8 Manufacturing Information:

ALLEN MEDICAL SYSTEMS, INC.
100 DISCOVERY WAY
ACTON, MA 01720 USA
800-433-5774 (NORTH AMERICA)
978-266-4200 (INTERNATIONAL)

1.9 EU Importer Information:

Baxter Medical Systems GmbH + Co. KG
Carl-Zeiss Straße 7-9
07318 Saalfeld/Saale
Germany


1.10 Authorised Australian sponsor:

Welch Allyn Australia Pty. Ltd.
1 Baxter Drive
Old Toongabbie, NSW 2146
Australia

INSTRUCTIONS FOR USE

2 System

2.1 System components Identification:

2.2 Product Code and Description:

A-80400 - Bariat Nissen Straps / Thigh Retractor, US

A-80400E - Bariat Nissen Straps/Thigh Retractor, EU

A-80400UK - Bariatric Nissen Straps/Thigh Retractor, UK

2.3 List of Accessories and Consumable Components Table:

The following list are accessories and components that may be used with this device.

Name of Accessory	Product Number
Accessory Cart	A-30010 (US), A-30013E (EU), A-30010UK (UK), A-30010J (Japan)

Name of Consumable	Product Number
Not Applicable	Not Applicable

Note: Consult the corresponding IFU for the products mentioned in the above table.

INSTRUCTIONS FOR USE

2.4 Indication for use:

The Nissen strap is used in a variety of surgical procedures including, but not limited to gynecology, urology, laparoscopy, general and robotic surgery during Reverse Trendelenburg positioning. These devices are capable of being used with a broad patient population as determined appropriate by the caregiver or institution.

2.5 Intended use:

The Nissen strap is designed to position and support the patient's upper leg in a variety of surgical procedures including, but not limited to gynecology, urology, laparoscopy, general and robotic surgery during Reverse Trendelenburg positioning. These devices are intended to be used by healthcare professionals within the Operating Room setting.

3 Equipment Setup and Use:

3.1 Prior to use:

- a. Inspect the product looking for any visible damage or sharp edges that could be caused by a drop or impact during storage.
- b. Make sure product has been properly cleaned and disinfected and wiped dry prior to each use.


3.2 Setup:

- a. Position the Strap Mount on the table rail adjacent to the patient's thigh. Tighten.
- b. Pull the strap under and around the patient's thigh. The pad must be between the thigh and the strap encircling the thigh.

INSTRUCTIONS FOR USE

- c. Position the buckle mount 150mm (6") from the strap mount on the table rail toward the head end of the table. Tighten.

Strap Section & Mount Buckle Buckle Section & Mount

- d. Pull the strap through the buckle until it is taut, and the patient is held on the table. Lock by pressing the buckle latch to the closed position.
- e. Repeat steps a. through d. for the other strap.

3.3 Device controls and indicators:

Controls and indicators of this device are described within the Setup instructions.

3.4 Storage, Handling and Removal Instructions:

3.4.1 Storage and Handling:

The product should be stored in a clean and safe environment to prevent product damage. See storage Specifications under Product Specification section.

3.4.2 Removal Instruction:

- a. Loosen the strap from the buckle.
- b. Loosen the buckle section mount knob and remove mount.
- c. Release the patient thigh and remove the strap section mount.

INSTRUCTIONS FOR USE

3.5 Troubleshooting Guide:

This device does not have a troubleshooting guide. For technical support user of the device shall first contact Hill-Rom Technical Support.

3.6 Device Maintenance:

Periodic Maintenance is not required for this device.

4 Safety Precautions and General Information:

4.1 General Safety Warnings and Cautions:

WARNING:

- a. Do not use if product shows visible damage.
- b. Prior to using this device, please read the instructions for equipment set up and use. Familiarize yourself with the product before application on a patient.
- c. To prevent patient and/or user injury and/or equipment damage, examine the device and surgical-table side rails for potential damage or wear prior to use. Do not use the device if damage is visible, if parts are missing or if it does not function as expected.
- d. Do not attach strap to any type of table width extender as patient injury could result

CAUTION:

- a. Do not over-tighten straps around patient's legs as injury may result.
- b. Do not exceed safe working load shown in the product specification table.

INSTRUCTIONS FOR USE

4.2 Product Specifications:

Mechanical Specifications	Description
Product Dimensions	Pad Dimension 79 cm x 18 cm (31" x 7") (L x W)
Material	Aluminum, Conductive strap, Viscoelastic Foam, Fusion 3, Black, Lectrolite Fusion 3HP, Black, Lectrolite Duotone, Black, Vinyl-nitrile rubber
Safe Working Load on the device	1.000 lbs. (453.5 kg)
Overall Weight of Complete Device	3.6 lbs. (1.63 kg)
Storage Specifications	Description
Storage temperature	-29° C to +60° C
Storage Relative humidity range	15% to 85%
Operating temperature	This device is intended to be used in a controlled Operating Room environment.
Operating Relative humidity range	
Electrical Specifications	Description
Not Applicable.	Not Applicable.
Software Specifications	Description
Not Applicable.	Not Applicable.
Compatibility Specifications	Description
Operating Room Table Compatibility	The Bariatric Nissen Thigh Strap is compatible with the following surgical table rail styles: US, UK, EU

Note: Consult the corresponding IFU for the products mentioned in the above table.

4.3 Sterilization Instruction:

This device is not intended to be sterilized.

4.4 Cleaning and Disinfection Instruction:

WARNING:

INSTRUCTIONS FOR USE

- Do not use bleach or products that contain bleach to clean the device. Injury or equipment damage can occur.
- After each use, clean the device with alcohol-based wipes.
- Do not put the device into water. Equipment damage can occur.
- Use a cloth and a quaternary ammonium disinfecting/cleaning solution to clean and disinfect the device.
- Read and follow the manufacturer's recommendation for low-level disinfection.
- Read and follow the cleaning product's instructions. Use caution in areas where liquid can get into the mechanism.
- Wipe the device with a clean, dry cloth.
- Make sure that the device is dry before you store it or use it again.

CAUTION: DO NOT IMMERSE PADS IN ANY LIQUID

CAUTION: DO NOT USE BLEACH OR PHENOLICS ON PADS

5 List of Applicable Standards:

Sl. no	Standards	Description
1	EN 62366-1	Medical devices - Part 1: Application of usability engineering to medical devices
2	EN ISO 14971	Medical devices- Application of risk management to medical devices.
3	EN 1041	Information supplied by the manufacturer of medical devices
4	EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
5	EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
6	IEC 60601-2-46	Medical electrical equipment - Part 2-46: Particular requirements for the basic safety and essential performance of operating tables

INSTRUCTIONS FOR USE

7	ISTA	International Safe Transit Association standards for package testing
---	------	--

Hillrom™

Bariatric Nissen Strap/Thigh Retractors

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