Document Description:	CSB-ELI Security-Hard Coded Password	Document Number:	80029673
		Version:	A
	Welch Allyn		
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.			

Customer Service Bulletin					
Products:	Welch	_	BUR150c, Welch Allyn ELI 250c, and Welch Allyn ELI 380 Resting	Date:	16-Jun-2022
Subject:	Security Vulnerability of Hard Coded Password				
HW Version(s) Affected:	N/A	SW Version(s) Affected:	Welch Allyn ELI 150c device – v2.2.0 and prior BUR150c device – v2.1.2 and prior Welch Allyn ELI 250c device – v2.1.2 and prior Welch Allyn ELI 280 device - v2.3.1 and prior Welch Allyn ELI 380 device - v2.6.0 and prior		
Serial Numbers Affected:	N/A	Lot or Date Code Affected:	N/A		

Disposition:	Reactive Service Use			
Distribution:	□ Customer Care		⊠ Field Service	
	⋈ ASPs ⋈ Distributors	□ Customers	☐ Company Confidential	
Training Required:	⊠ Yes □ No			
Summary:	The specified products utilize hard-coded password information to control inbound authentication and/or outbound communication to external components. This software architecture has been identified as presenting a security vulnerability related to Improper Access Control at the device. No known public exploits specifically target these vulnerabilities. This vulnerability has a low attack complexity.			

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.

Mitigation:

Hillrom is releasing software updates for all impacted devices to address these vulnerabilities. New product versions that mitigate these vulnerabilities are available as follows:

- Welch Allyn ELI 380 Resting Electrocardiograph: available by Q4 2023
- **Welch Allyn ELI** 280/BUR280/MLBUR 280 Resting Electrocardiograph: Software version 2.4.1, available June 2022
- Welch Allyn ELI 150c/BUR 150c/MLBUR 150c Resting Electrocardiograph: available by Q4 2023

Hillrom recommends users upgrade to the latest product versions. Customers should contact Technical Support for assistance.

The following actions can be taken by facility IT resources to reduce risk associated with this vulnerability.

- Apply proper network and physical security controls.
- Ensure a unique encryption key is configured for ELI Link and cardiographs, per instructions in the product user manual.
- Locate control system networks and remote devices behind firewalls and isolate them from the business network.
- When remote access is required, use secure methods, such as Virtual Private Networks (VPNs), recognizing that VPNs may have vulnerabilities and should be updated to the most current version available. Also recognize that a VPN is only as secure as the connected devices.

Organizations observing any suspected malicious activity should follow their established internal procedures and report their finding to CISA or tracking and correlation against other incidents.

CISA Contact Information:

- Phone +1-888-282-0870
- Email <u>CISAservicedesk@cisa.dhs.gov</u>

Welch Allyn and ELI are registered trademarks of Baxter International Inc. or its subsidiaries.

Resolution:

This issue has been documented into our feedback system (B104663, B104649, B104661) for correction in a future software release.

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	16-Jun-2022	SLB

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.

THIS INFORMATION IS THE PROPERTY OF WELCH ALL THE MANUFACTURE OR SALE OF EQUIPMENT OF	YN, INC. AND AS SUCH S R DEVICES WITHOUT TH	SHALL NOT BE REPR IE EXPRESS WRITTEI	ODUCED, COPIE N PERMISSION C	ED, OR USED AS A BASIS OF WELCH ALLYN, INC.	
FCD-90159538 Rev A	Parent: WI-90159537	APR372702 Rev	A 09-Jun-2022		3