

Preparation and Application

PREPARATION



1 Remove from box and allow **Preveleak** Surgical Sealant to reach room temperature prior to use.



2 Open the outer pouch and pass inner sterile pouch into sterile field.



3 Remove the double-barrel syringe and the delivery tips from the inner pouch.

Hold the syringe with the cap end upward. Remove the syringe cap by turning it 90° counterclockwise and pulling upwards, using a gentle side-to-side motion.



4 To attach a delivery tip, orient the small alignment tab on the hub end of the delivery tip with the corresponding V-shaped alignment notch on the syringe.



5 Firmly push delivery tip toward the syringe and turn it 90° clockwise to lock the tip in place.

APPLICATION



1 Discard the first 0.25 mL of sealant immediately prior to use.



2A | APPLIED TO SUTURE LINE
Apply **Preveleak** Surgical Sealant to a clamped application site in a slow and steady manner, approximating the suture line with the delivery tip.



2B | APPLIED TO SYNTHETIC GRAFT



3 Leave the clamps in place for at least 60 seconds before restoring circulation or touching sealant.