Hillrom™

Lift Assist Beach Chair

Instructions for Use Product

No. A-91501

80028178
Version E

INSTRUCTIONS FOR USE

ENGLISH.....	3
中文简体.....	34
中文繁體.....	66
DANSK	99
NEDERLANDS.....	130
SUOMI.....	162
FRANÇAIS	193
DEUTSCH.....	225
ΕΛΛΗΝΙΚΑ.....	257
BAHASA INDONESIA	288
ITALIANO	322
日本語版.....	353
한국어	386
LATVISKI	417
LIETUVIŲ	448
NORSK	480
POLSKI	512
PORTUGUÊS.....	544
ROMÂNESC.....	575
SRPSKI	607
SLOVENSKY.....	640
SLOVENŠČINA.....	671
ESPAÑOL	702
SVENSKA.....	734
TIẾNG VIỆT	766

INSTRUCTIONS FOR USE

IMPORTANT NOTICES

Prior to using this or any other type of medical apparatus with a patient, it is recommended that you read the Instructions for Use and familiarize yourself with the product.

- Read and understand all warnings in this manual and on the device itself prior to use with a patient.
- The  symbol is intended to alert the user to important procedures or safety instructions regarding the use of this device.

The  symbol on the labels is intended to show when the IFU should be referenced for use.

- The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this device remains with the attending physician.
- Device function should be checked prior to each usage.
- This device should only be operated by trained personnel.
- All modifications, upgrades, or repairs must be performed by an authorized specialist.
- Keep this manual available for future reference.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority listed in this document.

INSTRUCTIONS FOR USE

Table of Contents

Allen® Lift Assist Beach Chair (A-91501)

1	General Information:	6
1.1	Copyright Notice:	6
1.2	Trademarks:	6
1.3	Contact Details:	7
1.4	Safety Considerations:	7
1.4.1	Safety hazard symbol notice:	7
1.4.2	Equipment misuse notice:	7
1.4.3	Notice to users and/or patients:	7
1.4.4	Safe disposal:	8
1.5	Operating the system:	8
1.5.1	Applicable Symbols:	8
1.5.2	Intended User and Patient Population:	9
1.5.3	Compliance with medical device regulations:	10
1.6	EMC considerations:	10
1.7	EC authorized representative:	10
1.8	Manufacturing Information:	10
1.9	EU Importer Information:	10
1.10	Authorised Australian sponsor:	10
2	System	11
2.1	System components Identification:	11
2.2	Product Code and Description:	13
2.3	List of Accessories and Consumable Components Table:	13
2.4	Indication for use:	14
2.5	Intended use:	14
2.6	Residual Risk:	14
3	Equipment Setup and Use:	15
3.1	Prior to use:	15
3.2	Setup:	15
3.3	Device controls and indicators:	27

INSTRUCTIONS FOR USE

3.4	Storage, Handling and Removal Instructions:.....	28
3.4.1	Storage and Handling:	28
3.4.2	Removal Instruction:.....	28
3.5	Troubleshooting Guide:	29
3.6	Device Maintenance:.....	29
4	Safety Precautions and General Information:	29
4.1	General Safety Warnings and Cautions:.....	29
4.2	Product Specifications:.....	30
4.3	Sterilization Instruction:	31
4.4	Cleaning and Disinfection Instruction:.....	31
5	List of Applicable Standards:	32

INSTRUCTIONS FOR USE

1 General Information:

Allen Medical Systems, Inc. is a subsidiary of Hill-Rom, Inc., (NYSE: HRC), a leading worldwide manufacturer and provider of medical technologies and related services for the health care industry. As an industry leader in patient positioning, our passion is improving patient outcomes and caregiver safety, while enhancing our customers' efficiency. Our inspiration comes from providing innovative solutions to address our customers' most pressing needs. We immerse ourselves in our customers' world, to better address these needs and the daily challenges of their environment. Whether developing a solution to address patient positioning challenges or creating a system to offer safe and effective surgical site access for the surgical team, we are committed to providing products of exceptional value and quality.

Allen products are backed by responsive and reliable service and complimentary on-site product demonstrations.

1.1 Copyright Notice:

Revision

© 2019 Allen Medical Systems Inc. ALL RIGHTS RESERVED.

No part of this text shall be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or by any information or retrieval system without written permission from Allen Medical Systems, Inc. (Allen Medical).

The information in this manual is confidential and may not be disclosed to third parties without the prior written consent of Allen Medical.

1.2 Trademarks:

Trademark information can be found at [Allenmedical.com/pages/terms-conditions](https://allenmedical.com/pages/terms-conditions).

Products may be covered by one or more patents. Please consult listing at [Hill-rom.com/patents](https://hill-rom.com/patents) for any patent(s).

INSTRUCTIONS FOR USE

1.3 Contact Details:

For ordering information please see catalog.

Allen Customer Service Contact Information:

 North America	 International
 (800) 433-5774	 +1 978 266 4200 ext.4286
 (978) 263-8846	 +1 978 266 4426

1.4 Safety Considerations:

1.4.1 Safety hazard symbol notice:

	DO NOT USE IF PRODUCT SHOWS VISIBLE DAMAGE OR MATERIAL DEGRADATION.

1.4.2 Equipment misuse notice:

Do not use the product if package is damaged or unintentionally opened before use. All modifications, upgrades, or repairs must be performed by an authorized specialist.

1.4.3 Notice to users and/or patients:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Note: Refer to the surgical table manufacturer's user guide for instructions on use. Always refer to the surgical table manufacturer's weight limits.

	NEVER EXCEED THE WEIGHT CAPACITY OF THE OPERATING ROOM TABLE

INSTRUCTIONS FOR USE

1.4.4 Safe disposal:


Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories. If in doubt, the user of the device shall first contact Hill-Rom Technical Support for guidance on safe disposal protocols.

1.5 Operating the system:

1.5.1 Applicable Symbols:

Symbol used	Description	Reference
	Indicates the device is a medical device	MDR 2017/745
	Indicates the medical device manufacturer	EN ISO 15223-1
	Indicates the manufacturer's serial number. The device serial number is encoded as 1YYWWSSSSSS . YY indicates the year of manufacture. i.e. 18WWSSSSSS where 18 represents the year 2018. WW indicates the number of the manufacturing week per a standard shop calendar. (Leading zeros included.) SSSSSS is a sequential unique number.	EN ISO 15223-1
	Indicates the medical device Global Trade Item Number	21 CFR 830 MDR 2017/745
	Indicates the manufacturer's lot code using the Julian Date yyddd, where yy indicates the last two digits of the year and ddd indicates the day of the year. i.e. April 4 th , 2019 would be represented as 19094.	EN ISO 15223-1

INSTRUCTIONS FOR USE

Symbol used	Description	Reference
	Indicates the date when the medical device was manufactured	EN ISO 15223-1
	Indicates the manufacturer's catalogue number	EN ISO 15223-1
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions.	EN ISO 15223-1
	Indicates the device do not contain natural rubber or dry natural rubber latex	EN ISO 15223-1
	Indicates the authorized representative in the European Community	EN ISO 15223-1
	Indicates the Medical Device complies to REGULATION (EU) 2017/745	MDR 2017/745
	Indicates a Warning	IEC 60601-1
	Indicates the need for the user to consult the instruction for use	EN ISO 15223-1

1.5.2 Intended User and Patient Population:

Intended User: Surgeons, Nurses, Doctors, Physicians and OR healthcare professionals involved in the device intended procedure. Not intended for Lay persons.

Intended Populations:

This device is intended to be used with patients that do not exceed the weight in the safe working load field specified in the product specification section 4.2.

INSTRUCTIONS FOR USE

1.5.3 Compliance with medical device regulations:

This Product is a non-invasive, Class I Medical Device. This system is CE-marked according to Annex VIII, Rule 1, of the Medical Device Regulations (REGULATION (EU) 2017/745)

1.6 EMC considerations:

This is not an electromechanical device. Therefore, EMC Declarations are not applicable.

1.7 EC authorized representative:

HILL-ROM SAS
B.P. 14 - Z.I. DU TALHOUE
56330 PLUVIGNER
FRANCE
TEL: +33 (0)2 97 50 92 12

1.8 Manufacturing Information:

ALLEN MEDICAL SYSTEMS, INC.
100 DISCOVERY WAY
ACTON, MA 01720 USA
800-433-5774 (NORTH AMERICA)
978-266-4200 (INTERNATIONAL)

1.9 EU Importer Information:

Baxter Medical Systems GmbH + Co. KG
Carl-Zeiss Straße 7-9
07318 Saalfeld/Saale
Germany

1.10 Authorised Australian sponsor:


Welch Allyn Australia Pty. Ltd.
1 Baxter Drive
Old Toongabbie, NSW 2146
Australia

INSTRUCTIONS FOR USE

2 System



2.1 System components Identification:

Allen® Lift-Assist™ Beach Chair (Front)

INSTRUCTIONS FOR USE

Allen® Lift-Assist™ Beach Chair (back)

INSTRUCTIONS FOR USE

2.2 Product Code and Description:

A-91501 - Allen Lift Assist Beach Chair

A-91502 - Allen Intubation Plate and Pad

A-91503 - Allen Lift Assist Beach Chair Clamp

2.3 List of Accessories and Consumable Components Table:

The following list are accessories and components that may be used with this device.

Name of Accessory	Product Number
Lift-Assist Beach Chair Reusable Pad Set	A-90040
Intubation Plate and Pad	A-91502
Lift-Assist Beach Chair Clamp	A-91503
Easy lock Socket	A-40022
LPS Arm Support	A-90002
Easy lock Socket	A-40022
Counter Traction Strap	A-90015

Name of Consumable	Product Number
Universal head Restraint Disposables	A-90023

Note: Consult the corresponding IFU for the products mentioned in the above table.

CAUTION: Reuse of disposables will lead to patient cross contamination and/or device failure.

INSTRUCTIONS FOR USE

2.4 Indication for use:

The Beach Chair Shoulder Positioner is used in a variety of surgical procedures including, but not limited to orthopedic surgeries like, rotator cuff, total shoulder, reverse total shoulder, SLAP repair and other shoulder surgery. These devices are capable of being used with a broad patient population as determined appropriate by the caregiver or institution.

2.5 Intended use:

The Beach Chair Shoulder Positioner is designed to position and support the patient's head and torso in a variety of surgical procedures including, but not limited to orthopedic surgeries like, rotator cuff, total shoulder, reverse total shoulder, SLAP repair and other shoulder surgery. These devices are intended to be used by healthcare professionals within the Operating Room setting.

2.6 Residual Risk:

This product complies with relevant performance, safety standards. However, device or patient harm from misuse or device damage hazards cannot be completely excluded.

INSTRUCTIONS FOR USE

3 Equipment Setup and Use:

3.1 Prior to use:

- a. Inspect the product looking for any visible damage or sharp edges that could be caused by a drop or impact during storage.
- b. Make sure product has been properly cleaned and disinfected and wiped dry prior to each use.
- c. Ensure surgical table is in its final operation location with respect to the room.
- d. Ensure the table is locked.
- e. Remove X-ray tops and back section pads from the surgical table.
- f. Remove or lower the back (Preferred) or leg section of the table by 90°.
- g. Loosen both rail clamp mount knobs.

3.2 Setup:

1. Unlock shoulder panel locks and remove shoulder panels as shown below by swinging panels back and unhooking from head positioner tubes. This will make beach chair lighter and easier to manage.

INSTRUCTIONS FOR USE


2. Place top jaw of the rail clamp onto the torso (seat) section of the table rail as shown.

3. Snap hinged bottom jaw onto the rail as shown.


4. Slide rail clamp to the end of the table rail. Preferably to rail stops.


5. Repeat steps 2-4 with second rail clamp on opposite side of table.

6. Loosen width extension knob by turning to the left.

INSTRUCTIONS FOR USE


7. Fully extend width of the support bar.

8. Back-out (loosen) rail clamp knobs completely.

9. Lift Beach Chair by placing dominant hand on the red release handle and opposite hand on the appropriate side as shown. **The non-dominant hand should be placed well behind the rail clamp mount.**

INSTRUCTIONS FOR USE


10. Align rail clamp mounts above rail clamps.


11. Rest rail clamp mount on edge of rail clamp.

12. Adjust width of support bar as needed until second rail clamp mount is correctly aligned with second rail clamp.

13. Lower rail clamp mount onto rail clamp. Rail clamp mount should slide over rail clamp. Do NOT force. Rail clamp mount is not installed properly until green is visible through indicator slot and pin is flush with top (both circled).

INSTRUCTIONS FOR USE


14. Ensure there is No Gap between the Beach Chair and the OR table pad. Ensure blue foam roll is hanging below the Beach Chair.

15. Tighten both rail clamp mount knobs by turning clockwise.


16. Tighten width extension knob by turning it clockwise (to the right).

INSTRUCTIONS FOR USE

17. Replace shoulder panels by hooking the latches around the head positioner tubes and swinging the panel forward as shown. Engage shoulder panel locks by turning it outward.


18. Adjust lateral position of Beach Chair by gently sliding it near the base along support bar until desired position is achieved.

INSTRUCTIONS FOR USE

19. Turn Beach Chair release handle to the right or left. Lower Beach Chair slowly until it lies horizontal, flush with the OR Table pad. Use head positioner handle to assist in lowering.

20. Transfer patient onto table in supine position. Ensure patient's hip is located as shown in image (ALIGNED WITH EDGE OF O.R. TABLE).


21. Ensure patient's head is placed on the intubation pad as shown below. Perform necessary anesthetization and intubation.


INSTRUCTIONS FOR USE


22. Open head positioner tube release as shown below.

23. Extend head positioner tubes by grasping the handle at the top of the tubes and pulling gently outward. Be sure to keep intubation pad in place and monitor the patient's head during this process.


24. Slide patient head positioner onto head positioner tubes.

INSTRUCTIONS FOR USE

25. Lift patient's head and remove intubation pad. Slide head positioner forward and place the patient's head in the device.


Ensure that patient's head and neck are always fully supported.

26. Secure patient's head into the head positioner as dictated by the device user instructions.


27. Ensure head positioner tubes are extended at $\geq 8''$ (20.3cm) above head positioner.

INSTRUCTIONS FOR USE

28. Lock head positioner tubes by engaging head positioner tube release as shown below.

29. Place non-operative arm in restraint pad as shown.

30. There are 4 simple steps to achieving the modified Fowler Position (the “Beach Chair” position).

- a. Place OR Table into 5°-10° of Trendelenburg. **Monitor the patient's head and neck during adjustments.**

INSTRUCTIONS FOR USE

- b. While holding head positioner tube handle, unlock Beach Chair release handle (shown in step 20) and raise Beach Chair and patient approximately 45°. **Monitor the patient's head and neck during adjustments.**

- c. Flex OR Table to achieve flexion in the knee (10°-30° of flexion is desirable). **Monitor the patient's head and neck during adjustments.**

- d. Raise Beach Chair until desired angle is achieved. **Monitor the patient's head and neck during adjustments.**

INSTRUCTIONS FOR USE


31. Unlock and remove shoulder panels from operative side of patient as shown below.

Note: Be sure to support the patient and arm to avoid unwanted movement.

32. Apply counter traction device to non-operative side of patient using a strap or brace. Be sure to follow the device user instructions.

33. Use headrest lock if needed. Slide headrest lock below head positioner and tighten knob.

INSTRUCTIONS FOR USE

3.3 Device controls and indicators:

This label on the device indicates the rail clamp knob locks (clockwise) and unlock (counter clockwise).

This label on the device indicates "Do not lean on panels".

This label on the device indicates unlocking the side panels.

This label on the device indicates locking the side panels.

INSTRUCTIONS FOR USE

3.4 Storage, Handling and Removal Instructions:

3.4.1 Storage and Handling:

The product should be stored in a clean and safe environment to prevent product damage. See storage Specifications under Product Specification section 4.2.

3.4.2 Removal Instruction:

1. Remove counter traction device, arm supports or other restrain devices from patient.
2. Ensure head positioner lock is unlocked and free floating.

3. Reverse steps bb to dd to return patient to horizontal position.

Warning: Monitor the patient's head and neck during adjustments

4. Remove patient's head from the head positioner. While supporting the patient's head, remove head positioner from head positioner tubes and replace intubation pad before placing patient's head down.

5. Remove patient from device and remove Beach Chair from OR table by reversing steps 11-18.

INSTRUCTIONS FOR USE

3.5 Troubleshooting Guide:

This device does not have a troubleshooting guide. For technical support user of the device shall first contact Hill-Rom Technical Support.

3.6 Device Maintenance:

Make sure that all labels are installed and can be read. Replace labels as necessary by using a plastic scraper to remove the label. Use an alcohol wipe to remove any adhesive residue.

Contact Allen Medical Systems, Inc. if you need to repair or replace the device contact us using the information from the contact details section (1.3).

4 Safety Precautions and General Information:

4.1 General Safety Warnings and Cautions:

WARNING:

- a. Do not use if product shows visible damage.
- b. Prior to using this device, please read the instructions for equipment set up and use. Familiarize yourself with the product before application on a patient.
- c. To prevent patient and/or user injury and/or equipment damage, examine the device and surgical-table side rails for potential damage or wear prior to use. Do not use the device if damage is visible, if parts are missing or if it does not function as expected.
- d. Be sure to support the patient and arm to avoid unwanted movement.
- e. Monitor the patient's head and neck during adjustments.
- f. Whenever inclining the patient, monitor the head closely to ensure that the head positioner slides freely without obstruction to prevent undue pressure on the patient's head and neck.
- g. Ensure that patient's head and neck are always fully supported.
- h. Do not exceed 500 lbs. (226 kg) patient weight.

CAUTION:

- a. Do not exceed safe working load shown in the product specification table
 - b. Do not apply excessive force to shoulder panel (61 lbs (27 kg) weight limit).

INSTRUCTIONS FOR USE

4.2 Product Specifications:

Mechanical Specifications	Description
Product Dimensions	38" X 26" X 13 ^{3/4} " (96.5 cm x 66 cm x 35 cm)
Material	Stainless Steel, Aluminum, Fiberboard, Delrin, Ultra-High Molecular Weight Polyethylene, Nylon, Velcro®, Foam (High Density & Memory), Fusion 3 (or equivalent), Lectrolite (or equivalent), Glass-Filled Nylon, Polyurethane
Safe Working Load on the device	500 lbs. (226 kg)
Overall Weight of Complete Device	35 lbs. (15.9 kg)
Storage Specifications	Description
Storage temperature	-29° C to +60° C
Storage Relative humidity range	15% to 85%
Operating temperature	This device is intended to be used in a controlled Operating Room environment.
Operating Relative humidity range	
Electrical Specifications	Description
Not Applicable.	Not Applicable.
Software Specifications	Description
Not Applicable.	Not Applicable.
Compatibility Specifications	Description
The Lift Assist Beach Chair is compatible with:	Lift-Assist Beach Chair Clamps • A-91503 (US) • A-91503 (EU) • A-91503 (UK) • A-91503 (JP) Intubation Plate and Pad • A-91502
Operating Room Table Compatibility	The Lift Assist Beach Chair is compatible with the following surgical table rail styles: US, UK, EU, DEN, JP, SWISS. 3/8" x 1 1/8" (US Rail Size) 10 x 25 mm (EU Rail Size) 6.35 x 31.75 mm (UK Rail Size) 6.35 x 38.1 mm (Denyer Rail Size) 5. 9 x 32 mm (Japanese Rail Size)

Note: Consult the corresponding IFU for the products mentioned in the above table.

INSTRUCTIONS FOR USE

4.3 Sterilization Instruction:

This device is not intended to be sterilized. Equipment damage may occur.

4.4 Cleaning and Disinfection Instruction:

WARNING:

- Do not use bleach or products that contain bleach to clean the device. Injury or equipment damage can occur.
- After each use, clean the device with alcohol-based wipes.
- Do not put the device into water. Equipment damage can occur.
- Use a cloth and a quaternary ammonium disinfecting/cleaning solution to clean and disinfect the device.
- Read and follow the manufacturer's recommendation for low-level disinfection.
- Read and follow the cleaning product's instructions. Use caution in areas where liquid can get into the mechanism.
- Wipe the device with a clean, dry cloth.
- Make sure that the device is dry before you store it or use it again.

CAUTION: DO NOT IMMERSE PADS IN ANY LIQUID

CAUTION: DO NOT USE BLEACH OR PHENOLICS ON PADS

INSTRUCTIONS FOR USE

5 List of Applicable Standards:

Sl. no	Standards	Description
1.	EN 62366-1	Medical devices - Part 1: Application of usability engineering to medical devices
2.	EN ISO 14971	Medical devices- Application of risk management to medical devices.
3.	EN 1041	Information supplied by the manufacturer of medical devices
4.	EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
5.	EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
6.	IEC 60601-2-46	Medical electrical equipment - Part 2-46: Particular requirements for the basic safety and essential performance of operating tables
7.	ISTA	International Safe Transit Association standards for package testing

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Allen Medical Systems, Inc. är ett dotterbolag till Hill-Rom, Inc.