Hillrom

WELCH ALLYN pocket LED diagnostic set

Instructions for use

℞ ONLY

REF 775111, DIR 80029234 Ver. B 2022-08

928 & 928+ Series | # **901026** DIAGNOSTIC SET
228 & 228+ Series | # **901080** OTOSCOPE, POCKET LED
128 & 128+ Series | # **901082** OPHTHALMOSCOPE, POCKET LED

Baxter, Hillrom, and Welch Allyn are trademarks of Baxter International, Inc., or its subsidiaries.

CaviWipe is a trademark of Metrex Research Corp.

Welch Allyn, Inc. ("Welch Allyn") assumes no responsibility for any injury to anyone that may result from (i) failure to properly use the product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual, or (ii) any illegal or improper use of the product.


Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
hillrom.com

Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

EC REP and EU IMPORTER
Welch Allyn Limited
Navan Business Park, Dublin Road,
Navan, Co. Meath, C15 AW22
Ireland

Authorized Australian Sponsor
Welch Allyn Australia Pty Limited
1 Baxter Drive
Old Toongabbie NSW 2146
Australia

Authorized Representative for Kazakhstan
Too Orthodox Pharm
Uly Dala Avenue 7/4, apt 136, Nur-Sultan 010000
Kazakhstan

ENGLISH

Intended usd
The WELCH ALLYN pocket LED diagnostic set is intended for use by clinicians and medically qualified personnel for pediatric and adult patients of any gender in both ambulatory and acute care settings.
The WELCH ALLYN pocket LED otoscope and associated accessories are intended to be used for transmitting light into the ear canal to assess the health of the ear and support diagnosis of ear pathologies.
Patients undergoing routine health assessments or minor procedures involving the ear.
The WELCH ALLYN pocket LED ophthalmoscope is intended for viewing the fundus, supporting the diagnosis of ocular structure conditions, and viewing ocular abnormalities that may be indicative of conditions elsewhere in the body. The device also is intended for illuminating the cornea and iris to detect foreign bodies in the cornea and irregularities of the pupil.
The ophthalmoscope is intended for use by clinicians and medically qualified personnel for pediatric and adult patients of any gender in both ambulatory and acute care settings.

Indications for use
WELCH ALLYN otoscopes are hand-held, battery-powered devices that are intended to be used by trained clinicians during a patient health assessment to inspect the external ear, ear canal and tympanic membrane under illumination and magnification. Secondary uses of the otoscopes are to provide general illumination of the throat and/or nasal cavities.
The ophthalmoscope is used in routine ambulatory settings for routine exams and in acute care physical assessments. The ophthalmoscope is used to view ocular structures and is used to support the diagnosis of ocular structure conditions and for viewing ocular abnormalities that may be indicative of conditions elsewhere in the body.
Mode of Operation: Continuous

Symbols
For information on the origin of these symbols, see the Welch Allyn symbols glossary: welchallyn.com/symbolsglossary

GTIN — Global Trade Item Number

Consult Instructions for Use

Do not dispose of as unsorted municipal waste. Separate collection of batteries.

℞ ONLY — By prescription or order of use by a physician or dentist

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Diagnostic set parts

Diagram number	Description	Diagram number	Description
1a	Otoscope head	9	Lens selection dial
1b	Speculum	10	Diopter display window
2	Ophthalmoscope head	11	Neutral density filter switch (Plus model only)
3	2 AA 1.5 Volt non-rechargeable alkaline batteries (batteries not included)	12	Aperture selection dial
4	Handle	13	Lens anchor
5	On/Off switch	14	Lens
6	Handle bumper	15	Lens bumper with bottom slot
7	Brow bumper	16	Lens tab
8	Viewing window		

⚠ WARNINGS
The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death.
WARNING Patient injury risk. No modification of this equipment is allowed. Modifying this device can affect patient and operator safety and can compromise product performance and accuracy, and void the product warranty.
WARNING Patient injury risk. Use only Welch Allyn approved accessories according to the manufacturer's instructions for use. Using unapproved accessories with the device can affect patient and operator safety and can compromise product performance and accuracy, and void the product warranty.
WARNING Patient injury risk. The cleaning and disinfecting procedures and materials described must be conducted by persons trained and familiar with medical device cleaning and disinfecting.

⚠ CAUTIONS
The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.
CAUTION Failure to follow these instructions may cause damage to the device.
CAUTION Device damage risk. Do not sterilize the device.
CAUTION Device damage risk. Do not immerse the device in any solution.
CAUTION Device damage risk. Only authorized Welch Allyn Service Centers should perform service.
CAUTION These devices may be affected by electrostatic discharge (ESD) events ≥ 15kV. We recommend a dry environment with floors made of wood, concrete or other natural materials. If floor coverings are synthetic, relative humidity should be ≥ 30%.

Residual risk
This product complies with relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential patient or user harm from the following:
- Harm or device damage associated with electro-magnetic hazards,
- Harm from mechanical hazards,
- Harm from device, function, or parameter unavailability,
- Harm from misuse error, such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction.

Handle instructions
Attach handle bumper
1. Align the handle bumper ❻ and handle ❹ curves.
2. Slide the handle bumper onto the handle.

Insert the batteries
1. Obtain 2 AA non-rechargeable alkaline batteries.
2. Insert the batteries with the positive (+) terminals ❸ towards the head.

Attach the head to the handle
1. Align the threads on the head (❶ or ❷) and handle ❹.
2. Rotate the head clockwise until the head and handle are firmly joined.

Turn the handle on or off
1. Grip handle ❹.
2. To turn the handle ON, slide the switch ❺ towards the bottom of the handle.
3. To turn the handle OFF, slide the switch towards the top of the handle.

Otoscope
Depth of field measured from the tip of a 4mm speculum is 11 to 25mm.
Magnification of viewing lens is 3,64 +/- 5%.

⚠ WARNINGS
WARNING Patient injury risk. Do not shine otoscope in patient eyes. Shining light in the eye may cause eye injury.
WARNING Patient injury risk. Do not use otoscope (❶b without attaching a speculum (❶b. Using the otoscope without attaching a speculum may cause ear injury or cross-contamination.
WARNING Cross-contamination risk. Do not re-use disposable ear tips or disposable specula. Reusing disposable ear tips or disposable specula may result in contamination from one patient to another patient.

Attach the lens bumper
1. Rotate the lens (❾) clockwise or counter-clockwise to release the lens from the lens anchor (⓭).
2. Slide the lens bumper slot (⓯) onto the bottom lens tab (⓰).
3. Secure the lens and lens bumper by rotating the lens and lens bumper under the anchor.
4. Secure the top half of the lens bumper to the lens by pressing the lens bumper firmly against the edge of the lens.

Ophthalmoscope

⚠ WARNINGS
WARNING Patient injury risk. The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity for more than 9 minutes exceeds the safety guideline.

Ophthalmoscope operation
1. Turn the aperture selection dial (❷, with continuous rotation, to select small, medium, large, fixation, slit, cobalt blue, and red free apertures.
2. Move the lens selection dial (❾) in either direction to focus. There are 19 total lenses, ranging from -20 to +20 diopters.
3. Move the two-position filter switch (⓫ (Plus model only) to engage and disengage the neutral density filter. The filter reduces the light output flux by 48%, +/- 5%.
4. Use the diopter display window (⓱) to view direct diopter values (positive in green, negative in red).

Maintenance, inspection, and storage
Cleaning and disinfection

⚠ **WARNING:** Patient injury risk. The cleaning and disinfecting procedures and materials described must be conducted by persons trained and familiar with medical device cleaning and disinfecting.
Welch Allyn has validated the instructions for preparing the Pocket LED Diagnostic Set for re-use. The user must ensure that the cleaning and disinfection as actually performed by the user's personnel, with the user's equipment and materials, achieves the desired result. This may require validation and routine monitoring of the user's actual process.

Clean
1. Clean all exterior surfaces with a commercial healthcare wipe containing 70% isopropyl alcohol or a CaviWipe.
2. Inspect the device for any visible soil.
3. Repeat Steps 1 and 2 if visible soil remains.

Use a cotton swab with 70% isopropyl alcohol to remove residual film on instrument windows.

To clean the Pocket LED:
1. Rinse the Pocket LED with warm tap water for one minute.
2. Soak the Pocket LED in 70% isopropyl alcohol for 20 minutes.
3. Rinse Pocket LED with tap water and air dry.

Disinfect
1. Disinfect the device using a new 70% isopropyl alcohol wipe or a CaviWipe.
2. Follow the wipe manufacturer's instructions for appropriate use, contact times, precautions, and warnings.

Inspection
Routinely inspect the device for wear or other damage. Do not use if you see signs of damage, if the instrument malfunctions, appears not to be working properly, or if you notice a change in performance. Contact Hillrom's Technical Support department for assistance.

Storage
Store the device according to facility guidelines to keep device clean, dry, and ready for service.

Disposal
Users must adhere to all federal, state, regional, and/or local laws and regulations as they pertain to the safe disposal of medical devices and accessories. If in doubt, the user of the device should first contact Hillrom Technical Support for guidance on safe disposal protocols.

PATENT/PATENTS
hillrom.com/patents
May be covered by one or more patents. See above Internet address. The Hill-Rom companies are the proprietors of European, US, and other patents and pending patent applications.

Hillrom Technical Support
For information about any Hillrom product, contact Hillrom Technical Support at:
hillrom.com/en-us/about-us/locations

Notice to users and/or patients in EU
Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Specifications
Environmental
Transport: -40 °C – +70 °C, 95% R.H. max.
Storage: -10 °C – +55 °C, 95% R.H. max.
Operating: 10 °C – 40 °C, 90% R.H. max., 700hPa – 1060hPa altitude

Standards and compliance
The device complies with the following standards:
IEC 60601-1, 60601-1-2, 60601-1-6, Otoscopo e Oftalmoscopio; ISO 15004-1, 15004-2, 10942 Group B, Ophthalmoscopic; ISO 10993-1, Otoscope
Country-specific standards are included in the applicable Declaration of Conformity.

Lot code
For lot code YY-JJJ, YY= Last 2 digits of the year, and JJJ=consecutive day of the Julian year.

Guidance and manufacturer's declaration
For information about electromagnetic compatibility (EMC), see the Welch Allyn website:
welchallyn.com/emc-pocket
A printed copy of EMC Emissions and immunity information can be ordered from Welch Allyn for delivery within 7 calendar days.

Accessories

Reorder number	Description
05928-U	Soft Case
410521	Hard Case
52432-U	2.75mm Single Use Pediatric Tip, Case
52434-U	4.25mm Single Use Adult Tip, Case
106081	Accessory Kit
21103-U	Tube of assorted Ear Tips
24400-U	Set of Reusable Ear Specula (2.5mm, 3mm, 4mm, 5mm)

Limited warranty
Welch Allyn warrants that: 1) the Pocket LED Ophthalmoscope and Otoscope will be free of defects in workmanship or material for 3 years from the date of purchase, and 2) the Pocket Plus LED Ophthalmoscope and Otoscope will be free of defects in workmanship or material for 5 years from the date of purchase. Welch Allyn will correct any of these defects at its expense by repairing or replacing the product. This warranty does not cover damage caused by accident, repairs performed by a provider who is not authorized by Welch Allyn, or failure to follow Welch Allyn's maintenance instructions.
THIS LIMITED WARRANTY IS PROVIDED INSTEAD OF ANY OTHER IMPLIED PRODUCT QUALITY WARRANTIES AND IT IS THE SOLE WARRANTY OF THE PRODUCTS.

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