



Hillrom™

# ACCELERATE WOUND HEALING AND ELEVATE YOUR STANDARD OF CARE

## AS YOUR PATIENTS WORK TOWARD HEALING, ARE COMPLEX WOUNDS WORKING AGAINST THEM?

Advanced pressure injuries are all too common in hospitals and long-term care facilities — as are complex skin conditions like flaps, grafts, burns and deep tissue injuries. These serious complications can have significant impacts on clinical and financial outcomes.<sup>1</sup>

Pressure injuries are associated with:<sup>1</sup>



**57% LONGER  
LENGTH OF STAY**

↑ 3X higher mortality rate

↑ 22% higher readmission  
rate within 30 days

## For Better Healing, Rethink Your Surface

It's time to elevate your standard of care — and the clinical and financial outcomes that come along with it. We developed the Envella® bed, our most advanced wound care surface to date, to help your patients heal more quickly and comfortably.

Hillrom is a part of Baxter

**Baxter**



## Elevate Your Standard of Care

The Envella® Bed makes complex wound care simpler.



- 1** Easy-to-use touch controls  
Integrated scale and auto-leveling bead bath  
Head-of-bed angle indicator and alert  
Bed exit with alert silence
- 2** Central brake and alert
- 3** Superior performance for all major skin risk factors with Air Fluidized Therapy
- 4** Comfortable and patient-adjustable transition  
Weight-based pressure redistribution
- 5** Easy to remove, self-standing siderails
- 6** Side transfer for easy patient ingress/egress

### TECHNICAL SPECIFICATIONS

Bed	Low Position	22" (56 cm)	Air Fluidized Therapy	Bead Temperature Control Range	86° - 100°
	High Position	31.5" (80 cm)		Minimum Width of the Bead Section	19" (48.3 cm)
	Overall Length	92.5" (235 cm)		Minimum Length of the Bead Section	45" (114 cm)
	Overall Width	40" (101.6 cm)	Caster and Braking System	Caster Size	6"
	Head Elevation Range	0° - 60°		Four Wheel Braking	Yes
	Bed Frame Ground Clearance	6" (15.2 cm)		Break Not Set Alert	Yes
	Patient Weight Range	70-350 lbs (32-159 kg)	Power Requirement	Power	120V, 60 Hz
	Maximum Patient Height	75" (191 cm)		Max Current	10 A

Elevate your standard of care today. We can help. Get started at [1-800-638-2546](tel:1-800-638-2546) or visit [rental.hillrom.com](http://rental.hillrom.com).

[hillrom.com](http://hillrom.com)

<sup>1</sup> Braden B. "Cost of Pressure Ulcer Prevention. Is it really cheaper than treatment?" NPUAP. 2012.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.