

THE CHALLENGING MED-SURG ENVIRONMENT REQUIRES A BALANCE OF KEEPING PATIENTS SAFE AND SATISFIED.

That's why Hillrom's family of surfaces is designed with innovative technology to help manage pressure, shear and moisture to deliver optimal wound prevention and healing.

The Hillrom™ pro+ surface provides advanced technology that delivers superior pressure redistribution and microclimate management.

Pressure injuries are prevalent and costly.

IN THE AVERAGE 100-BED FACILITY,

170

PATIENTS WILL ACQUIRE A PRESSURE INJURY.

Results: Your facility could bear costs of \$1.8 million or more per year.^{1,2,3}



HILLROM™ PRO+ SURFACE

- Focused Advanced Microclimate® therapy is proven to effectively remove excess heat and moisture away from the patient's skin, especially through the sacral area, where pressure injuries are more common.4
- Customized non-powered pressure redistribution adjusts to patient weight and movement for optimal pressure redistribution, shear management, and patient comfort
- Cross contamination risk may be reduced through the internal pump design.
- Internal pump design can help improve caregiver efficiency and therapy delivery.
- Patient comfort levels can increase as compared to similar low air loss surfaces.5.
- Quiet operation
- Sloped heel section can off load pressure in vulnerable heel area.
- Intuitive status indicator provides quick identification that Focused Advanced Microclimate® Therapy is being delivered.

TECHNICAL SPECIFICATIONS

Width	35.5" (rental/purchase), 38.5" (purchase)
Length	80" or 87.5"
Height	8"
Patient Weight	70 – 500 lbs
Seams	Welded
Surface Weight	50 lbs
Frame Compatibility	Step-deck and flat-deck
Warranty	Top Cover: 1 Year Blower Unit: 1 Year Internal Components: 5 Years

Flammability Standards

Cal TB-129, State of California, Bureau of Home Furnishings and Thermal Insulation Technical Bulletin Number 129, entitled "Flammability Test Procedure for Mattresses for Use in Public Buildings," dated October 1992; 16CFR1632 (1-1-14 Edition); NFPA 101-2012, Life Safety Code; BFD IX-11, Mattress Fire Test, dated September 1, 1993 (Boston); 16CFR1633 (1-1-14 Edition); CAN/CGSB-4.2 No. 27.7-2013, Textile Test Methods, Combustion Resistance of Mattresses - Cigarette Test

PROVEN SURFACE PERFORMANCE


Dry Flux* Measured in watts per square meter (W/sqm)



Hillrom™ pro+ surface maximizes heat removal

Evaporative Capacity*


Measured in grams of H₂O removed per square meter per hour (g/m²-hr)



Hillrom™ maximizes moisture removal

Peak Sacral Pressure**

Measured in millimeters mercury (mmHg)



Hillrom™ pro+ surface minimizes peak sacral pressure

*30 (use degree symbol) head-of-bed elevation ** Data weighted for head-of-bed elevation angle frequency.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

¹²⁰¹⁴ AHA Data.

²2013 IPUP™ Survey Results. Data on file.

³ Jackson, et al. "Pressure Ulcer Prevention in High Risk Postoperative Cardiovascular Patients" Crit Care Nurse. 2011;31:44-53.

⁴ Black J, et al. Pressure ulcer incidence and progression in critically ill subjects: Influence of a low-air loss mattress verses a powered air redistribution mattress. J Wound Cont Nurs. 2012;59(5)267-275.

⁵ Hillrom Internal Test Results. Data on file.