



CENTRELLA SMART+ BEDCONNECTED TO THE CENTER OF PATIENT CARF

STANDARD FEATURES AVAILABLE ON BASE MODEL **Designed to be**:

SAFE

- Enhanced 3-mode bed exit with Silence and Auto Reset
- Brake Not Set Alert
- Foley Position Limit & Drainage Bag Holder
- Integrated Scale
- Head Angle <30 Limit Control
- Motion Activated Nightlight
- Emergency CPR
- Trendelenburg/Reverse Trendelenburg
- Integrated Restraint Holders
- Falls Guarantee Program*

CONNECTED

- WiFi Module
- Smart Bed Ready (**SideCom** Required)

SIMPLE

- Advanced **SlideGuard** Technology
- One-button Positioning; Patient Stand Assist,
 FlexAfoot powered bed extension, dining chair,
 and Boost feature
- Graphical Caregiver Interface (GCI)
 Touchscreen
- Battery Backup
- Device Storage Area
- Integrated Urinal Holders
- Four Corner Brake & Steer
- Built-in Line and Cord Management
- Four Corner IV Pole Sockets
- Lockout Controls
- Built-in Help Menus

^{*} Customer must meet eligibility and all other terms and conditions



ONE SCALABLE CENTRELLA PLATFORM

BASE FRAME	PACKAGE 255 OPTION	PACKAGE 355 OPTION
Standard Features (listed on front page)	Advanced Patient Pendant	Centrella max Powered Air Surface
Scalability	USB Charging	SafeView+ System with IllumiGuide Siderail Handgrip
	Contact-Free, Continuous Monitoring Capable	Second Graphical Caregiver Interface (GCI) Touchscreen
	Verbal Safety Prompts	

TECHNICAL SPECIFICATIONS		
Overall frame width	40" or 43" (Comfort Edition) (101.6 cm or 109.22 cm)	
Overall frame length with foot extension	86.7" - 98.7" (220 - 251 cm) (retracted - extended)	
Bed height maximum	32.5" (83 cm)	
Bed height minimum	14.5" (37 cm)	
Head section inclination	0° to 65°	
Thigh section inclination	0° to 30°	
Foot section inclination	0° to 25°	
Trendelenburg/Reverse Trendelenburg position (max)	15°	
Caster size type	5" (13 cm) Central Locking System	
Safe working load	650 lbs (295 kg) (includes patient weight, surface and accessories)	
Patient weight range	70 – 500 lbs (32 – 227 kg)	
Centrella core Surface 35.5"W x 75.5" - 87.5"L x 6"H (90 cm x 192-222 cm x 15 cm)		
Centrella pro Surface 35.5" or 38.5"W x 75.5" - 87.5"L x 7"H	[90 or 98 cm x 192-222 cm x 17.78 cm]	
Centrella pro+ Surface 35.5" or 38.5"W x 75.5" - 87.5"L x 9"	H (90 or 98 cm x 192-222 cm x 23 cm)	
Centrella max Surface 35.5" or 38.5"W x 75.5" - 87.5"L x 8"H	[90 or 98 cm x 192-222 cm x 20.32 cm]	

A LA CARTE OPTIONS

SAFE

Innovative Support Surfaces:

- Centrella max Powered Air Surface (Advanced Microclimate Technology with Integrated Pump, Enhanced Turn Assist, Patient Comfort Controls, and optional Continuous Lateral Rotation)
- Centrella pro+ Surface (Advanced Microclimate Technology with Integrated Pump)
- Centrella pro Surface (Targeted and Customized Non-powered Pressure Redistribution)
- **Centrella** core Surface (Non-powered Pressure Redistribution)

SafeView+ Lighting System (with IllumiGuide Siderail Handgrip)
Contact-Free, Continuous Monitoring
Verbal Safety Prompts
Patient Pendant
Experience POD Device

Voalte Patient Safety

SIMPLE	
Comfort Edition 40" width	
USB charger	
X-ray Cassette Sleeve	
IntelliDrive Powered Transport	
Obstacle Detect System	
Second Graphical Caregiver Interface (GCI) Touchscreen	
Accessory Outlet	
Additional Line Management	

CONNECTED

Cable-Free Nurse Call Connection (Hillrom ReadyConnect)

Active Nurse Call Button

Universal lighting controls, TV controls

SmartCare Remote Management Wireless Bed Connectivity (Visualization and Clinical Application Integration options)

Baxter.com

Baxter, Advanced Microclimate, Centrella, Experience POD, FlexAfoot, Graphical Caregiver Interface (GCI), Hillrom, IllumiGuide, IntelliDrive, NaviCare, Obstacle Detect, ReadyConnect, SafeView, SideCom, SlideGuard, and Voalte are trademarks of Baxter International, Inc. or its subsidiaries.