Baxter Welch Allyn RetinaVue

WELCH ALLYN RETINAVUE 700 PEDESTAL

Convert the handheld **Welch Allyn RetinaVue** 700 Imager into a tabletop format with ease using the **RetinaVue** 700 Pedestal.

- Multiple adjustment points for flexible patient positioning
- Quick release enables "take it and go" portability for mobile exams
- Aids in stabilizing across all patient populations
- Minimal to no training required across staff
- Ideal for high-volume clinics



PRODUCT SPECIFICATIONS	
Weight	8.4lbs (3.8kg)
Height	5in (12.8cm)
Length	14.5in (36.8cm)
Width	11in (27.9cm)
Raising Angle	70°

ORDERING INFORMATION

728134

Pedestal Base for **RetinaVue** 700 Imager

FOR MORE INFORMATION, PLEASE VISIT **RETINAVUE.COM**.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Hillrom, RetinaVue and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC147-230077 V1 06/23