Document Description:	Welch Allyn Service Tool v1.10.0 Software Release	Document Number:	80027878
	Notes	Version:	A
	Welch Allyn ™		
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.			

Customer Service Bulletin

Product:	Welch Allyn Service Tool		Date:	2021-01-26
Subject:	Software Release Notes v1.10.0			
HW Version(s) Affected:	All	SW Version(s) Affected:	All	
Serial Numbers Affected:	All	Lot or Date Code Affected:	N/A	

Classification:	As Needed			
Distribution:	⊠ Custom	ner Care	□ Product Service □	⊠ Field Service
	⊠ ASPs	□ Distributors	□ Customers	☐ Company

Training Required:	⊠ Yes □ No		
Summary:	The below summarizes the changes between the previous installer release (WAST v1.9.0.0) and the current release (WAST v1.10.0). This includes changes made to the Core and Plugins, and all new device Plugins.		
	Items addressed are listed below:		
	1. Removed the Axeda Agent from the Installer package as this is no longer supported. When installing one must install the Service Tool then install the PartnerConnect agent		
	2. Updated the Calibration Certificate Template to address the Australia Service Center new address.		
	3. Removed Windows 7 support as this is considered end of support from Microsoft.		
	4. Updated application with security changes to minimize security risk.		

Document Description: Welch Allyn Service Tool v1.10.0 Software Release	Document Number: 80027878	
Notes	Version: A	
Welch Allyn ™		
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.		

Hardware/Software Compatibility

The below summarizes the Windows OS systems compatibility for the Welch Allyn Service Tool v1.10.0

Device Description	Versions
Tested Operating Systems	Win 10 x86, Win 10 x64

Service Action Needed: All Welch Allyn Service Tool users should update to this latest version in order to get the most recent changes.

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
Α	N/A	Initial Release	2021-01-26	KMG

THIS INFORMATION IS THE PROPERTY OF WELCH ALLYN, INC. AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT OR DEVICES WITHOUT THE EXPRESS WRITTEN PERMISSION OF WELCH ALLYN, INC.