

Document Number: 80029669

Document Version: A

Printed or electronic versions of this document not accessed directly from the designated Hillrom Controlled Quality Information System are For Reference Only.

Customer Service Bulletin

Product:	Welch Allyn® Connex Vital Signs Monitor (CSM)			
Subject:	Masimo LNCS® SpO2 Patient Cable			
HW Version(s) Affected:	ALL	SW Version(s) Affected:	ALL	
Serial Numbers Affected:	ALL	Lot or Date Code Affected:	ALL	

Classification:	Information Only		
Distribution:	⊠ Customer Care	☑ Product Service	⊠ Field Service
	□ ASPs □ Distributors	□ Customers	☐ Company Confidential

Summary:

The Welch Allyn® Connex Spot Monitor operates with a Masimo LNCS® Series 20-pin SpO2 Patient Cable. Due to global supply constraints, Baxter will be shipping a 10-foot Masimo LNCS® Series SpO2 patient cables (RED -LNC-10) with Connex Spot Monitors, in place of the Masimo 4-foot LNCS® Series SpO2 patient cables (RED LNC-04).



Once supply constraints improve and stock again is available for the 4-foot cables, we will return shipping the Connex Spot Monitor with those again (there is no specific timeline at this time). Due to this, some customers may receive mixed shipments of 4 and 10-foot LNCS® Series cables. Unfortunately, Baxter is unable to offer the 10-foot Masimo LNCS® Series Patient Cables to customers that are affected by this at no charge.

END OF BULLETIN

THIS INFORMATION IS THE PROPERTY OF HILLROM HOLDINGS AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS
FOR THE MANUFACTURE OR SALE OF EQUIPMENT, DEVICES, OR SOFTWARE WITHOUT THE EXPRESS WRITTEN PERMISSION OF HILLROM.



Document Number: 80029669

Document Version: A

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
Α	N/A	Initial Release	2022-03-10	TP