

Document Number: 20017180

Document Version: A

Printed or electronic versions of this document not accessed directly from the designated Hillrom Controlled Quality Information System are For Reference Only.

Customer Service Bulletin

Product:	CSM – CONNEX SPOT MONITOR CVSM – CONNEX VITAL SIGNS MONITOR CIWS – CONNEX INTEGRATED WALL SYSTEM						
Subject:	MASIMO SP02 MODULE REPLACEMENT						
HW Version(s) Affected:	CVSM - P3,P5,P6 CSM - ALL	SW Version(s) Affected:	ALL				
Serial Numbers Affected:	ALL	Lot or Date Code Affected:	ALL				
Γ							
Classification:	Information Only						
Distribution:	□ Customer Care	☐ Product Service	☐ Field Service				
		□ Customers	☐ Company Confidential				
Summary:							
THE MASIMO SP02 MODULE MODEL MX5 HAS BEEN REPLACED WITH MASIMO SP02 MODULE MODEL							
MX7 DUE TO COMPONENT OBSOLESCENCE. THERE IS NO CHANGE TO THE FORM, FIT, OR FUNCTION							
OF THE SP02 MODULE.							
01 1112 31 02 WIOD							
Bulletin Detail 1	L:						
CSM: SERVICE KITS 105931 and 105950 HAVE BEEN UPDATED WITH THE REPLACEMENT MX7 MODULE.							
CVSM: SERVICE KIT 107282 HAS BEEN UPDATED WITH THE REPLACEMENT MX7 MODULE.							
CIWS: SERVICE KIT 107284 HAS BEEN UPDATED WITH THE REPLACEMENT MX7 MODULE.							

THIS INFORMATION IS THE PROPERTY OF HILLROM HOLDINGS AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT, DEVICES, OR SOFTWARE WITHOUT THE EXPRESS WRITTEN PERMISSION OF HILLROM.



Document Number: 20017180

Document Version: A

Table 1		Compatible Vers	sions			
Subsystem Software Hardware Welch Allyn Part Number						
Description	Version	Version	v elem ming in a une i (unisser			
Main Board		P3, P5, P6	407574-06, 407574-08			
Braun 6000	5.25	N/A				
Software						
Braun 4000	1.00.03	N/A				
Software Lammar Comm's	1.00.00	P2,P3	407524			
module	A0003	1 2,1 3	407324			
Newmar Comm's	V2	P2,P3	410389-04			
module						
Lammar Radio	3.00.04	MAC address	406187			
Software		starting				
Newmar Radio	2.00.01	with:00:1A:FA MAC address	413023			
Software	2.00.01	starting	413023			
Software		with:00:17:23				
Masimo SPO2	7.10.1.2	0x701	107282, 107284			
Nellcor SPO2	1.00.13,	P3,P4	405712, 412896, 107020, 406273-01			
	1.00.15	- 7	, , , , , , , , , , , , , , , , , , , ,			
EARLY Sense	01.06.007	P3	725692			
SureTemp Plus	2.08	P3	405701			
NIBP	2.05	Rev A.	405672			
ECG	1.00.05	P3	412726			
CO2 Oridion	1.10.02	P3	408259			
CO2 Medtronic	1.14	P3	412897			
CO2 Medtronic	1.16	P3	420502			
Scale Software	1.0.26	N/A	N/A			
WA Connectivity	2.4	N/A	N/A			
SDK	5002	DT/A	NY/A			
NCE Connex CS	5.0.0.3	N/A N/A	N/A N/A			
Configuration Tool	1.8.4	N/A	N/A			
DCP Server	3.0.0.3	N/A	N/A			
Data Catcher	4.0.0.20	N/A	N/A			
Connex VM	2.0.1	N/A	N/A			
Welch Allyn Service Tool	1.10.x.x	N/A	N/A			
Welch Allyn Service Hub	5.0.x	N/A	N/A			

THIS INFORMATION IS THE PROPERTY OF HILLROM HOLDINGS AND AS SUCH SHALL NOT BE REPRODUCED, COPIED, OR USED AS A BASIS FOR THE MANUFACTURE OR SALE OF EQUIPMENT, DEVICES, OR SOFTWARE WITHOUT THE EXPRESS WRITTEN PERMISSION OF HILLROM.



Document Number: 20017180

Document Version: A

END OF BULLETIN

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
Α	N/A	Initial Release	2022-03-16	CGH