


Document Description: CSB – Discontinued Masimo LNOP Sensors	Document Number: 80021708 Version: B
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

Customer Service Bulletin

Product: Masimo LNOP Sensors & Cables	Date: 2016-08-10
Subject: CSB – Discontinued Masimo LNOP Sensors & Cables	
HW Version(s) Affected: See list below	SW Version(s) Affected: NA
Serial Numbers Affected: NA	Lot or Date Code Affected: NA


Classification: Informational Only
Distribution: <input checked="" type="checkbox"/> Customer Care <input checked="" type="checkbox"/> Product Service <input checked="" type="checkbox"/> Field Service <input checked="" type="checkbox"/> ASPs <input checked="" type="checkbox"/> Distributors <input checked="" type="checkbox"/> Customers <input type="checkbox"/> Company Confidential

Training Required: Yes No

Summary: Masimo LNOP sensors will no longer be available from Welch Allyn after current stock is depleted. Customers with Spot Vital Signs, Spot LXi Vital Signs, VSM300 Vital Signs, Micropaq, Propaq CS, or Propaq Encore with a Masimo SpO2 option using LNOP sensors and cables can continue using these devices with the LNCS line of sensors and cables. Compatible replacement LNCS sensors and cables can be found on the product page accessories tab at <https://www.welchallyn.com>. When replacing the LNOP sensors with LNCS sensors a replacement LNCS cable is required to connect the sensor with the monitor.

Discontinued Materials:
These material will no longer be available from Welch Allyn once current stock is exhausted.

Masimo PN	WelchAllyn PN	Description	Where used
1001	LNOP-ADT	Adhesive Finger Sensor - Adult, case of 20	VSM300/Spot/Spot LXi
1001	008-0655-00	LNOP Adult SpO2 Sensors, Disposable, box of 20	VSM300/Spot/Spot LXi
1002	008-0694-00	LNOP Neonate Sensors, Disposable, box of 20	VSM300/Spot LXi

Document Description: CSB – Discontinued Masimo LNOP Sensors	Document Number: 80021708 Version: B
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

1005	PC-08	8-foot Cable with DB-9 connector for LNOP	VSM300/Spot LXi
1005	008-0824-00	8-FOOT Cable W/Sensor Connector (LNOP)	VSM300/Spot LXi
1006	PC-12	12-foot Cable with Sensor Connector (LNOP)	Spot LXi
1025	LNOP-PDT	Adhesive Finger Sensor - Pediatric, case of 20	VSM300/Spot/Spot LXi
1025	008-0654-00	LNOP Pediatric/Slender Digit SpO2 Sensors, Disposable, box of 20	VSM300/Spot/Spot LXi
1173	PC-04	4-foot Cable with Sensor Connector for LNOP	VSM300/Spot LXi
1173	008-0954-00	LNOP Interface Cable, SPO2, 4 foot (14 pin connector)	VSM300/Spot LXi
1269	LNOP-DCI	LNOP-DCI Adult Reusable Sensor	VSM300/Spot/Spot LXi
1269	008-0648-00	LNOP-DCI Adult Reusable Sensor	VSM300/Spot/Spot LXi
1276	LNOP-DCIP	Reusable Finger Sensor - Pediatric	VSM300/Spot/Spot LXi
1276	008-0656-00	LNOP-DCIP Pediatric Reusable Sensor	VSM300/Spot/Spot LXi
1794	008-0953-00	LNOP-TCI Ear SPO2 Sensor, Reusable	VSM300/Spot/Spot LXi
1800	LNOP INF-L	Adhesive Finger Sensor - Infant, case of 20	VSM300/Spot/Spot LXi
1800	008-0951-00	LNOP-Inf-L Neonate Long, Disposable, box of 20	VSM300/Spot/Spot LXi
2288	PC-04-WA	4-foot Cable with DB-9 connector for LNOP	Spot
2289	PC-08-WA	8-foot Cable with DB-9 connector for LNOP	Spot
1002/1798	LNOP-NEO	Adhesive Finger Sensor - Neonate, case of 20	VSM300/Spot LXi

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	D*	GCS
B	Header	Rev A to Rev B Admin change to allow PDF to be added to SAP DIR	D*	GCS