



Hillrom™



TS 7500 OR Table with floor-mounted column and SQ14 X-TRA tabletop

SIEMENS INTEGRATED TABLETOPS

SPECIALTY HILLROM™ TABLETOPS FOR SIEMENS ARTIS PHENO® DEVICE

Integrate with Hillrom tabletops to:

- Achieve imaging speed equivalent to the Siemens ARTIS pheno® device.
- Operate tabletop via Siemens control or Sensor Control FloatLine for traditional control style.
- Facilitate collision protection between tabletop and Siemens ARTIS pheno® device.
- Utilize large 360 radiolucent imaging window and carbon fiber components.
- Support capabilities with DynaCT.



INTEGRATED TABLETOPS WITH THE SIEMENS ARTIS PHENO® DEVICE

CARBON FLOATLINE™ TABLETOP



Neuro components attach to the tabletop for precise imaging and additional length.



The flat tabletop floats into position, even in combination with reverse Trendelenburg and Trendelenburg.



The Sensor Control FloatLine remote offers the traditional control style of a floating tabletop.

SQ14 X-TRA TABLETOP



The articulating tabletop allows the patient to be properly positioned for spine procedures (left), laterally for VAT procedures (center), and seated for cranial or neuro procedures (right).

TECHNICAL DETAILS

	Carbon FloatLine	SQ14 X-TRA
Material no	# 1616201	# 1581502
Dimensions (L x W) Without siderails With siderails	90.2" x 20.5" (2290 mm x 520 mm) 90.2" x 22" (2290 mm x 560 mm)	94.3" x 21.3" (2395 mm x 540 mm) 94.3" x 23.6" (2395 mm x 600 mm)
Height range	19.3" – 40.9" (490 mm -1040 mm)	19.3" – 40.9" (490 mm -1040 mm)
X-ray area	Head End: 70.9" x 20.5" (1800 mm x 520 mm) with Headrest: 82.7" x 20.5" (2100 mm x 520 mm)	Head end: 36.2" (920 mm) Foot end: 59.5" (1510 mm)
Longitudinal slide	31.5" (800mm)	13.8" (350 mm)
Lateral slide	±4.9" (±125 mm)	±3.5" (±88 mm)
Trendelenburg / reverse Trendelenburg	± 45°	± 45°
Tilt	± 30°	± 30°
Maximum patient weight	Without lateral movement: 496 lbs (225 kg) In all positions positions: 396.8 lbs (180 kg)	352.7 lbs (160 kg)

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

ARTIS pheno is a registered trademark of Siemens Healthcare GmbH

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.