



Hillrom™

## Yellofins Apex™ Stirrups

### PRODUCT SPECIFICATIONS

Patient weight capacity	500 lb
Lithotomy adjustment (from sterile field)	Yes
Lithotomy range (from sterile field)	+80° to -32°
Abduction range (from sterile field)	+23° to -9°
Maximum medial/lateral shift	±20°
Product weight	18.5 lb



[Learn how at hillrom.com.](https://www.hillrom.com)

### ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at [hillrom.com](https://www.hillrom.com).

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

130 E. Randolph St. Suite 1000, Chicago, IL 60601

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Allen Medical Systems, Inc. ALL RIGHTS RESERVED. APR60601 rev 1 26-FEB-2020 ENG - US



Hillrom™



Yellofins Apex™ Stirrups

# YOU ARE ALWAYS RAISING THE BAR ON SAFETY AND EFFICIENCY IN THE OPERATING ROOM.

## SO ARE WE.

The Yellofins stirrups line from Hillrom has long been the gold standard in lithotomy positioning. Now we're taking the gold standard to a new level with the Yellofins Apex™ stirrups—so you can achieve new levels of performance and safety in the operating room.

We're elevating the gold standard to help you elevate care.



### SAFE

With a dual-rod design to prevent medial leg drop and automatic position locking technology to safely position your patient – Apex helps you protect your patients and your surgical team.



### INTUITIVE

The single-point boot release handle allows for ease of adjustment, taking the heavy lifting out of your workflow.



### INFECTION CONTROL

Welded pads and removable silicone straps that are easy to clean and disinfect, allowing you to control potential infections.



### EFFICIENT

Accommodate a wide variety of patients with the stirrups' 500 lb weight capacity.

## Next-Level Innovation for Next-Level Care

In developing the Yellofins Apex stirrups, we listened to surgical teams like yours who told us what's needed to support the future of surgical care. Starting from the ground up, we uncovered every opportunity to help you improve safety and efficiency.

