



Hillrom™



# ADVANCING CONNECTED CARE™ WITH VITAL SIGNS CONNECTIVITY

## Hillrom and Epic

Help improve staff productivity, patient outcomes and risk reduction with accurate, automatic vitals sent directly to your EMR. See how a connected solution can help you improve clinical workflows and reduce errors.

### Unconnected Workflow

1

Take a patient's vital signs

2

Document vitals on paper

3

Open patient chart in Visit Navigator®

4

Type each vital sign recorded on paper into its respective field in the patient's chart

### Connected Workflow

1

Take patient's vital signs

2

Click the "Vitals" button in Visit Navigator

3

Click "Import" on the External Results Entry

4

Validated vitals will appear in Visit Navigator

## Connectivity Counts

Help save time and improve data accuracy by streamlining your clinical workflows.

- Reduce time previously spent manually entering data into the EMR
- Improve accuracy by minimizing manual processes and associated errors
- Work more efficiently by freeing up your time for other patient priorities



Connect the vital signs device to your computer via a USB cord or Bluetooth®.

Take patient vitals. No protected health information (PHI) is entered on device.


Vitals data is transferred over the USB cable to the EMR and is updated in the patient's chart.

## Your Workflow Supported by the Citrix® Environment

### Thin Client

The Welch Allyn Virtual Channel Client (VCC) software enables a remote thin client to utilize USB-connected devices. You will need to purchase a license and installation services to support this feature.

Part numbers: 103945, 106371 and PS-CS-13



Contact your Hillrom representative to learn more.

[hillrom.com](http://hillrom.com)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Bluetooth is a registered trademark of Bluetooth SIG, Inc. Citrix is a registered trademark of Citrix Systems, Inc.

© 2020 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR65904 Rev 1 80023917 Ver A 28-AUG-2020 ENG-US