



Welch Allyn
KleenSpec® Disposable Vaginal Specula
and Cordless Illumination System



An advancement in women's health
technology that's nothing short of brilliant

Welch Allyn KleenSpec® Disposable Vaginal Specula
and Cordless Illumination System—

now in every size
for every type of patient.



590 Series KleenSpec® Disposable Vaginal Specula

Enhance Patient Comfort

- The extra small speculum allows you to perform a comfortable exam with virginal/pediatric patients as well as post-menopausal and post-hysterectomy women—no more referrals or uncomfortable exams using metal specula or the wrong size
- Smooth, molded acrylic designed to prevent pinching or binding of tissue for greater patient comfort

Help Reduce Cross-Contamination

- Single-use speculum helps to reduce cross-contamination for improved patient safety while eliminating the time and expense of reprocessing metal specs

Improve Staff Productivity

- Features a wider handle design for improved maneuverability and better balance during exams
- A strong grip and a yoke that's more flexible than low-cost substitutes

The 590 Series KleenSpec Vaginal Speculum is compatible with the Cordless Illumination System (#80010) as well as the Corded Illumination System (#78900).



Look for the distinctive Welch Allyn blue, green, white and yellow colored handles

“ I had been using the corded version of the Welch Allyn Illumination System since 1980 and was thrilled when the cordless system was introduced. We have been using the Welch Allyn cordless Illumination system with disposable vaginal speculums for all of our patient exams since 2008 and couldn't be happier. This light is brighter than any I have ever seen and shows true tissue color. If you forget to charge the light, it is easy to grab a replacement from another room or use a spare. The light is small, easy and quick to use which truly makes us more efficient. With Manhattan space at a premium, not having floor space taken up by a light source is a great benefit to us. I would recommend this system to any clinicians looking to improve the quality of their care.”

Shelley Kolton, MD / Downtown Women OB/Gyn Associates / New York City, NY



KleenSpec Cordless Illumination System

- LED technology as opposed to halogen, for a whiter, brighter light and true tissue color during an exam
- No bulb replacement necessary for a reduced cost of ownership
- Features a cordless design so you'll never have to worry about cords getting in your way, breaking or needing to be cleaned, speeding the exam and giving you a full range of motion
- Push-button on/off; the illumination system easily fits into every size 590 Series speculum
- Newly designed light pipe provides a homogeneous spot and eliminates glare back so all of the light is projected forward
- Uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps and has a continuous on time of 80 minutes
- A lithium ion battery provides a recharge time of five hours

Ordering Information

KleenSpec Vaginal Specula—Premium 590 Series

- 590XS KleenSpec 590 Series Premium Disposable Vaginal Specula, **X-Small (24/box; 4 boxes/case)**
- 59000 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Small (24/box; 4 boxes/case)**
- 59001 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Medium (24/box; 4 boxes/case)**
- 59004 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Large (18/box; 4 boxes/case)**

KleenSpec Vaginal Specula—Premium 590 Series with Smoke Tube

- 59005 KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Small (12/box; 4 boxes/case)**
- 59006 KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Medium (12/box; 4 boxes/case)**

KleenSpec Cordless Illumination System

- 80010 Complete KleenSpec Cordless Illuminator with Charging Station
- 80015 Same as above, for export purposes—specify country and voltage; conforms to IEC 60601-1
- 80000 KleenSpec Cordless Illuminator
- 74010 Charging Station for Cordless Illuminator
- 59010 Disposable Sheath for Cordless Illumination System **(100/box; 5 boxes/case)**

KleenSpec Corded Illumination System

- 78900 KleenSpec Illumination System (with cord)—4.7 V Halogen Light Handle with on/off switch and transformer. 110-130 V, 60 Hz for North America. UL Listed, CSA Certified
- 78950 Same as above, for export purposes—specify country and voltage; conforms to IEC 60601-1
- 08800-U 4.6 V to 4.7 V Halogen Replacement Lamp for 78900 Illumination System
- 07800-U 6 V Halogen Replacement Lamp for 78010 Illumination System
- 58010 Disposable Sheath for Corded Illumination System **(25/box; 5 boxes/case)**

KleenSpec Vaginal Specula—Premium 580 Series w/Sheath

- 58000S KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Small (18/box; 4 boxes/case)**
- 58001S KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Medium (18/box; 4 boxes/case)**
- 58004S KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Large (12/box; 4 boxes/case)**

KleenSpec Storage Systems

- 78820 KleenSpec Corded Illumination Storage System
- 58014 KleenSpec Vaginal Speculum Dispenser

Green Series™ Exam Light IV—Another LED innovation from Welch Allyn

- 48810 GS Exam Light IV with Table/Wall Mount
- 48950 Mobile Stand for GS Exam Light IV
- 52640 Disposable Sheath for Exam Light III and GS Exam Light IV **(25/box, 5 boxes/case)**



#78900



#58001S
#78900



#48810



Explore the entire Welch Allyn women's health line, including the KleenSpec Series and Cordless Illumination System online at welchallyn.com/womenshealth or call us at **800.535.6663**.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

