



FCS 500

HYBRID MONITOR HOLDER







Hillrom offers a variety of arm systems that allow for flexible arrangement and positioning of a various number of monitors. Options range from one large monitor up to eight smaller monitors to fit your OR needs.



	FCS 500-SI	FCS 500-TP	FCS 500-PH
Monitor Compatibility	Siemens	Third Party (excluding Siemens or Philips)	Philips
Support Arm Lengths	72" (1840 mm), 83" (2105 mm), 93" (2370 mm), 100" (2540 mm)	72" (1840 mm), 83" (2105 mm), 93" (2370 mm), 100" (2540 mm)	63" (1600 mm), 73" (1865 mm), 84" (2130 mm), 91" (2300 mm)
Vertically Articulating Arm Type	Spring	Spring	Motorized
Vertically Articulating Arm Length	1,240 mm	1,240 mm	1,000 mm
Maximum Reach	100" (2540 mm)	100" (2540 mm)	90" (2300 mm)
Monitor Configurations	1 - 8 monitors (see back for configuration types)	1 - 8 monitors (see back for configuration types)	1 monitor up to 56"
Ability to Mount Monitors to Backside	2 small monitors, up to 26"	2 small monitors, up to 26"	2 small monitors, up to 26"
Arm Brake Options <i>*Friction required for monitor bearing brake</i>	Friction or Electro-Pneumatic (gas-assisted spring arm)	Friction or Electro-Pneumatic (gas-assisted spring arm)	Friction, Electro-Pneumatic, or Electro-Magnetic (motorized articulating arm)



FCS 500 CONFIGURATIONS

	 Type 1	 Type 2	 Type 3	 Type 4	 Type 5	 Type 6
# of Monitors	3 Monitors	4 Monitors	2 x 2 Monitors	2 x 3 Monitors	2 x 4 Monitors	56" Monitor
Spring Arm Load Capacity	44-110 lb (20-50 kg)	88-154 lb (40-70 kg)	88-154 lb (40-70 kg)	132-198 lb (60-90 kg)	176-242 lb (80-110 kg)	176-242 lb (80-110 kg)
Monitor Carrier Load Capacity	max. 70 lb (32 kg)	max. 79 lb (36 kg)	max. 79 lb (36 kg)	max. 136 lb (62 kg)	max. 136 lb (62 kg)	max. 158 lb (72 kg)

For more information, please contact your Hillrom sales representative.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.