



Affinity® 4 Birthing Bed

Item Description: Electric Specialty Patient Bed, Birthing

Model Number: P3700

JSN Number: M7000

Warranty: 1 Year From Delivery Date

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

V-Cut or Straight Cut Surface

Removable Foot Section

Battery Backup

Powered Foot Section

Four Wheel Dual Locking Casters

OneStep® Siderail Release

Plastic Deck Panels

Automatic Nightlight

Lock-out Controls

Manual Trend-Like Positioning

Instant CPR

Instant Labor Grips

Automatic Seat Tilt (15°)

Safe Working Load 500 lbs (227 kgs)

Maximum Patient Weight 470 lb (213 kgs)

Siderail Controls

IV Pole

Hinged Foot Section Mattress

Central Brake and Steer

Fluid Basin

Liftoff Bed Weight: 480 lbs (217 kgs)

Stow and Go® Bed Weight: 511 lbs (232 kgs)

Product Accessory Options

☐ Calf Supports (AF-CALF)

☐ Push Handles with IV Hook (AF-PH)

☐ ISS Transfer Pole Assembly

☐ ISS Offset Bar

☐ Standard IV Pole

☐ Vertical Oxygen Tank Holder

Affinity® 4 Birthing Bed

Item Description: Electric Specialty Patient Bed, Birthing

Model Number: P3700

JSN Number: M7000

Proposal Number: _____



Affinity® 4 Birthing Bed Technical Specifications

Bed length	92" (234 cm)
Bed width siderails up	39" (99 cm)
Bed width siderails down	36" (91 cm)
Bed height low (with mattress)	24" (61 cm)
Bed height high (with mattress)	39" (99 cm)
Bed maximum height of seat section in Trend-like position	40" (102 cm)
Mattress length	78" (198 cm)
Mattress width	35" (89 cm)
Mattress thickness (head/seat section)	5" (12.7 cm)
Mattress thickness (foot section)	3.5" (8.9 cm)

Maximum head elevation	63°
Maximum seat elevation	15°
Maximum trend-like positioning	8°
Maximum lift capacity bed	500 lbs (227 kg)
Maximum lift capacity foot section	400 lbs (181 kg)
Caster Size	6" (15 cm)
Bed Weight – Lift Off	480 lbs (218 kg)
Stow and Go	511 lbs (232 kg)
Electrical Information	100V 5A 50Hz 120V 4A 60Hz 230V 2A 50Hz

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
203025 rev 2 22-MAY-2018 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom