



## DELIVERING SAFETY, COMFORT AND CONVENIENCE

The labor and delivery experience can be unpredictable for everyone involved. To help manage this, Hillrom offers a solution that can adapt and deliver for mom, baby and you.

Every mom is different, every labor is different. That's why Affinity® 4 Birthing Bed is designed to help you adapt by delivering safety, comfort and convenience, whatever the situation.

### DELIVERING SAFETY

- Infinite Patient Positioning
- Patient Transfers
- Patient Transports

### DELIVERING COMFORT

- Automatic Pelvic Tilt
- Wide Surface
- Lumbar Support

### DELIVERING CONVENIENCE

- Battery Backup
- Clinical Access
- Infinite Positioning Calf Supports

## DELIVERING SAFETY


The Affinity® 4 Birthing Bed is designed to keep mom, baby and you safe.

## DELIVERING COMFORT

Keeping mom comfortable throughout the birthing experience should be paramount.

## DELIVERING CONVENIENCE

The Affinity® 4 Birthing Bed is designed to make every movement fast and easy.



### Maximum Patient Weight

500 lb bed patient weight limit.  
400 lb foot section patient weight limit.

### Patient Transport

Ergonomic push handles and light-weight birthing bed frame.



### Gas Assist Bed Tilt

Achieve 8° gas assisted bed tilt on demand.

### 15° Pelvic Tilt

Tilted seat designed to cradle mom to increase her feeling of security.



### Patient Transfer

Zero-transfer gap for a seamless transfer of your patient.

### Clinical Access

Available in both traditional U-Cut or our exclusive V-Cut mattress to improve perineal access.



### “Dual-Locking” Casters

Two points of contact prevents casters from rolling or swiveling.



### Lumbar & Seat Support

Customize patient comfort with the push of a button.

### Battery BackUp

Enables function of the bed during transport or loss of power.



### Stow & Go® Foot Section

Exclusive, ergonomic design; foot section slides and stores under the bed.

### Light-Weight Foot Section

Handles and angular stand to keep mattress off the floor.



### Wide Sleep Deck

35" surface width gives mom room to care for baby in bed.

### Infinite Positioning Calf Supports

Ergonomically designed, provides comfort and infinite positioning.





## TECHNICAL SPECIFICATIONS

Product Number	P3700
Full Bed Length	90" (229 cm) P3700A-P3700C 92" (234 cm) P3700D-P3700E
Bed Width Siderails Up Siderails Down	39" (99 cm) 36" (91 cm)
Bed Height Low (with mattress) High (with mattress) Maximum Height of Seat Section in Trend-like Position	24" (61 cm) 39" (99 cm) 40" (102 cm)
Mattress Length Width Thickness (head section) Thickness (foot section)	78" (198 cm) 35" (89 cm) 5" (12.7 cm) 3.5" (8.9 cm)
Critical Angles Maximum Head Elevation Maximum Seat Elevation Maximum Trend-like Positioning	63° 15° 8°
Maximum Lift Capacity Bed Foot Section Caster Size Bed Weight	500 lbs (227 kg) 400 lbs (181 kg) 6" (15 cm) 480 lbs (218 kg)

## STANDARD FEATURES

- V-Cut or Straight Cut Surfaces
- Removable Foot Section
- Battery Backup
- Powered Foot Section
- Four Wheel Dual-locking Casters
- OneStep® Siderail Release
- Plastic Blow-Molded Head Panel
- Automatic Nightlight
- Lock-out Controls
- Manual Trend-like Positioning
- Instant CPR
- Instant Labor Grips
- Automatic Seat Tilt (15°)
- Siderail Controls
- IV Pole
- Plastic Blow Molded Head Panel
- Hinged Foot Section Mattress
- Central Brake and Steer
- Fluid Basin

## OPTIONAL FEATURES

- Prevention Surface
- Air Package
- Nurse Call, Lighting, Entertainment & UTV
- Contemporary Headboard
- Freedom Hill™ Headboard
- Liberty Hill™ Headboard
- Calf Supports
- Push Handles with IV Hooks



## ACCESSORIES

- Universal Instrument Tray
- Labor Bar
- Drainage Bag Hook Kit
- Permanent IV Pole
- Line Managers
- Foot Rest Bar (for lift off beds)
- Pendant Holder
- Pendant Kit
- Full Leg Supports
- Full Leg Supports, Telescoping
- Drapes - Medline/Proxima
- Anesthesia Screen
- Oxygen Tank Holder

For more information, please contact your Hillrom sales representative at 1-800-445-3730.

[hillrom.com](http://hillrom.com)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.