

Baxter

Barton Convertible Chair Series

Early and safe patient mobilization is critical to supporting patient recovery in the acute care and post acute care settings. The **Barton** Convertible Chair can accommodate adult patients of all sizes and multiple options to support a wide range of patient weights, abilities, and morphology.



KEY FEATURES

- **Barton** Convertible Chairs designed to easily convert between a flat and upright mobile chair that allows caregivers to transfer patients.
- Mobile chair featuring unlimited number of reclined and tilt positions.
- Tilt-in-Space capability
- Trendelenburg position capability
- Customized features allow the chair to be set to the patient: adjustable head, leg and arm rests, side supports and height.
- Heavy-duty central locking and steering casters support easy mobility.



OPTIONS

- I-400 (capital only)[†]
- I-660 (rental only)
- I-700 (both rental and capital)[†]

[†]The I-400 and I-700 capital purchase comes with quick release belt.

Barton Convertible Chair



How it works—Positioning & Transfer System (PTS)

- Place the chair at bedside in the horizontal position.
- Attach the sheet to the Positioning & Transfer System (PTS).
- Transfer the patient by turning a handle.
- Incline the chair and the patient is ready to be transported.
- Reverse the process to transfer the patient from chair to bed with ease.

TECHNICAL SPECIFICATIONS

Seat Width	
I-400*	21" (54.6 cm)
I-660*	27" (68.8 cm)
I-700*	25.5" (65 cm)
Overall Width	
I-400*	30" (76 cm)
I-660*	35" (89 cm)
I-700*	35" (89 cm)
Overall Length	
Prone Position	72.5" - 78.5" (184 cm - 199 cm)
Wheelchair Position	45.5" - 62.5" (115 cm - 159 cm)
Back angle	0-70°
Tilt angle	0-17°
Caster Diameter	5.9" (15 cm)
Maximum weight capacity	
I-400*	400 lbs (181 kgs)
I-660*	660 lbs (300 kgs)
I-700*	700 lbs (318 kgs)

* I-400 (capital only); I-660 (rental only); I-700 (capital and rental)



Designed to quickly and easily adjust to Tilt-in-Space positions.



Fully-adjustable headrest with side support.



Heavy-duty central locking/steering casters.



The optional **Barton** PTS helps transfer your patient with ease.



Manual height adjustability.

Distributed by:

Baxter Healthcare Corporation
One Baxter Parkway / Deerfield, Illinois 60015
Baxter.com / T 800.445.3730

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

Baxter.com

Baxter is a trademark of Baxter International Inc. or its subsidiaries.
Barton is a registered trademark of Human Care HC AB.
US-CS402-250002 (v1.0) 07/2025