



**Corporate Headquarters**

4341 State Street Road, P.O. Box 220  
Skaneateles Falls, NY 13153-0220 USA  
Phone: 800.535.6663 Fax: 315.685.3361  
[www.welchallyn.com](http://www.welchallyn.com)



November 21, 2010

To Whom It May Concern:

The measurement time specification for the Braun PRO 4000 Thermometer may be further detailed as follows:

Measurement time:

- i. Initial start-up time: 3-4 seconds
- ii. Measurement time: 2-3 seconds

Sincerely,

A handwritten signature in cursive script that reads "Alisa R. Salibra".

Alisa R. Salibra  
Global Category Manager, Thermometry  
Welch Allyn