

Welch Allyn, Inc
 4341 State Street Road
 P.O. Box 220
 Skaneateles Falls, NY 13153-0220
 USA
 Telephone: 315-685-4100



December 14, 2015

To: Welch Allyn US Distributors and Customers

Re: Pro4000 and Pro6000 Probe Cover Compatibility

From: Jamie Wood, US/C Product Manager

With the introduction of the new Braun Pro6000 Ear Thermometer, we have also launched a new Probe cover for the Pro6000 device. This probe cover is also backwards compatible to the Braun Pro4000 device when proper loading technique is used (see diagram below):

The Braun Pro4000 Probe covers are only compatible with the Braun Pro4000 device and will not work on the Pro6000.

PRO 4000 Probe Cover	PRO 6000 Probe Cover
Flat Flange	Raised Flange
Compatibility: Pro 4000 <u>only</u>	Compatibility: Pro6000 and Pro4000
Part #: 05075-005	Part #: 06000-005
<p>The image shows the packaging for the PRO 4000 Probe Cover. The box is white with black and blue text. It features the BRAUN logo and 'ThermoScan PC 200' branding. Below the box, a white probe cover with a flat flange is shown. The box text includes '200 x' and 'Use only Braun certified probe covers for accurate readings'.</p>	<p>The image shows the packaging for the PRO 6000 Probe Cover. The box is white with blue and black text. It features the BRAUN logo and 'ThermoScan PC 200' branding. Below the box, a white probe cover with a raised flange is shown. The box text includes '200 x' and 'Use only Braun certified probe covers for accurate readings'.</p>

Follow these loading instructions for proper loading of probe covers on Braun Pro4000:

BRAUN

PRO 4000
Proper Probe
Cover Installation

WelchAllyn

1. 2. 3. 4.

This product is manufactured by Kaz USA, Inc. under license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany. ©2014 Kaz USA, Inc.