

# SIMPLE, SECURE AND CONNECTED WORKFLOWS WITH BAXTER CARDIOLOGY

## THE BARDYDx **CAM** PATCH AND **CARDIO SERVER** ECG MANAGEMENT SYSTEM

We're helping you improve clinical workflows and reduce manual documentation for cardiac monitoring data by sending it straight to your EMR. See how the latest **Cardio Server** ECG Management System supports seamless connectivity with the **CAM** Patch.

### STANDARD LT ECG WORKFLOW

- 1 Order patch
- 2 Manually document patient demographics and apply patch
- 3 Patient wears patch then mails it back
- 4 Data is captured and a report is generated
- 5 Report is posted
- 6 Staff logs in and downloads report
- 7 Add physician notes to report
- 8 Upload to the EMR
- 9 Link the report to the patient's chart

### CONNECTED **CAM** PATCH WORKFLOW WITH **CARDIO SERVER**

- 1 Order **CAM** Patch
- 2 Patient wears **CAM** Patch for up to 14 days then the clinic directly uploads data to the cloud\*
- 3 Bardy IDTF captures data and generates report that is posted in **BDxConnect** Portal and automatically sent to **Cardio Server** with discrete data
- 4 Physician opens report in **Cardio Server** or **BDxConnect** Portal to review, edit and sign. Final report automatically sent to the EMR


\*Option to mail CAM Patch back to BardyDx is available

# CONNECTIVITY COUNTS

## HELP SAVE TIME BY STREAMLINING YOUR CLINICAL WORKFLOWS

- Reduce workflow steps previously spent manually entering data
- Automatically associate results with the EMR orders
- Free up your time for other patient priorities

## WORKFLOW DIAGRAM



Contact your **Cardio Server** representative to learn more about how the connected Baxter cardiology solutions can help you simplify your workflows, save time and improve data analysis.

**Rx only.** For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

### Baxter.com

Baxter International Inc.

Baxter, BardyDx, BDxConnect, CAM and Cardio Server are trademarks of Baxter International Inc. or its subsidiaries.

HL7 is a trademark of Health Level Seven, Inc.

US-FLC198-240006 (v1.0) 03/2024