

Centrella® Smart+ Bed

Meets the changing needs of healthcare facilities and delivers elevated care through optimized patient safety, enhanced patient satisfaction and advanced caregiver-focused technology.



Standard Features Available on Base Model

OPTIMIZED PATIENT SAFETY

Enhanced 3-Mode Bed Exit with Silence and Auto Reset

Motion Activated Nightlight

One-button Patient Stand Assist

Advanced SlideGuard® Technology

Foley Position Limit & Drainage Bag Holder

Head Angle <30° Limit Control

Trendelenburg/Reverse Trendelenburg

Dampened Emergency CPR

Integrated Patient Helper Bracket

ENHANCED PATIENT SATISFACTION

Device Storage Area

Integrated Urinal Holders

FlexAfoot™ Powered Bed Extension

ADVANCED CAREGIVER-FOCUSED TECHNOLOGY

Graphical Caregiver Interface (GCI)® Touchscreen (right side)

Brake Not Set Alert

One-button Chair

Integrated Scale (kg/lbs)

Battery Backup

Four Corner Brake and Steer

Four Corner IV Pole Sockets

Integrated Restraint Holders

Lockout Controls

Built-in Line and Cord Management

WiFi Module

Smart Bed Ready (SideCom® Required)



Hillrom™

One Scalable Centrella® Platform

BASE FRAME	PACKAGE 255 OPTION	PACKAGE 355 OPTION
Standard Features (listed on front page)	Advanced Patient Pendant	Centrella® max Powered Air Surface
Scalability	USB Charging	SafeView®+ System with IllumiGuide™ Siderail Handgrip
	Contact-Free, Continuous Monitoring Capable	Second Graphical Caregiver Interface (GCI)® Touchscreen
	Verbal Safety Prompts	

TECHNICAL SPECIFICATIONS

Overall frame width	40" or 43" (Comfort Edition) (101.6 cm or 109.22 cm)
Overall frame length with foot extension	86.7" - 98.7" (220 - 251 cm) (retracted - extended)
Bed height maximum	32.5" (83 cm)
Bed height minimum	14.5" (37 cm)
Head section inclination	0° to 65°
Thigh section inclination	0° to 25°
Foot section inclination	0° to 50°
Trendelenburg/Reverse Trendelenburg position (max)	15°
Caster size type	5" (13 cm) Central Locking System
Safe working load	650 lbs (295 kg) (includes patient weight, surface and accessories)
Patient weight range	70 - 500 lbs (32 - 227 kg)
Centrella® core Surface 36"W x 76" - 88"L x 6"H (91 cm x 193-223 cm x 15 cm)	
Centrella® pro Surface 36" or 40"W x 76" - 88"L x 7"H (91.4 or 101.6 cm x 193-223 cm x 17.78 cm)	
Centrella® pro+ Surface 36" or 40"W x 76" - 88"L x 9"H (91.4 or 101.6 cm x 193-223 cm x 23 cm)	
Centrella® max Surface 36" or 40"W x 76" - 88"L x 8"H (91.4 or 101.6 cm x 193-223 cm x 20.32 cm)	

IEC 60601-2-52 Compliance

hillrom.com

A LA CARTE OPTIONS

OPTIMIZED PATIENT SAFETY FEATURES	Centrella® max Powered Air Surface (includes Blower, Comfort Controls, Turn Assist and Focused Advanced Microclimate® Management Therapy)
	Centrella® pro+ surface (includes integrated blower, and Focused Advanced Microclimate® Management Therapy)
	Centrella® core Surface
	PATIENT SAFETY APPLICATIONS
	NaviCare® Patient Safety
	Contact-Free, Continuous Monitoring
ENHANCED PATIENT SATISFACTION FEATURES	WatchCare® Incontinence Management
	Comfort Edition 40" Width
	Patient Pendant (up to 2)
	Experience POD™ Device
	Active Nurse Call Button
ADVANCED CAREGIVER-FOCUSED TECHNOLOGY FEATURES	Universal TV and Lighting Controls
	SafeView®+ System (with IllumiGuide™ Siderail Handgrip)
	5th Wheel Steering
	IntelliDrive® Powered Transport
	Second Graphical Caregiver Interface (GCI)® Touchscreen
	Accessory Outlet
	Additional Line Management
Obstacle Detect® System	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 203601 rev 7 28-JUN-2021 ENG - US