



Compella™ Bariatric Bed

Item Description: Advancing Bariatric Care

Model Number: P7800

JSN Number: _____

Warranty: 1 Year From Delivery Date

Installation: _____

Proposal Number: _____

Facility: _____

City, State: _____

Estimated Delivery: _____

Delivery Cost: _____

Taxes: _____

Product Standard Features

EXCLUSIVE Powered Width Expansion from 40" – 50"

FlexAfoot™ Bed Length Adjustment Feature

Single-Mode Bed Exit

Digital HOB Angle Indicators and Alert

Integrated Caregiver and Patient Controls

In-Bed Scale

Battery Backup

Product Accessory Options

IntelliDrive® XL Powered Transport
(Standard on Rental)

3-Mode Bed Exit (Standard on Rental)

Patient-Helper Trapeze

Patient Pendant

Smart Bed Ready – SideCom® System Ready
(Standard on Rental)

Nurse Call – Capital Purchase Only

Low Air Loss Surface with Turn Assist

Low Air Loss Surface with CLRT


Compella™ Bariatric Bed

Item Description: Advancing Bariatric Care

Model Number: P7800

JSN Number: N/A

Proposal Number: _____



Technical Specifications

Bed	
Low position	Sleep deck to floor 18.5" (47 cm)
High position	Sleep deck to floor 26.8" (68 cm)
Transport position	Sleep deck to floor 17" (43 cm)
Minimum overall width (width extenders retracted)	43" (109 cm)
Minimum overall length (foot extension retracted)	91" (231 cm)
Maximum Patient weight capacity	1,000 lbs (453.6 kgs)
Electrical Information	120/230 V 12/6A 50/60 Hz

Bed Angles

Head section	0-50°
Knee section	0-30°
Max Tilt/Reverse Tilt	0-9°

Surface

Surface width	40" to 50" (102 cm to 127 cm)
Surface length - inflated	80" to 88" (203 cm to 224 cm)
Surface Height - inflated (foam)	8" (20 cm)
Surface Height - inflated (low air loss)	8.5" (21.5 cm)

Caster & Braking System

Caster size	6" (15 cm)
Four-wheel braking	Yes
Brake Not Set Alert	Yes
Brakes lock both rolling & swiveling	Yes

Scale System

Accuracy	1% of patient weight (with sleep deck flat)
Weigh in any position	Yes
Scale weight capacity	1,000 lbs (453.6 kgs)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
207196 rev 2 21-MAY-2018 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom