



SIMPLIFY CLINICAL WORKFLOWS AND SECURE DEVICE ACCESS.

FAST, SECURE ACCESS TO YOUR **WELCH ALLYN** VITAL SIGNS DEVICES WITH IMPRIVATA SINGLE SIGN-ON

Balancing the benefits of vital signs monitor connectivity to your EMR with the need for strong user authentication doesn't need to mean sacrificing workflow simplicity for security. **Welch Allyn** vital signs monitors with integration with Imprivata bring fast, secure badge-tap access to the bedside.



EXTENDING BADGE ACCESS THAT CLINICIANS ARE ALREADY FAMILIAR WITH TO YOUR CONNEX DEVICES CAN HELP YOU PROTECT PATIENT DATA WHILE SAVING VALUABLE TIME.



STREAMLINE CLINICAL WORKFLOWS

Securely log into devices without manually entering your username and password so you can spend more time focused on patient care.** With the **Welch Allyn** dual function barcode scanner & RFID reader, clinicians have the ease of using one device to sign on to the vital signs monitor and scan their patient's wristband for a completely connected workflow.




CONTROL DEVICE ACCESS

A simple badge scan or tap provides a quick authentication at the bedside granting device access only to authorized users.

YOUR CONNECTED CLINICAL WORKFLOW WITH SINGLE SIGN-ON.

The process is simple. Start by opening a session with Imprivata — either at the EMR workstation, on the vital signs monitor or elsewhere. This authentication is then extended to your **Connex** vital signs device. Access is granted in a single step, allowing you to take patient vitals and send data to your EMR.



CONTACT YOUR BAXTER REPRESENTATIVE OR VISIT HILLROM.COM TO LEARN MORE.

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator’s Manual or Instructions for Use.

* RFID support on **Connex** Spot Monitor software version 1.32 or higher and **Connex** Vital Signs Monitor 2.42 or higher.
You must first sign in with Imprivata in accordance with your facility’s policies. Then you can simply tap or scan your badge on select **Connex devices to log into the device.
*** Consult with your Baxter representative to validate usability with your specific RFID card type.

Baxter.com | Baxter International Inc.

Baxter, Connex and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.
Imprivata is a trademark of Imprivata, Inc.
Masimo and SET are trademarks of Masimo Corporation.