

MOUNTING AND MOBILITY OPTIONS

VITAL SIGNS DATA WITHIN YOUR
REACH — ACROSS CARE SETTINGS
AND WORKFLOWS.



Mobile stands support clinical workflows that require a shared device across multiple rooms.

SPECIFICATION	STANDARD BASKET STAND	WORK SURFACE AND ACCESSORY MANAGEMENT STAND	ACCESSORY POWER MANAGEMENT STAND
	Simple portability with a small footprint	Organized storage and cable management, plus a compact mobile work surface footprint	Mobile work surface and integrated battery to power through a busy day
Part Number*	7000-MS3	7000-MWS	7000-APM
Height	49" (124.46 cm) total 38" (96.52 cm) to top of basket handle	47" (119.38 cm) total 39.25" (99.7 cm) to work surface	47" (119.38 cm) total 41" (104.14 cm) to work surface
Base Diameter	23" (58.42 cm)	23" (58.42 cm)	23" (58.42 cm)
Work Surface Area	N/A	11.62" (29.51 cm) wide 11.25" (28.56 cm) deep	11.5" (29.21 cm) wide 10" (25.4 cm) deep
Accessory Basket	8.0" (20.32 cm) wide 5.0" (12.70 cm) deep 7.0" (17.78 cm) high	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" (19.05 cm) high 4 accessory bays	11.75" (29.85 cm) wide 8.25" (20.96 cm) deep 7.5" (19.05 cm) high 4 accessory bays
Weight	13.8 kg / 30.5 lb (fully loaded with device and accessories)	15.4 kg / 34 lb (fully loaded with device and accessories)	18.4 kg / 40 lb (fully loaded with device and accessories)
Power Accessory Functionality	N/A	N/A	9 cell Lithium Ion Battery (Field replaceable, integrated power supply, integrated work light, battery life indicator display, two USB ports)

*Connex Spot Monitor sold separately.



DESKTOP STAND
7000-DST



WALL MOUNT WITH BASKET
7000-GCX



INTEGRATED WALL SYSTEM
777-PM2WCX-US**
777-PM3WCX-US**

Wall mount options are ideal for care settings that need a dedicated device in every exam room, while desktop stands can be moved as needed.

SPECIFICATION	DESKTOP STAND	WALL MOUNT WITH BASKET	INTEGRATED WALL SYSTEM
	Compact and portable tabletop mount suitable for placement on any work surface	Easy-to-access storage and cable management	Complete set of physical exam tools in one integrated package
Part Number(s)*	7000-DST	7000-GCX	777-PM2WCX-US, 777-PM3WCX-US
Dimensions	Height 9.5" (24.13 cm) Width 7.25" (18.42 cm) Depth 4.75" (12.07 cm)	Wall Channel: 19" (48.26 cm) Pivot Arm: 8" (20.32 cm) Wall Mount Utility Basket: 12.75" (32.39 cm) wide x 7.5" (19.05 cm) deep x 5.9" (14.99 cm) high	Height: 12" (30.48 cm) Length: 34" (86.36 cm) Depth: 2" (5.08 cm) <i>Depth dimension is for wall board only</i>

Contact your Baxter Sales Representative or visit hillrom.com/CSM to learn more.

Baxter.com

Baxter International Inc.

*Connex Spot Monitor sold separately.

777-PM2WCX-US model includes mounting bracket for **Connex Spot Monitor or **Spot Vital Signs** 4400 Device with **PanOptic** Plus LED Ophthalmoscope, **MacroView** Plus LED Otoscope for iExaminer, Ear Specula Dispenser.

***777-PM2WCX-US model includes mounting bracket for **Connex** Spot Monitor or **Spot Vital Signs** 4400 Device with **PanOptic** Plus LED Ophthalmoscope, **MacroView** Plus LED Otoscope for iExaminer, Ear Specula Dispenser.

Baxter, Connex, Hillrom, KleenSpec, MacroView, PanOptic, Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC45-230027 V1 07/23