

Part Number Configurations

| Connex® Vital Signs Monitor | | | | | |
|-----------------------------|------------------|-------------------------------|----------------------------------|--------------------|------------------|
| 6 | X | X | X | X | X |
| Model | Radio | Pulse Oximetry | Respiration | Thermometry | Printer |
| 6 | 7 = Not Included | N = Nellcor® SpO ₂ | C = CO ₂ | T = SureTemp® Plus | P = Internal |
| | 8 = Included | M = Masimo® SpO ₂ | R = RRa | E = Braun 6000 | X = Not Included |
| | | | S = EarlySense® F = ES Safety | D = Both | |
| | | | X = Not Included | X = None | |

| Connex® Integrated Wall System | | | | | | |
|--------------------------------|------------------|------------------------------|-------------------|---------------------------------|-----------------|--|
| 8 | X | X | X | X | X | X |
| Model | Radio | Pulse Oximetry | Oral Temp | Storage | Ear Temp | Heads |
| 8 | 4 = Not Included | N = Nellcor SpO ₂ | T = SureTemp Plus | V = SCS (Standard Cuff Storage) | E = Braun 6000 | C = MacroView™ 23820 + Coaxial Ophthalmoscope 11720 |
| | 5 = Included | M = Masimo SpO ₂ | X = No Oral Temp | | X = No Ear Temp | P = MacroView 23820 + PanOptic™ Ophthalmoscope 11810 |
| | | X = No SpO ₂ | | | | X = No Heads |

For Connex devices, power cords are indicated as “-X” at the end of the part number as a 1-digit alphanumeric code (e.g., 84MTVEC-B or 67MXTP-B).

| Power Cord Key | |
|----------------|-----------------------|
| B | North America |
| 2 | Europe |
| 4 | United Kingdom |
| 6 | Australia/New Zealand |
| 7 | South Africa |

| Model | 67 | 68 | 84 | 85 |
|---|----|----|----|----|
| Wireless Connectivity | | | | |
| Upgradeable | ✓ | | ✓ | |
| Integrated Radio | | ✓ | | ✓ |
| Wired Connectivity | | | | |
| USB | ✓ | ✓ | ✓ | ✓ |
| Ethernet | ✓ | ✓ | ✓ | ✓ |
| Nurse Call | ✓ | ✓ | ✓ | ✓ |
| Number of Accessories Capable of Running | 4 | 4 | 4 | 4 |