



Baxter

ECG SOLUTIONS FOR LONG-TERM CARE



CP 150 RESTING ECG

Simple to use with a 7" touchscreen, full-size printout and the ability to save and send PDF files from the device

DIAGNOSTIC CARDIOLOGY SUITE ECG

Designed for EMR-connected workflows, the Diagnostic Cardiology Suite includes ECG and spirometry capabilities to help you gather key data for diagnosis

ELI 150c RESTING ECG

Compact and lightweight, the **ELI 150c** Electrocardiograph provides comprehensive functionality in an ultra-portable device



SIMPLE

Reduce manual data entry and leverage simple-to-use devices to deliver accurate results



SECURE

Control user logins, assign user roles and access audit trails to secure patient data with Diagnostic Cardiology Suite



CONNECTED

Capture patient data and send immediate results to your EMR automatically

WITH BAXTER ECG SOLUTIONS, YOU CAN:

- Perform ECG exams real-time and in-house, while symptoms are occurring, so conditions can be quickly identified and addressed
- Reduce money spent on hospital transport and outsourced ECG costs

When you send residents to a hospital for ECG procedures, you may be sending more out the door than you think. Transportation costs, logistics management, resident stress and many other factors can adversely impact the level of care you provide — and your bottom line.

THAT'S WHY BAXTER ECG SOLUTIONS ARE DESIGNED TO HELP YOU KEEP THESE VALUABLE PROCEDURES IN-HOUSE.

CP 150

CP150A-1ENB	CP 150 Interpretive Electrocardiograph
CP150AS-1ENB	CP 150 Interpretive Electrocardiograph with Spirometry

DIAGNOSTIC CARDIOLOGY SUITE

CC-RXX-WAXX	DCS Resting ECG with WAM
CC-RXX-WADX	DCS Resting ECG w WAM and DICOM
CC-RXX-AAXX	DCS Resting ECG w AM12
CC-RXX-AADX	DCS Resting ECG w AM12 & DICOM

ELI 150C

ELI150c-AAA-AACAD	ELI 150c, AM12 , AHA banana end, interpretation, 40 patient storage, LAN/USB and DICOM
ELI150c-BAA-AACAD	ELI 150c, WAM, AHA banana end, interpretation, 40 patient storage, LAN/USB and DICOM

Contact your Baxter representative or visit hillrom.com to learn more.

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Baxter.com

Baxter International Inc.

Baxter, AM12, CP 150 and ELI are trademarks of Baxter International, Inc. or its subsidiaries.
US-FLC50-220015 (v1.0) 01/2025