



ECG MONITORING ELECTRODES

ELECTRODE DETAILS

Electrode Style	Adult, Foam, Single-Use Electrode
Reference	108070
Adhesive Type	Wet Gel
Quantity	10 Electrodes/Pouch (300 Electrodes/Case)
Application	General Patient Monitoring (for 72 Hours)
Shelf Life	36 Months (in Accordance with Storage Guidelines)

TECHNICAL SPECIFICATIONS

Shape	Teardrop
Size (Excluding Grip)	1.5 in (38 mm)
Sensor Diameter (Eyelet)	0.35 in (9.0 mm)
Substrate Thickness (Adapter Excluded)	0.06 in (1.5 mm)
Total Product Surface Area	1.7763 in ² (45.117 mm ²)
Gel Area	0.361 in ² (9.1694 mm ²)
Adhesive Area	1.4153 in ² (35.9474 mm ²)
Integrated Lead Wire (Length/Color)	N/A

MATERIALS

Substrate Material	Polyethylene Foam
Adhesive	Medical Grade Adhesive
Gel Type	Conductive Hydrogel
Foam (Sponge) Material	Polyurethane Reticulated Foam
Release Liner	PET
Sensor Polymer	Ag/AgCl Plated ABS
Adaptor/Connector (Stud)	Stainless Steel
Integrated Lead Wire Jacketing	N/A
Integrated Lead Wire Cord	N/A

BIOCOMPATIBILITY

ISO 10993	Passed
Not Made with Natural Rubber Latex	Yes

ELECTRICAL PERFORMANCE (ANSI/AAMI EC 12)

ACZ impedance @10 Hz (Before Defib Simulation)	58 - 93 Ohm
DC Offset Voltage (Before Defib Simulation)	0.2 - 0.5 mV
SDR Slope @ 30 Sec Int. (Remaining Potential After Defib)	0.2 - 0.3 mV/sec
ACZ impedance repeat (After Defib Simulation)	48 Ohm
Combined Offset Instability and Inner Noise (COIIN)	18 - 20 μ V
Bias Current Tolerance (DC Offset Voltage After DC Loading)	< 1.0 mV

FEATURES

MR Conditional	No
X-Ray Translucence	No
Integrated Abrader	No
Repositionability	No

PACKAGING

Product Packaging Material	Paper/PE/Foil/PE
Resealable Pouch	No
Product Packaging Size (L x W)	11.42 in x 3.74 in (29.0 cm x 9.5 cm)
Department Packaging Box (L x W x H)	N/A
Transport Packaging Case (L x W x H)	12.6 in x 8.27 in x 6.3 in (32.0 cm x 21.0 cm x 16.0 cm)

ENVIRONMENTAL

Halogenated Hydrocarbon Content (e.g. PVC)	None
Phthalate Derivatives Content (e.g. DEHP)	None
RoHS Compliant	Yes
REACH Compliant	Yes

VISIT [HILLROM.COM](https://www.hillrom.com) OR CONTACT YOUR BAXTER SALES REPRESENTATIVE TO LEARN MORE.

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

[Baxter.com](https://www.baxter.com)

Baxter International Inc.

One Baxter Parkway / Deerfield, Illinois 60015

Baxter is a trademark of Baxter International Inc.

US-FLC40-220007 (v1.0) 05/2024