

**Baxter**

# TAKE YOUR ED ECG WORKFLOW TO NEW HEIGHTS

## AUTO-RECONCILIATION FROM BAXTER CARDIOLOGY



Get STAT ECGs to your physicians for review without delay — with auto-reconciliation from Baxter Cardio Server and **ELI** ECG devices.

Is your team spending time manually reconciling studies to orders for tests performed before the order has been placed?

Automated order reconciliation from Baxter Cardiology's end-to-end solution greatly reduces manual workflow steps used in the Emergency Department and beyond, saving clinicians valuable time.



## AUTO-RECONCILIATION WITH BAXTER CARDIOLOGY:



Spend less time manually reconciling ECG studies.



Automate reconciliation even if an order is entered after Cardio Server receives the exam.



Improve billing efficiencies with simplified STAT/order reconciliation process.

# SIMPLIFY YOUR ED ECG WORKFLOW

1

Enter correct mandatory patient demographic information.

2

Acquire ECG on Baxter **ELI 380** or **ELI 280**.

3

Send ECG to Cardio Server.

4

Cardio Server will auto-reconcile to existing and/or future orders.



## TAKE ADVANTAGE OF AUTO-RECONCILIATION TODAY

Cardio Server uses mandatory fields below to auto-reconcile exams that are not already associated with an order. Simply configure your **ELI 380** mandatory fields to match Cardio Server.

1

Patient ID

2

Date of Exam

3

First and Last Name

4

Date of Birth  
*\*Optional*

Contact your Baxter representative for assistance.

[Baxter.com](http://Baxter.com) | Baxter International Inc.

Baxter and ELI are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC64-230018 V1 08/2023