

Baxter

TAKE YOUR ECG WORKFLOW TO NEW HEIGHTS

BAXTER CARDIOLOGY ENHANCEMENTS WITH AUTO-RECONCILIATION AND PATIENT DEMOGRAPHIC QUERY



Get STAT ECGs to your physicians for review without delay — with auto-reconciliation from Baxter Cardio Server and ELI ECG devices.

Is your team spending time manually reconciling studies to orders for tests performed before the order has been placed?

Automated order reconciliation from Baxter Cardiology's end-to-end solution greatly reduces manual workflow steps used in the Emergency Department and beyond, saving clinicians valuable time.



AUTO-RECONCILIATION WITH BAXTER CARDIOLOGY:



Spend less time manually reconciling ECG studies.



Automate reconciliation even if an order is entered after Cardio Server receives the exam.



Improve billing efficiencies with simplified STAT/order reconciliation process.

SIMPLIFY YOUR ECG WORKFLOW

1

Scan patient barcode or enter mandatory patient demographic information.

2

Acquire ECG on Baxter **ELI 380** or **ELI 280**.

3

Send ECG to Cardio Server.

4

Cardio Server will auto-reconcile to existing and/or future orders.



TAKE ADVANTAGE OF AUTO-RECONCILIATION AND PATIENT DEMOGRAPHIC QUERY TODAY

Automatically match exams to orders in Cardio Server from Baxter **ELI 280** and **380** carts.

UPG380-A-PDQ-A:

Patient Demographic Query upgrade for **ELI 380**

UPG280-A-PDQ-A:

Patient Demographic Query upgrade for **ELI 280**

Contact your Baxter representative to learn more.

Baxter.com | Baxter International Inc.

Baxter and ELI are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC64-230018 (v2.0) 08/2024