



Hillrom™

Excel Care® Bariatric Bed
Rental offering

THE FLEXIBILITY YOU NEED

From lower acuity care to the ICU, hospitals are adjusting to accommodate an increasing population of bariatric patients.

Because of increased lengths of stay, injuries related to mobility issues and higher readmission rates can be costly. Therefore, it's critical to have access to the right equipment to provide appropriate bariatric care for a range of conditions.^{1,2}



The Excel Care® Bariatric Bed helps you deliver exceptional care with:

- Flexible patient support
- Simplified care
- Support patient independence

FLEXIBLE PATIENT SUPPORT

The Excel Care® Bariatric Bed provides flexible support to help you manage safety issues and patient outcomes.

- Pressure redistribution surface helps treat and prevent pressure injuries.
- Continuous air flow keeps patient skin cool and dry.
- Turn Assist helps make it easier to perform linen changes and wound checks.
- Air surface remains inflated during transport.
- Continuous Lateral Rotation Therapy (CLRT) supports protocols for helping prevent ventilator associated complications.
- Low bed height and seat deflation help make transition to and from bed easy and safe.





SIMPLIFIED CARE

From repositioning patients to performing other routine care, the **Excel Care® Bed** helps simplify the way patients are handled with easy-to-use integrated features.

- Expandable width frame and surface adjust from 40" to 50" for patient comfort.
- Adjustable bed length helps minimize patient migration and the need for patient repositioning.
- Automated features such as an integrated scale to help offer multiple options to manage your patient — in the comfort of their own bed.



SUPPORT PATIENT INDEPENDENCE

The **Excel Care® Bariatric Bed** helps safeguard patient dignity with accessible controls that help support patient independence.

- With multiple grip points and single post patient helper trapeze, patients can assist the caregiver during repositioning or make themselves more comfortable.
- Easy-to-reach pendant offers easy access and a convenient way for patients to comfortably adjust themselves in bed.
- 1,000 lb patient weight capacity

TECHNICAL SPECIFICATIONS

Bed	
Low Position - Sleep Deck to Floor	17" (43.2 cm)
High Position - Sleep Deck to Floor	29" (73.7 cm)
Maximum Overall Width	44"-54" (111.8-137.2 cm)
Maximum Overall Length	91"-100" (254 cm)
Patient Weight Capacity	1,000 lbs (453.6 kgs)
Bed Angles	
Head Section	0°-48°
Knee Section	0°-20°
Max Trend/Reverse Trend	12°/10°
Sleep Surfaces & Air Supply Unit	
Foam	80"-88"L x 40"-50"W x 8.0"H (223.5 x 101.6-127 x 20.3 cm)
Low Air Loss with Turn Assist	80"-88"L x 40"-50"W x 8.5"H (223.5 x 101.6-127 x 21.6 cm)
Low Air Loss with CLRT	80"-88"L x 40"-50"W x 8.5"H (223.5 x 101.6-127 x 21.6 cm)
Air Supply Unit Dimensions	12.1"H x 14.5"W x 6.0"D (30.7 x 36.8 x 15.2 cm)
Caster & Braking System	
Caster Diameter	6" (15.2 cm)
Four Wheel Braking	Yes
Steer and Lock Both Rolling and Swivelling	Yes
Scale System	
Accuracy	1% of patient weight (with sleep deck flat)
Weigh in Any Position (except Trend/Rev Trend)	Yes
Scale Weight Capacity	1,000 lbs (453.6 kgs)
Safety Listing	UL 60601-1
Power Requirements	120v, 60 cycle AC
Fuse, Maximum Current	5A, 6A
Fire Safety	
Meets the Requirements of Foam Mattress also Meets	16CFR1632 CAL129, CAL603, NFPA101, BFDIX-11

STANDARD FEATURES - RENTAL

Expandable Width and Length Frame and Surface

Patient Pendant

Battery Backup

In-bed Scale

Head and Trend Angle Indicators

OPTIONAL FEATURES

Patient-helper Trapeze

Low Air Loss Surface with Turn Assist

Low Air Loss Surface with CLRT



Hillrom™

For more information, please contact your Hillrom sales representative at 800-638-2546 or order online at rental.hillrom.com.

hillrom.com

¹ Rice S. Retrofitting hospitals for obese patients. Modern Healthcare. 2014 Feb. modernhealthcare.com/article/20140208/MAGAZINE/302089980.
² McGinley L. Best Practices for Safe Handling of the Morbidly Obese Patient. Bariatric Nursing and Surgical Patient Care. Volume 3, Number 4, 2008.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.