

Excel Care[®] ES Bariatric Bed

Skills Checklist



Staff Signature: _____

Trainer's Signature: _____

Date: _____

Perform the following.

COMPLETED

Bed Operations

1. Demonstrate how to prepare bed for transport.
2. Demonstrate how to zero the scale.
3. Demonstrate how to weigh a patient.
4. Demonstrate how to adjust a patient's weight.
5. Demonstrate how to switch between lbs and kgs.
6. Demonstrate how to activate battery backup.
7. Demonstrate how to expand and contract sleep deck frame.
8. Demonstrate use of caregiver controls and lockout controls.
9. Demonstrate how to activate CPR for both frame and surface.
10. Demonstrate cardiac chair position.
11. Demonstrate how to adjust bed length for both frame and surface.
12. Locate drainage bag holders.
13. Demonstrate the brake and steer mechanism.
14. Demonstrate Trendelenburg and Reverse Trendelenburg.
15. Locate head-of-bed angle indicator.
16. Demonstrate operation of siderails.
17. Demonstrate patient helper trapeze movement.

Standard Surface Functions

18. Demonstrate how to navigate air supply controls.
19. Locate on/off switch, power cord attachment and hose attachments.
20. Demonstrate how to deflate/inflate side bolsters
21. Demonstrate how to deflate/inflate foot bolsters

Excel Care® ES Bariatric Bed

Skills Checklist

Perform the following.

COMPLETED

Optional Low Air Loss Surface Functions

- | | COMPLETED |
|--|-----------|
| 22. Demonstrate patient setup. | |
| 23. Demonstrate how to adjust comfort settings. | |
| 24. Demonstrate fowler boost adjustments. | |
| 25. Demonstrate seat deflate mode. | |
| 26. Demonstrate max inflate mode. | |
| 27. Demonstrate Turn Assist. | |
| 28. Demonstrate Continuous Lateral Rotation Therapy. | |

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
183340 rev 3 25-NOV-2015 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom