

Explore Welch Allyn Physical Exam Sets and Systems featuring our new, LED wide-view PanOptic Plus Ophthalmoscopes and MacroView Plus Otoscopes

## WALL SYSTEMS AND TRANSFORMERS

|                                |               |               |                         | OPHTHALMOSCOPES                              |                   |              | OTOSCOPES & ACCESSORIES                  |                           |   |                       | VITAL SIGNS     |                              |  |                    |
|--------------------------------|---------------|---------------|-------------------------|--|-------------------|--------------|--|---------------------------|---|-----------------------|-----------------|------------------------------|--|--------------------|
|                                |               | Wall<br>Board | 777 Wall<br>Transformer | PanOptic Plus<br>with Quick Eye<br>Alignment | PanOptic<br>Basic | Coaxial      | MacroView<br>Plus, Digitally-<br>Capable | <b>MacroView</b><br>Basic | iExaminer<br>SmartBracket for<br>MacroView Plus | Ear Spec<br>Dispenser | Wall<br>Aneroid | SureTemp Plus<br>Thermometer | Connex Spot<br>Monitor/4400<br>Vital Signs Device* | 5-Year<br>Warranty |
|                                |               |               |                         | 118-3-US                                     | 118-2-US          | 11720-L      | 238-3                                    | 238-2                     | 11840-IEX-US                                    |                       |                 |                              | Vital digits bevice                                |                    |
|                                | 777-PM3WCX-US | ✓             | ✓                       | ✓  |                   |              | ✓  |                           | ✓   | ✓                     |                 |                              | ✓  | ✓                  |
|                                | 777-PM3WAS-US | ✓             | ✓                       | ✓  |                   |              | ✓  |                           | ✓   | ✓                     | $\checkmark$    | ✓                            |  | ✓                  |
|                                | 777-PM3WAX-US | ✓             | ✓                       | ✓  |                   |              | ✓  |                           | ✓   | ✓                     | $\checkmark$    |                              |  | ✓                  |
| ed<br>st                       | 777-PM3WXS-US | ✓             | ✓                       | ✓  |                   |              | ✓  |                           | ✓   | ✓                     |                 | ✓                            |  | ✓                  |
| 777 Integrated<br>Wall Systems | 777-PM3WXX-US | ✓             | $\checkmark$            | ✓  |                   |              | ✓  |                           | ✓   | ✓                     |                 |                              |  | ✓                  |
|                                | 777-PM2WCX-US | ✓             | $\checkmark$            |  | $\checkmark$      |              |  | ✓                         |   | ✓                     |                 |                              | ✓  |                    |
| \                              | 777-PM2WAS-US | ✓             | $\checkmark$            |  | $\checkmark$      |              |  | ✓                         |   | ✓                     | $\checkmark$    | ✓                            |  |                    |
|                                | 777-PM2WAX-US | ✓             | $\checkmark$            |  | $\checkmark$      |              |  | ✓                         |   | ✓                     | $\checkmark$    |                              |  |                    |
|                                | 777-PM2WXS-US | ✓             | $\checkmark$            |  | $\checkmark$      |              |  | ✓                         |   | ✓                     |                 | ✓                            |  |                    |
| Transformers                   | 777-PM2WXX-US | ✓             | $\checkmark$            |  | $\checkmark$      |              |  | ✓                         |   | ✓                     |                 |                              |  |                    |
|                                | 777-SM2WXX    | ✓             | $\checkmark$            |  |                   | $\checkmark$ |  | ✓                         |   | ✓                     |                 |                              |  |                    |
|                                | 777-SM2WAS    | ✓             | $\checkmark$            |  |                   | $\checkmark$ |  | ✓                         |   | ✓                     | $\checkmark$    | ✓                            |  |                    |
|                                | 777-SM2WAX    | ✓             | $\checkmark$            |  |                   | $\checkmark$ |  | ✓                         |   | ✓                     | $\checkmark$    |                              |  |                    |
|                                | 777-PM3XXX-US |               | ✓                       | ✓  |                   |              | ✓  |                           | ✓   |                       |                 |                              |  | ✓                  |
|                                | 777-PM2XXX-US |               | ✓                       |  | $\checkmark$      |              |  | ✓                         |   |                       |                 |                              |  |                    |
|                                | 777-SM2XXX    |               | ✓                       |  |                   | <b>√</b>     |  | <b>√</b>                  |   |                       |                 |                              |  |                    |

 $^*$ Choosing these models will ensure the necessary mounting bracket ships with your wall board.

## **DIAGNOSTIC SETS**

|                | OPHTHALMOSCOPES   |          |   | OTOSCOPES & ACCESSORIES                  |       |   |                          | POWER HANDLES (sets include one power handle) |                   |              |              |                      |   |
|----------------|---|----------|---|--|-------|---|--------------------------|---|-------------------|--------------|--------------|----------------------|---|
|                | PanOptic Plus with Quick Eye Alignment PanOptic Basic Coaxial |          | <b>MacroView</b> Plus, Digitally- Capable | ıs, Digitally-   Macroview   Traditional |       | iExaminer<br>SmartBracket for<br>MacroView Plus | Lithium-lon<br>Plus with | Lithium-lon                                   | Nickel<br>Cadmium | Hard Case    | Soft Case    | 5-Year<br>Warranty** |   |
|                | 118-3-US  | 118-2-US | 11720-L                                   | 238-3                                    | 238-2 | 250-2   | 11840-IEX-US             | built-in USB                                  |                   |              |              |                      |   |
| 71-PM3LXE-US   | ✓   |          |   | ✓  |       |   | ✓                        | ✓   |                   |              | ✓            |                      | ✓ |
| 71-XM3LXE      |   |          |   | ✓  |       |   | ✓                        | ✓   |                   |              | ✓            |                      | ✓ |
| 71-PX3LXERT-US | ✓   |          | $\checkmark$                              |  |       |   |                          | ✓   |                   |              | ✓            |                      | ✓ |
| 71-PX3LXE-US   | ✓   |          |   |  |       |   |                          | ✓   |                   |              | ✓            |                      | ✓ |
| 71-PM3LXES-US  | ✓   |          |   | ✓  |       |   | ✓                        | ✓   |                   |              |              | ✓                    | ✓ |
| 71-PM2LXU-US   |   | ✓        |   |  | ✓     |   |                          |   | $\checkmark$      |              | ✓            |                      |   |
| 71-SM2LXU      |   |          | $\checkmark$                              |  | ✓     |   |                          |   | $\checkmark$      |              | ✓            |                      |   |
| 71-SS2CXX      |   |          | $\checkmark$                              |  |       | ✓   |                          |   |                   | $\checkmark$ | ✓            |                      |   |
| 71-SM2CXX      |   |          | ✓   |  | ✓     |   |                          |   |                   | $\checkmark$ | $\checkmark$ |                      |   |

Welch Allyn Green Series 777 Integrated Wall Systems and Transformers can help you stay efficient during your examinations. Select from wide-view and standard ophthalmoscopes and otoscopes, digital or manual blood pressure devices and integrated or stand-alone oral/axillary/rectal thermometry with a Welch Allyn SureTemp Plus Thermometer.



Welch Allyn Diagnostic Sets are the perfect solution to keep your exam tools secure for whenever you need them. Try our wide-view optics and fast charge power supplies built into our Plus models.



Welch Allyn Connex Spot Monitor and/or Spot Vital Signs Device must be ordered separately.

\*\*5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

## **DESK SETS**

All Welch Allyn Desk Sets shown below ship with a Universal Desk Charger and power supply.

Welch Allyn Desk Sets are compact and perfect for smaller offices to support eye and ear examinations. Try our wide-view optics our Plus models.



#71-PM3LDE-US

|              | ОРНТН  | ALMOSCO | OPES         | отоѕсс  | PES & ACC    | CESSORIES  | POWER HANDLES (sets include two of the same handle) |                   |                      |  |
|--------------|--|---------|--------------|---|--------------|--|---|-------------------|----------------------|--|
|              | PanOptic Plus with Quick Eye Alignment  PanOptic Basic Coaxial |         | Coaxial      | MacroView Plus, Digitally- Capable  MacroView Basic |              | iExaminer<br>SmartBracket<br>for MacroView<br>Plus | Lithium-lon   | Nickel<br>Cadmium | 5-Year<br>Warranty** |  |
|              | 118-3-US 118-2-US  |         | 11720-L      | 238-3   | 238-2        | 11840-IEX-US                                       |   |                   |                      |  |
| 71-PM3LDE-US | ✓  |         |              | $\checkmark$  |              | ✓  | ✓   |                   | ✓                    |  |
| 71-PM2LDX-US |  | ✓       |              |   | $\checkmark$ |  | $\checkmark$  |                   |                      |  |
| 71-SM2LDX    |  |         | $\checkmark$ |   | $\checkmark$ |  | $\checkmark$  |                   |                      |  |
| 71-SM2NDX    |  |         | $\checkmark$ |   | $\checkmark$ |  |   | ✓                 |                      |  |
| 71-XXXLDE    |  |         |              |   |              |  | ✓   |                   |                      |  |

<sup>\*\*5-</sup>year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

## Welch Allyn **Ophthalmoscopes**







Welch Allyn

Otoscopes

DuraShock

image capture

Digitally-capable for

| Total Maria |
|-------------|
| 250.2       |





| Model                                 | Traditional | PanOptic<br>Basic | PanOptic<br>Plus |
|---------------------------------------|-------------|-------------------|------------------|
| Warranty                              | 1 year      | 1 year            | 5 year           |
| Larger field of view                  | _           | 20X               | 20X              |
| Designed with<br><b>DuraShock</b>     | _           | Yes               | Yes              |
| <b>Quick Eye</b> alignment technology | _           | _                 | Yes              |

|                      | 230-2       | 230-2              | 230-3                    |
|----------------------|-------------|--------------------|--------------------------|
| Model                | Traditional | MacroView<br>Basic | <b>MacroView</b><br>Plus |
| Warranty             | 1 year      | 1 year             | 5 year                   |
| Larger field of view | _           | 3X                 | 3X                       |
| Designed with        | _           | Yes                | Yes                      |

Connect with your Baxter representative today to discuss the best Welch Allyn Physical Exam & Diagnostic solution for your clinical needs.



Hillrom...

hillrom.com

Yes