







Hillrom™

HILLROM
1069 STATE ROUTE 46 EAST
BATESVILLE, IN 47006
PHONE: (812) 934-7777
FAX: (812) 934-8189
www.hillrom.com



PERSPECTIVE

TOP VIEW



LEGEND:

SIDE VIEW

- ① MAIN UNIT: 13.5 KG
- ② CONFERENCE UNIT: 8 KG
- ③ 4K UNIT OR 4K PLUS UNIT: 5.5 KG

SPECIFICATIONS

MAIN UNIT

DIMENSIONS: 133 MM X 430 MM X 450 MM (5.2" X 16.9" X 17.7").

UNIT WEIGHT: 13.5 KG.

CONFERENCE UNIT:

DIMENSIONS: 44 MM X 430 MM X 450 MM (1.7" X 16.9" X 17.7").

UNIT WEIGHT: 8 KG.

4K UNIT:


DIMENSIONS: 44 MM X 430 MM X 450 MM (1.7" X 16.9" X 17.7").

UNIT WEIGHT: 5.5 KG.

© [[2021]] HILLROM SERVICES, INC. ALL RIGHTS RESERVED.

NOTES:

1. INSTALLATION TO BE COMPLETED IN ACCORDANCE WITH MANUFACTURER'S SPECIFICATIONS.
2. DO NOT SCALE DRAWING.
3. THIS DRAWING IS INTENDED FOR USE BY ARCHITECTS, ENGINEERS, CONTRACTORS, CONSULTANTS AND DESIGN PROFESSIONALS FOR PLANNING PURPOSES ONLY. THIS DRAWING MAY NOT BE USED FOR CONSTRUCTION.
4. ALL INFORMATION CONTAINED HEREIN WAS CURRENT AT THE TIME OF DEVELOPMENT BUT MUST BE REVIEWED AND APPROVED BY THE PRODUCT MANUFACTURER TO BE CONSIDERED ACCURATE.
5. CONTRACTOR'S NOTE: FOR PRODUCT AND COMPANY INFORMATION VISIT www.CADdetails.com/info AND ENTER REFERENCE NUMBER 4794-126



HELION DETAILS

HELION VIDEO MANAGEMENT SYSTEM (MAIN, CONFERENCE & 4K UNITS)