


**Baxter**



## HELUX PRO SURGICAL LIGHT

THE FUTURE OF  
CONNECTED ILLUMINATION

---

Empower clinicians with a surgical lighting system that helps keep teams focused on the surgical site, while the light heads do the rest. It's time to view your OR in a new light with the **Helux Pro** Surgical Light – connected, reliable, superior illumination.



### CONSISTENT ILLUMINATION

From shadow elimination to deep cavity lighting, the **Helux Pro** provides consistent lighting conditions so your team can focus on the patient in front of them.



### ENHANCED PATIENT SAFETY

Helps reduce cross-contamination with sterile light controls and the high intensity indicator brings team awareness to surgical light safety.



### SEAMLESS ORI CAPABILITY

The surgical light and 4k camera can seamlessly integrate with ORI systems, enabling images to be routed by 3rd party equipment, providing a truly connected OR experience.



### OPTIMIZED LAMINAR AIRFLOW

The **Helux Pro** Surgical Light is developed to optimize laminar air flow, providing consistent operating conditions.

# View your OR in a new light with the Helux Pro Surgical Light



## TECHNICAL SPECIFICATIONS

Illumination level Ec in 1 m (lx)	160,000 lux
Focusable light field size (d10) in 1 m / 39.37"	6.3-9.8" (16-25 cm)
Color temperature (K)	3,500   4,000   4,500   5,000
Color rendering index (Ra)	max. 97
Dimming range (%)	< 10% Endo   30-100%
Depth of illumination (L1+L2) at 20% Ec / 2nd Edition	76" with ALC Plus (193 cm)
Shadowing with 1 shadower without tube	92%
Shadowing with 1 shadower with tube	91%
Shadowing with 2 shadowers without tube	68%
IP Rating	54

For more information please contact your local Baxter representative.

[hillrom.com](https://www.hillrom.com)

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

Baxter, Helux Pro and Hillrom are trademarks of Baxter International, Inc. or its subsidiaries. US-CS391-230003 v1 MAY-2023