



Hillrom™

Hillrom™ Stretcher Portfolio

Go Beyond Mobility

Procedural Stretcher

TECHNICAL SPECIFICATIONS

Overall Length	83" (210.8 cm)
Maximum Width (Siderails Up)	32" (813 mm) or 36" (914 mm)
High/Low Position	34.25" (87 cm)/ 20.7" (52.6 cm)
Siderail Length	47" (119.4 cm)
Siderail Height Above Sleep Deck	14.5" (368 mm)
Maximum Trend/Reverse Trend	18°
Maximum Head Elevation	90° (80° with Auto Contour™)



Customizable to be as simple or versatile as you need, our procedural stretcher accommodates all patient conditions as well as your needs in an uncertain environment.

SIMPLE & SAFE MANEUVERING

Steering Plus™ System	+
OneStep® Tuck-Away Siderail	+
BackSaver® Fowler Feature	+
Head & Foot Brake/Steer	+
Four Sides Brake/Steer	O
Dual-Sided hi/lo Hydraulic Control Pedals	+
3-Sided hi/lo Hydraulic Control Pedals	O
¾" IV Pole Sockets	+
Integrated Dual-Display Scale	O
Radiolucent Deck	O
Upright Chest X-Ray Cassette Holder	O

SAFE & COMFORTABLE EXPERIENCE

26" Width Mattress	+
30" Width Mattress	O
Auto Contour™ BackSaver® & Knee Gatch	O
Enhanced Safety Siderail	O
3" Comfortline® Foam Mattress	+
4" Comfortline® Foam Mattress	O
5" Comfortline® Foam Mattress	O
5" AccuMax Quantum™ Mattress (Elastimax or Hypolex)	O

DURABILITY TO WITHSTAND DEMANDS

Safe Working Load	700 lbs (317.5 kg)
Full Perimeter Bumpering	+
Smooth Deck Design	+

+ = Included; O = Optional

Transport Stretcher

TECHNICAL SPECIFICATIONS

Overall Length	83" (210.8 cm)
Maximum Width (Siderails Up)	32" (81.3 cm) or 36" (91.4 cm)
High/Low Position	37" (94 cm)/ 20.5" (52.1 cm)
Siderail Length	47" (119.4 cm)
Siderail Height Above Sleep Deck	14" (35.6 cm)
Maximum Trend/Reverse Trend	18°
Maximum Head Elevation	90°

SIMPLE & SAFE MANEUVERING

Steering Plus™ System	O
Easy Release Fold-Down Siderail	+
Head & Foot Brake/Steer	+
Four Sides Brake/Steer	O
Dual-Sided hi/lo Hydraulic Control Pedals	+
3-Sided hi/lo Hydraulic Control Pedals	O
¾" IV Pole Sockets	+
Radiolucent Deck	O
Upright Chest X-Ray Cassette Holder	O

Our transport stretcher is an economical solution to help accommodate all patient conditions, and meet your needs in this uncertain environment.

SAFE & COMFORTABLE EXPERIENCE

26" Width Mattress	+
30" Width Mattress	O
3" Comfortline® Foam Mattress	+
4" Comfortline® Foam Mattress	O
5" Comfortline® Foam Mattress	O
5" AccuMax Quantum™ Mattress (Elastimax or Hypolex)	O

DURABILITY TO WITHSTAND DEMANDS

Safe Working Load	700 lbs (317 kg)
Full Perimeter Bumpering	+
Smooth Deck Design	+

+ = Included; O = Optional



Trauma Stretcher

TECHNICAL SPECIFICATIONS

Overall Length	83" (210.8 cm)
Maximum Width (Siderails Up)	32" (81.3 cm)
High/Low Position	37.5" (95.3 cm)/ 24.25" (61.6 cm)
Siderail Length	47" (119.4 cm)
Siderail Height Above Sleep Deck	11" (27.9 cm)
Maximum Trend/Reverse Trend	18°
Maximum Head Elevation	90°

SIMPLE & SAFE MANEUVERING

Steering Plus™ System	+
OneStep® Tuck-Away Siderail	+
Head & Foot Brake/Steer	+
Four Sides Brake/Steer	O
Dual-Sided hi/lo Hydraulic Control Pedals	+
3-Sided hi/lo Hydraulic Control Pedals	O
¾" IV Pole Sockets	+
Integrated Dual-Display Scale	O
Radiolucent Deck	+
Full Length X-Ray Cassette Holder	+

SAFE & COMFORTABLE EXPERIENCE

26" Width Mattress	+
3" Comfortline® Foam Mattress	+
4" Comfortline® Foam Mattress	O
5" Comfortline® Foam Mattress	O
5" AccuMax Quantum™ Mattress (Elastimax or Hypolex)	O

DURABILITY TO WITHSTAND DEMANDS

Safe Working Load	700 lbs (317 kg)
Full Perimeter Bumpering	+
Smooth Deck Design	+

+ = Included; O = Optional

Our trauma stretcher provides versatility to help you care for your high acuity patient. With the built in x-ray cassette platform, capture quality images right on the stretcher, without moving your patient.



Surgical Stretcher

TECHNICAL SPECIFICATIONS

Overall Length	92" (233.7 cm)
Maximum Width (Siderails Up)	32" (81.3 cm)
High/Low Position	33.25" (84 cm)/ 21.5" (54.6 cm)
Siderail Length	47" (119.4 cm)
Siderail Height Above Sleep Deck	14.5" (36.8 cm)
Maximum Trend/Reverse Trend	18°
Maximum Head Elevation	70°

SIMPLE & SAFE MANEUVERING

Steering Plus™ System	+
OneStep® Tuck-Away Siderail	+
BackSaver® Fowler Feature	+
Head & Foot Brake/Steer	+
Four Sides Brake/Steer	O
Dual-Sided hi/lo Hydraulic Control Pedals	+
¾" IV Pole Sockets	+

SAFE & COMFORTABLE EXPERIENCE

26" Width Mattress	+
4" Surgical Surface	+
4" Flat or Concave Head Cushion	O
Manual Knee Gatch	O

DURABILITY TO WITHSTAND DEMANDS

Safe Working Load	700 lbs (317 kg)
Full Perimeter Bumpering	+
Smooth Deck Design	+

+ = Included; O = Optional

The Hillrom surgical stretcher provides flexibility to comfortably accommodate all patient conditions as well as your needs. Its simple-to-use articulating head section gives surgeons optimum positioning control to perform surgical procedures.



OB-GYN Stretcher

TECHNICAL SPECIFICATIONS

Overall Length	80.125" (203.5 cm)
Maximum Width (Siderails Up)	36" (91.4 cm)
High/Low Position	37" (94.0 cm) / 24.25" (61.6 cm)
Siderail Length	37" (94.0 cm)
Siderail Height Above Sleep Deck	13" (91.4 cm)
Maximum Trend/Reverse Trend	18°
Maximum Head Elevation	90°

SIMPLE & SAFE MANEUVERING

Steering Plus™ System	+
Easy Release Fold-Down Siderail	+
Head & Foot Brake/Steer	+
Four Sides Brake/Steer	O
Dual-Sided hi/lo Hydraulic Control Pedals	+
3-Sided hi/lo Hydraulic Control Pedals	O
¾" IV Pole Sockets	+
Radiolucent Deck	+
Upright Chest X-Ray Cassette Holder	O

SAFE & COMFORTABLE EXPERIENCE

30" Width Mattress	+
3" OB/GYN Mattress	+

DURABILITY TO WITHSTAND DEMANDS

Safe Working Load	700 lbs (317 kg)
Full Perimeter Bumpers	+
Smooth Deck Design	+

+ = Included; O = Optional

With its exclusive Stow and Go™ foot section, our OB/GYN Stretcher is multi-functional for general use to pelvic exams. With a 700 lb weight capacity, integrated catch basin, and infinite labor positioning, our stretcher can be your stand-in birthing bed when you need it.




Stretcher Accessory Compatibility

ACCESSORY	TRANSPORT	PROCEDURAL	TRAUMA	SURGICAL	OB-GYN
Chart Holder P361	X	X	X	X	X
Head/Foot Board P412T	X	X	X	X (Footboard)	X (Headboard)
Foot Extender Pad P929G1 (3"), P929G2 (4")	X	X	X	X	X
ISS Transfer Pole & Socket Adapter P158 (ISS Transfer Pole) P163 (Socket Adapter)	X	X	X	X	X
IV Transporter Device P491	X	X	X	X	X
One Button IV Pole P2217A	X	X	X	X	X
Liquid Oxygen Tank Holder P273	X	X	X	X	X
Vertical Oxygen Tank Holder P27601	X	X	X	X	X
Patient Tray P490B	X	X	X	X	X
Siderail Pads P931BT	X	X	X	X	X
Transport Strap P349	X	X	X	X	X
Ankle Stirrups P347AT	X	X	X	X	
Armboard P344CT	X	X	X		X
Head Restraint P449				X	
Scale Upgrade Kit P746A01 (26") & P746A02 (30")		X	X		
Stretcher Pillow P1425C	X	X	X		X
Integrated PACU Extenders P261ECP				X	
Wrist Rest P262A01 (Superior), P262A02 (Temporal)				X	
Telescoping Calf Supports P35745AT					X
Detachable Catch Basin P265					X
Paper Roll Dispenser P364AT01 (26"), P364AT02 (30")	X	X	X		
Upright Chest X-Ray Cassette Holder P279AT	X	X	X (built in)	X	X
Lateral X-Ray Cassette Holder P264			X		
Utility Tray P297B01 (26"), P297B02 (30")		X			

Stretcher Accessories




Patient Tray
P490B



Foot Extender Pad
P929G1 (3")
P929G2 (4")



Vertical Oxygen Tank Holder
P27601



Paper Roll Dispenser
P364AT01 (26")
P364AT02 (30")



One Button Adjustable IV Pole
P2217A



Head/Footboard
P412T



Liquid Oxygen Tank Holder
P273



Utility Tray
P297B01 (26")
P297B02 (30")




IV Transporter Device
P491




Convertible Footboard
P350CTM07



Armboard
P344CT




Transport Strap
P349




Siderail Pad
P931BT



Ankle Stirrups
P347AT



ISS Transfer Pole
P158
Socket Adapter
P163



Upright Chest X-Ray Cassette Holder
P279AT




Scale Upgrade Kit
P746A01 (26")
P746A02 (30")



Lateral X-Ray Cassette Holder
P264

Surgical Stretcher Accessories



Superior Wrist Rest
P262A01



Integrated PACU Extenders
P261ECP



Temporal Wrist Rest
P262A02



Head Restraint Strap
P449

OB/GYN Stretcher Accessories



Detachable Catch Basin
P265



Telescoping Calf Supports
P35745AT





Hillrom™

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your Hillrom sales representative at 1-800-445-3730.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR61301 rev 3 09-SEP-2021 ENG – US