



Hillrom™


Golvo® Mobile Lifts



A MOBILE LIFT THAT WORKS AS HARD AS YOU DO

Patient mobility requires a unique blend of power and flexibility. Our Golvo® mobile lifts combine the versatility of a mobile lift with overhead lifting functionality — helping you address the many lifting needs (and care environments) you encounter every day.

With its vertical lifting structure and soft, flexible lift strap, the Golvo lift is designed to lift and support patients in a variety of situations, including:

- Transfers to/from
 - Repositioning
 - Ambulation
 - Horizontal lifting
 - Lifting from the floor
 - Helping patients out of their personal vehicles
- 



MOBILE LIFT FLEXIBILITY MEETS OVERHEAD FUNCTIONALITY



**Golvo® 9000
Mobile Lift**

Golvo® 9000 LowBase Mobile Lift

LowBase option gives extra clearance in tight spaces (e.g., beneath low beds and stretchers).

Safe Working Load: 200 kg / 440 lbs

Please note: the base width may be adjusted with weight up to 140 kg / 308 lbs.



ADJUSTS TO YOUR ENVIRONMENT

The lift's adjustable-width base makes it easy to fit into tight spaces like chairs, toilets and low beds.



ADDED ARMREST SUPPORT

Retractable armrests give patients support and confidence during ambulation, and make it easier for caregivers to maneuver the lift.



FLEXIBLE LIFTSTRAP

The soft, flexible liftstrap helps reach patients where they are—even on the floor or inside their personal vehicles.

EASY TO USE

- Armrests come standard to support ambulation and easy maneuvering
- Lightweight, environmentally friendly lead or lithium-ion battery is easy to charge and easy to handle

VERSATILE

- Flexible lifting interval is ideal for ambulation and lifting from the floor
- Optional quick-release system makes changing the slingbar or installing a scale quick and easy
- Variety of slings and soft accessories help you customize the lift for each situation

FUNCTIONAL

- Flexible lift strap helps you reach patients in many positions, and gives patients and caregivers space and comfort
- Intelligent diagnostic system shows caregivers when the lift is safe to use and when it's time for service
- Parallel base-width adjustment supports lifts from the floor or from chairs
- Aluminum design is strong, durable and environmentally friendly

SUPPORT THAT GOES WHERE YOU GO

With its unique lift strap design and flexible lifting interval, the Golvo mobile lift is ideal for many different care environments—including rehabilitation clinics, ICUs, childcare centers and operating rooms.



Ambulation

The Golvo lift's retractable armrests, open design and generous lifting height help patients take those first steps safely and confidently.



Lifting from the Floor

Lifting to and from the floor is a common need in childcare and rehabilitation settings. Our flexible, comfortable slings support patients' heads throughout the entire lifting operation.



Lifting to/from Bed

The Golvo's flexible lift strap and parallel base-width adjustment make lifting to and from bed smooth and easy. For very low beds and stretchers, the LowBase option provides extra clearance.



Repositioning and Treatment

The RepoSheet® positioning aid attaches to help turn patients in bed with minimal strain on the body. The MultiStrap™ lifting accessories support wound care, hygiene procedures and treatment.



Lifting to/from Chairs and Commodes

The Golvo lift base adjusts electronically, opening in a parallel and symmetrical motion. This provides excellent access to chairs, and the flexible lift strap makes positioning simple.



Vehicle Extraction

The flexible strap extends into vehicles, making it easier to apply the sling loops onto the slingbar. Then the Golvo lift supports the patient's weight while caregivers guide them out of the vehicle.

Golvo® Mobile Lifts Technical Specifications

	Golvo 9000	Golvo 9000 LowBase
Product Number	2000045	2000049
Maximum Load (With Base Width Adjustment)	200 kg / 440 lbs (140 kg / 308 lbs)	200 kg / 440 lbs (140 kg / 308 lbs)
Weight With Lead Battery With Lithium Ion Battery	41.6 kg / 91.7 lbs 40.2 kg / 88.63 lbs	45 kg / 99.2 lbs 43.6 kg / 96.1 lbs
Height	1455 – 2090 mm / 57.3 – 82.3"	1455 - 2090 mm / 57.3 – 82.3"
Base Height Length Width Opening of the Base	105 mm / 4.1" 1185 mm / 46.6" 735 – 1020 mm / 28.9 – 40.2" Parallel	105 mm / 4.1" 1185 mm / 46.6" 735 – 1020 mm / 28.9 – 40.2" Parallel
Casters Front Rear	75 mm / 3" dual 75 mm / 3" dual	46 mm / 1.8" dual 75 mm / 3" dual
Lifting Speed (Without Load)	27 mm/sec (1.1"/sec) or 37 mm/sec (1.5"/sec)	27 mm/sec (1.1"/sec) or 37 mm/sec (1.5"/sec)
Maximum Lifting Height	1816 mm / 71.5"	1816 mm / 71.5"
Lifting Interval	1245 mm / 49"	1245 mm / 49"
Emergency Lifting	Electric	Electric
Emergency Lowering	Electric and manual	Electric and manual
Sound Level	52.8 dB(A)	52.8 dB(A)
Battery Option	Lead acid gel or lithium ion	Lead acid gel or lithium ion
Charger	Built-in	Built-in
Hand Control	Wired	Wired
Protection Class (including hand control)	IPX4	IPX4
Material	Aluminum	Aluminum
Slingbar Standard Optional	Universal 450 Slingbar Universal 350 Slingbar Universal 600 Slingbar Universal Twinbar 670	Universal 450 Slingbar Universal 350 Slingbar Universal 600 Slingbar Universal Twinbar 670
Flexlink or Flexible Liftstrap	Soft flexible liftstrap	Soft flexible liftstrap
Accessories	Quick-Release Hook LikoScale™ 350, 200*, 400* Horizontal lifting device: OctoStretch™ FlexoStretch™ LikoStretch Mod 600 IC™ Slingbar Cover Paddy	Quick-Release Hook LikoScale™ 350, 200*, 400* Horizontal lifting device: OctoStretch™ FlexoStretch™ LikoStretch Mod 600 IC™ Slingbar Cover Paddy

* Available in the United States and Canada only.

**IF YOU KNOW
LIKO, YOU KNOW
HILLROM.**

**WE'VE UNIFIED
OUR BRANDS
TO ADVANCE
CONNECTED
CARE. VISIT US AT
HILLROM.COM**



Hillrom™

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your Hillrom sales representative at 1-800-445-3730.

hillrom.com

Golvo™ mobile lift is intended for transferring patients (adult or children) between devices (e.g., within the room), floor lifting, horizontal lifting, supporting patient limbs, ambulating patient, bathing patient, toileting patient, weighing patient and transferring patients from a car. Intended for use in the following environments: healthcare, intensive care, emergency ward, rehabilitation, habilitation.

Class 1

Manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå, Sweden

These products are regulated health care products which, where required by applicable regulations, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documentation accompanying medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Liko R&D AB. ALL RIGHTS RESERVED. APR296701 rev 1 10-NOV-2021 ENG – US