

Hill-Rom® Bariatric Solutions

Transforming bariatric care.



Enhancing outcomes for
patients and their caregivers:

Hill-Rom.

Compella™ Bariatric Bed

Advancing bariatric care.

The Compella™ Bariatric Bed helps overcome the barriers to delivering efficient, safe and dignified bariatric care.



Streamlining Workflow

- IntelliDrive® XL powered transport
- Exclusive Powered width expansion from 40"-50"
- Adjustable bed length

Delivering Safe Patient Care

- 3-mode bed exit with silence feature
- CLRT and 30° head-of-bed alert
- Turn Assist
- Pressure redistribution air surface

Enhancing Patient Dignity

- Aesthetic design
- Side sitting comfort
- Accessible controls
- Patient weight capacity: 1,000 lbs

TotalCare® Bariatric Bed

Innovative bariatric care.

TotalCare® Bariatric Bed delivers innovative support with features to help mobilize bariatric patients without compromising dignity.



Mobilize with Confidence

- CLRT and P&V
- 30° / 45° head-of-bed alert
- FullChair® and Chair egress positioning

Deliver Efficient Care

- Adjustable bed length
- Integrated scale
- Battery backup
- IntelliDrive® powered transport

Protect Patients and Caregivers

- 40" pressure redistribution surface
- Advanced Microclimate® Technology
- Turn Assist
- Bed exit
- Patient weight capacity: 500 lbs

Excel Care® ES Bariatric Bed

Adaptable bariatric care.

The Excel Care® ES Bariatric Bed offers the flexibility to help deliver exceptional care for a range of bariatric patient needs.



Flexible Patient Support

- Pressure redistribution air surface
- Surface remains inflated during transport
- Turn Assist
- CLRT
- Head-of-bed angle indicator

Simplified Care

- Expandable frame & surface (40"-50")
- Adjustable bed length (80"-88")
- Integrated scale
- Battery backup

Supports Patient Independence

- Single post patient-helper trapeze
- Seat deflate
- Patient weight capacity: 1,000 lbs

Tri-Flex II™ Bariatric Bed*

Dependable bariatric care.

The Tri-Flex II™ Bariatric Bed is a reliable and dependable choice for the care of bariatric patients.



**Rental only, 37"-54" capital*

Dependable Care

- Expandable frame (37"-48")
- Integrated scale
- Battery backup
- Easy to use caregiver controls

Reliability

- Low air loss and alternating pressure surface (Synergy® Air Elite)
- Low deck height of 14"

Patient Comfort

- Patient-helper trapeze
- Patient weight capacity: 1,000 lbs

Sentinel II Bariatric Low Bed

Sentinel II Bariatric low bed provides an integrated solution with features to help prevent patient falls and help mobilize bariatric patients.



Efficiency

- Expandable frame (38"-48")
- Adjustable bed length (80"-84"-88")
- Integrated scale
- One-button chair positioning

Fall Prevention

- Low bed height of 10.88"
- Bed exit alarm
- Connection to nurse call

Patient Dignity

- Patient-helper trapeze
- Patient weight capacity: 850 lbs

Bariatric Accessories*

Hill-Rom's bariatric portfolio delivers a range of products that work together to help perform routine tasks more efficiently and provide a comfortable experience for both patients and their families.



Bariatric Expandable Walker

- Adjustable width and height to accommodate various patient sizes.
- Patient weight capacity: 1,000 lbs

Bariatric Wheelchair

- Swing-away elevating leg rests allow for easier patient ingress and egress.
- Patient weight capacity: 700 lbs

Extra Large Commode

- Easy-to-release drop arm allows you to make safer lateral transfers.
- Patient weight capacity: 1,000 lbs

Shower Transfer Bench

- Durable, molded plastic bench and backrest provides adjustable height of individual zones.
- Patient weight capacity: 600 lbs

**Bariatric Accessories only available as rental products.*

Procedural Stretcher

- A powered wheel deployment lets you easily step on the pedal and go.
- Ergonomic push handles with integrated IV pole cradle supports easy transport positioning.
- Available in 30" width
- Patient weight capacity: 700 lbs



Patient Transfer Devices

AIRPAL® Patient Transfer System

- Lateral transfer system that gently glides patients across a friction-reducing air surface for a smooth transition.
- Patient Weight Capacity: 1,200 lbs



Barton® Transfer Chair

- Converts between a flat and upright mobile chair that allows a single caregiver to easily transfer patients.
- Patient weight capacity:
660 lbs (rental);
700 lbs (capital)



Patient Lift Devices

Viking® XL Mobile Floor Lift UltraTwin™ Overhead Lift

- Lightweight, easy to maneuver
- Versatile—supports ambulation, lifting from bed or floor, etc.
- Patient weight capacity:
660 lbs (Mobile Lift);
1,100 lbs (Overhead Lift)



Furniture

Art of Care® Bariatric Recliner

- Three-position seating
- 30" seat
- Front and rear locking casters
- Accessory hooks, optional removable arm and tray table
- Powered Bariatric Recliner available
- Patient weight capacity: 660 lbs



Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world, Hill-Rom's people, products, and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2019 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

CTG711 rev 13 11-MAR-2019 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730

Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom.