

**Baxter**

PROMOTIONAL OFFER

# UPGRADE AND SAVE

TRADE UP YOUR LEGACY VITAL  
SIGNS DEVICE AND UPGRADE  
TO A NEW QUALIFYING DEVICE\*  
TO RECEIVE A \$200 REBATE

**OCTOBER 1-DECEMBER 31, 2024**

**\$200**

The **Spot Vital Signs 4200**  
has been end of manufactured  
and is scheduled to go **End of  
Service June 30, 2025**.

The **Spot Vital Signs LXi** went  
**End of Service June 30, 2023**.

Avoid unnecessary downtime by  
upgrading to a new Welch Allyn unit today.

**Contact your Baxter representative to learn more  
or click the button to redeem this offer.**

\*Connex Spot Monitor, Connex Vital Signs Monitor, Connex Integrated Wall System,  
Spot Vital Signs 4400 device.



# UPGRADE YOUR DEVICE TODAY

## COMPARISON OF OUR CURRENT SOLUTIONS VS. LEGACY DEVICES

FEATURE	SPOT VITAL SIGNS 4200 DEVICE <i>(END OF MANUFACTURE)</i>	SPOT VITAL SIGNS LXi DEVICE <i>(END OF SERVICE)</i>	SPOT VITAL SIGNS 4400 DEVICE	CONNEX SPOT MONITOR
Workflows	Spot	Spot	Spot	Spot, Office, and Interval Profiles
BP Averaging	No	No	Yes	Yes
Parameters	Step BP, Temp, SPO <sub>2</sub>	<b>SureBP</b> , Temp, SPO <sub>2</sub>	<b>SureBP</b> , Temp, SPO <sub>2</sub>	<b>SureBP</b> , Temp, SPO <sub>2</sub> , RRp
SPO <sub>2</sub> Options	<b>Nellcor, Masimo</b>	<b>Nellcor, Masimo</b>	<b>Nonin</b>	<b>Nellcor, Masimo, Nonin</b>
Temperature Options	<b>SureTemp</b>	<b>SureTemp Plus, Braun PRO 4000</b>	<b>SureTemp Plus</b>	<b>SureTemp Plus, Braun PRO 6000</b>
Custom Parameters/Modifiers	None	Basic (height, weight, pain, RR, BMI)	Basic (height, weight, pain, RR, BMI)	Advanced (20 customizable fields)
Early Warning Scores	No	No	No	Yes
Wireless Connectivity	No	<b>Wi-Fi®</b> technology via Accessory	No	<b>Wi-Fi</b> technology, <b>Bluetooth®</b> technology
Wired Connectivity	Serial Cable	USB, Serial Cable	USB	USB or Ethernet
Barcode Scanner	No	Yes	No	Yes
Scale Connectivity	No	Yes via USB Dongle	No	Yes via USB Dongle
Display	5" Monochrome LCD	5.8" Monochrome LCD	7" Color Touchscreen	7" Color Touchscreen
Memory	1 Record	50 records	50 records	400 records
Battery Type	Li Ion	Li Ion	Li Ion	Li Ion
Charge Time	12 Hours	10 Hours	4 Hours	4 Hours
Run Time	8 Hours	8 Hours	8 Hours	8 Hours or 17 Hours with APM
Mounting Options	Basic Stand, Wall Mount, IV Pole Mount	Basic Stand, Wall Mount	Desktop, Basic Stand, 777 Wall Mount	Desktop, Basic Stand, APM Stand, 777 or GCX Wall Mount

**Rx only.** For safe and proper use of the products mentioned herein, please refer to the Operator's Manual or Instructions for Use.

Offer available on eligible purchases made October 1–December 31, 2024. To be eligible for offer, customer must purchase qualifying Welch Allyn products identified in this document and provide proof of Welch Allyn purchase as indicated on the redemption form found at the links above. All information must be submitted via the online redemption process by December 31, 2024. All customers must submit a completed W-9 for the vital signs rebate offer only. Promotion valid in U.S. only to primary care medical facilities, Long Term Care Facilities, and Ambulatory Surgery Centers. Distributor stocking orders do not qualify. Please allow 12 weeks for receipt of vital signs rebate check based on terms above. Promotion cannot be combined with any other promotion, offer, coupon, enhanced discount or special pricing arrangement; some exclusions may apply.

Baxter, Connex, Spot Vital Signs, SureTemp, and Welch Allyn are trademarks of Baxter International Inc or its subsidiaries.

Bluetooth is a trademark of Bluetooth SIG, Inc.

Braun is a trademark of Braun GmbH.

Wi-Fi is a trademark of Wi-Fi Alliance.

US-FLC157-240041 (V1.0) 09/2024