



Welch Allyn  
**SureTemp Plus**  
Electronic Thermometers



## WELCH ALLYN SURETEMP PLUS 690 ELECTRONIC THERMOMETER

- Removable probe well for easy cleaning
- Large, easy-to-read LCD
- C/F button for toggling between Fahrenheit and Celsius
- Recall button displays last temperature taken
- Rubberized grips
- Convenient storage for 25 probe covers
- Battery life indicator
- Mode button



### OPTIONAL WALL MOUNT

The SureTemp Plus 690 with optional wall mount has storage for an additional 50 probe covers.

## WELCH ALLYN SURETEMP PLUS 692 ELECTRONIC THERMOMETER

SAME FEATURES AS THE 690 MODEL, PLUS:

- NEW** Backlit display helps in low-light environments
- Pulse timer button/monitor mode

**NEW FEATURE** Easier to switch to Monitor Mode for measuring a patient's thermal steady state (approx. 3 to 5 minutes)



### STANDARD WALL MOUNT

The SureTemp Plus 692 with standard wall mount has storage for an additional 50 probe covers.

Interchangeable rectal probe and well reduce risk of cross-contamination

### IDEAL CARE SETTINGS:

- Family Practice
- Pediatrics
- Long-term Care
- Internal Medicine
- OB/Gyn



### IDEAL CARE SETTINGS:

- Urgent Care
- MedSurge
- Maternity
- Emergency Medicine
- Oncology
- Cardiology

# SETTING THE STANDARD OF CARE

Enhance a routine task with versatile, easy and accurate technology featuring oral, axillary and rectal modes to meet all your temperature taking needs.

## 1 EASY

An easy-to-read LCD display with flashing icons, and recall button displays last temperature taken. SureTemp Plus 692 model features a 60-second timer to support pulse and respiration readings.

## 3 FAST AND ACCURATE

The SureTemp Thermometer captures quick readings - approximately 4-6 seconds oral, 10-13 second pediatric axillary and rectal readings and 12-15 second adult axillary readings. Toggle between Celsius and Fahrenheit, or monitor mode that helps measure thermal steady state.

## 2 SAFE AND SECURE

Soft, single-use disposable covers help minimize cross-contamination risk, and interchangeable probe wells are easy to remove and clean as you switch between modes. The SureTemp Plus 692 model is designed with two security features to help protect against theft.



User-selectable icon and/or words for identification of thermometer mode

Battery life indicator

## TECHNICAL SPECIFICATIONS

Calibration Accuracy	±0.2 °F (±0.1 °C) (Monitor Mode)
Temperature Predict Time	<b>Oral:</b> Approximately 4-6 seconds <b>Adult Axillary (18 years and older):</b> Approximately 12 - 15 seconds <b>Pediatric Axillary (17 years and younger):</b> Approximately 10 - 13 seconds <b>Rectal:</b> Approximately 10 - 13 seconds
Unit Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)
Weight	<b>Model 690:</b> 12.6 oz (357 g), <b>Model 692:</b> 12.6 oz (357 g)
Ambient Temperature Range	50.0 °F (10 °C) to 104 °F (40.0 °C)
Patient Temperature Range	80.0 °F (26.7 °C) to 110 °F (43.3 °C)
Power Source	(3 ) 1.5 Vdc AA batteries
Battery Life	<b>Model 690:</b> Approximately 6,000 temperature, <b>Model 692:</b> Approximately 5,000 temperatures

## PRODUCT AND ACCESSORY ORDERING INFORMATION

### THERMOMETERS

#### 690 MODELS:

**01692-200:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 4ft cord

**01692-201:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 4ft cord

**01692-300:** SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord (must be locked to the wall)

**01692-301:** SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Rectal Probe and Rectal Probe Well with 9ft cord (must be locked to the wall)

**01692-700:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord, with Rolling Stand and Holder for Spare Probe & Well

#### 692 MODELS:

**01692-200:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 4ft cord

**01692-201:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 4ft cord

**01692-300:** SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord (must be locked to the wall)

**01692-301:** SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Rectal Probe and Rectal Probe Well with 9ft cord (must be locked to the wall)

**01692-700:** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord, with Rolling Stand and Holder for Spare Probe & Well

### ACCESSORIES

**21326-0000:** Wall Holder

**21330-0000:** Holder for Spare Probe and Well

**06138-000:** Calibration Key Assembly, 690/692

**02892-000:** Probe and Well Kit, 4ft, Rectal

**02892-100:** Probe and Well Kit, 9ft, Rectal

**02893-000:** Probe and Well Kit, 4ft, Oral

**02893-100:** Probe and Well Kit, 9ft, Oral

**02891-0000:** Probe Well, Blue

**02891-1000:** Probe Well, Red

**08273-000:** Thermometer Stand with Holder for Spare Probe and Well

**21333-0000:** Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well

**98025-0000:** Extended Warranty M692

### REPLACEMENT PROBE COVERS

**05031-101:** Case of 1,000 Probe Covers

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

**hillrom.com**