

Synergy[®] Air Elite Low Air Loss Therapy

A flexible choice for low air loss and alternating pressure.



Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Synergy® Air Elite Low Air Loss Therapy



The flexible choice for wound prevention and healing.

With both Low Air Loss and alternating pressure treatment, the Synergy® Air Elite Low Air Loss Therapy allows you, as a caregiver, to create customized therapy and prevention programs for your patients. This mattress replacement helps enable the cost-conscious hospital to implement therapy and prevention programs in accordance

with Joint Commission standards, without the expense of purchasing an entire new fleet of beds. The Synergy® Air Elite surface helps enhance the level of care you provide your patients while helping improve the patient safety standards in your facility.



Hill-Rom
A HILLENBRAND COMPANY

Key Features at a Glance

- Alternating pressure and low air loss mattress replacement system:
 - Ideally suited for patients in a variety of care settings
 - Polyurethane fitted cover reduces effects of shear and friction
 - Utilizes existing bed frames, reducing rental costs and storage issues
- Aids in the prevention and treatment of pressure ulcers
- Fully inflates in less than two minutes
- Air flow addresses maceration related problems and promotes a dry healing environment
- Available for rental or purchase

Synergy[®] Air Elite surface is available for rent on Hill-Rom[®] beds such as the Advanta[™] 2 frame and can be used in conjunction with other flat deck beds.*

* Confirm compatibility with frame manufacturer.



Additional Features and Benefits

Performance:

- Three therapy choices:
 - Low air loss
 - Alternating pressure/pressure relief
 - Pressure relief
- Flexibility to customize therapy for prevention and healing

Safety:

- Auto firm function provides uniform firmness to accommodate caregiver procedures
- Safety features:
 - Power failure triggers auto alarm
 - Main power switch can be shut off in an emergency
 - CPR hose disconnects for rapid deflation

Comfort and Convenience:

- Individual air cushions:
 - Connected to side panel with connector for easy detachment and cleaning
- Double insulated blower provides quiet operation



The average cost per hospitalization for patients who develop a Stage III & IV pressure ulcers has been reported to be \$43,180.*

Partnering for Wound Prevention and Treatment

Technical Specifications

Control Unit

Inflation time	under 2 minutes
Functions	alternating pressure, pressure relief, max inflate, low air loss
Weight	13.2 lbs
Alternating pressure	5-95 minutes
Max flow rate	1200 lpm
Rated voltage	110-120V
Max current	5 AMPS
Rated frequency	60 Hz
Protection type	Class II, Type BF
Mode of operation	continuous

Mattress Material

Latex-free comforter	polyurethane-coated taffeta nylon
Cell	polyurethane-laminated taffeta nylon
Base	high denier polyester
Dimensions (H" x W" x L")	
Inflated	10" x 36" x 85" 10" x 39" x 85" 10" x 48" x 85"
Suggested Medicare code	E0277
Product weight	33 lbs (Standard) 40 lbs (Bariatric)
Weight range	36" Width 600 lbs 39" Width 1,000 lbs 48" Width 1,000 lbs

* Center for Medicare & Medicaid Services Office of Public Affairs, April 14, 2008, Fact Sheet for: CMS PROPOSES ADDITIONS TO LIST OF HOSPITAL-ACQUIRED CONDITIONS FOR FISCAL YEAR 2009.

Hill-Rom is a leading global medical technology company with more than 7,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom's people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.®**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
153727 rev 3 28-OCT-2014 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom